



SHINING 3D



Aoralscan Elite

3.5.2

User Manual

Table of Contents

Overview

Foreword

Getting Started

- Hardware introduction
- Software functions

Hardware

Unpack the Package

Hardware Overview

Connect/Disconnect the Scanner

- Attach and Detach the scanner tip
- Connect and disconnect the scanner

Specifications

Environmental Requirements

About IntraoralScan

System Requirements

Register and Log in

Activation

- Activation for the first time
- Authorization renewal

Upgrade Software

Calibrate the Scanner

Create Order

Interface

- Preview
- ① Home and dental cloud connection
- ② Settings
- ③ Create/ Import Orders
- ④ Advertisement
- ⑤ Orders

Order Settings

Create Order

- Order Information
- Tooth Selection
- Import Orders

Patient Merge

- Patient Merge in Software
- Merge Patients of Imported Order

Scan

Scanning Preparations

- Intraoral Environment
- Scanner Preparation
- Scanning Position and Path
- Heat the Scanner Tip

Scanning Settings

Scan Interface Introduction

- Function preview
- ① Camera Window
- ② Guide
- ③ IMU menu/Intraoral scan/Edentulous scan
- ④ Progress Bar
- ⑤ Preview
- ⑥ Operations
- ⑦ Scan frames and time
- ⑧ Extra match
- ⑨ Tools and additional functions

Functions

Scan Workflow

Pre-op Scanning

- Create an Order
- Start Scanning
- Scanning Completed

Restoration Scan

- Create an Order
- Start Scanning

- Scanning Completed

Implant Scan

- Create an order
- Start Scanning
- Scanning Completed

Dynamic Bite

Removable Denture Scan

- Create an order
- Natural teeth bite
- Bite rim scanning
- Old denture scanning

Extra Match

- Create extra align
- Other functions

Multiple Bites

- Add multiple bites
- Pre-design
- Send

Coded Scanbody Scanning

- Create an order
- Import an order
- Scan

Scanning Tools and Additional Functions

- Scanning tools
- Additional functions

Send Order

Pre-design

Preview Edit

- Intraoral data
- Application list
- Fill Holes
- Fit View

Oral Health Report

- Create an oral report
- Edit the oral report

Coordinate Adjustment

Mark Teeth

Bite Detection

Extract Margin Line

Modified Model

- Button
- Buttons

Dynamic Bite

ConsulOS

Overview

- Features
- Quick guide

Create program

- Model adjustment
- Teeth confirm
- Program confirm
- Change program

Orthodontic simulation

- Generate report
- Animation
- Compare the programs
- Other Functions

Manual setup

Add brackets

AccuDesign

Introduction

Import the Model

- Create a New Project
- Import a Project

Coord Adjustment

Area Adjustment

Base Adjustment

Attachment Adjustment

- Articulator Adjustment
- Text Adjustment

Printing View

Export the Model

- MetronTrack

What is MetronTrack?

Characteristic

Segment

- Interface

Measurement Workflow

Measurement Guidance

- Measurement Items
- Measurement Template

Measurement Report

Model comparison

- Temp crown
- CreIBT

Introduction

Workflow

Settings

- Create a New Template

Segmentation

- Interface
- Functions
- Options

Arrangement

- Interface
- Operations
- Options

FA points

- Interface
- Operation
- Toolbar
- Options

Bracket

- [Interface](#)
- [Bracket selection](#)
- [Bracket editing](#)
- [Toolbar](#)
- [Options](#)

Guide

- [Interface](#)
- [Guide create](#)
- [Guide edit](#)
- [Options](#)
- [Export model with brackets](#)

Report

- [Functions](#)

- [CreSplint](#)

Splint Introduction

Characteristic

Settings

Data edit

- [Model list](#)
- [Path insertion](#)
- [Edge line](#)
- [Parameter](#)
- [Visualization tool](#)

- [AccuWare](#)
- [Tools](#)

Introduction to Dental Cloud

Use ScanBinder

Use DentalViewer

- [Care and Maintenance](#)

Pre-cleaning, Disinfection, and Sterilization

- [Scanner Body, Cradle and Calibrator Maintenance](#)
- [Scanner Tip Maintenance](#)

Scanner Storage




- [Contact Us](#)

Contact Us

Overview

Foreword

Safety Instruction

Signal	Meaning
	Note: This symbol is used to inform you of the additional information of the product.
	Caution: This symbol is used to inform you of incorrect operations that may damage the device or result in data loss. Any damages resulting from misuse are not covered by the warranty.
	Warning: This symbol is used to inform you of the potential risks that may result in serious personal injury and other safety incidents.

The Declaration of Intellectual Property and Disclaimer

Thank you for using the products of SHINING 3D TECH CO., LTD.(hereinafter referred to as the "SHINING 3D"). Before you use the products, please carefully read and understand this declaration. Once you use this product, it means that you fully accept this statement and promise to comply with the relevant regulations.

1. The contents of the Product Instruction and User Manual (hereinafter collectively referred to as the "Product Usage Documentation") are critical to your personal safety, legal rights, and liabilities. Before you use the products, Please ensure that you have carefully read the Product Usage Documentation, and use the product correctly in accordance with the requirements of the Product Usage Documentation. We also recommend that the products be operated by trained professional technicians.
2. Please inspect and/or maintain the product before use. If the product is damaged, deformed or in any other abnormal condition, stop using it immediately and contact the after-sales service personnel for maintenance. SHINING 3D will not be responsible for any problems caused by your failure to inspect or maintain the product in a timely manner.
3. SHINING 3D does not guarantee the applicability of the outcomes of your use of the products, and you are responsible for verifying the quality and functionality of the outcomes. You should check and verify thoroughly that any outcomes meet your requirements before using them, for which you bear full responsibility. If any damage arising from using the outcomes of any products, you shall bear the corresponding risk, and SHINING 3D shall not bear any responsibility.
4. SHINING 3D owns complete intellectual property rights for the contents of the for which you bear full responsibility. Without the written consent of SHINING 3D, it is not allowed to copy, transmit, publish, adapt, compile or translate any contents of the Product Usage Documentation in any form for any purpose.

5. The Product Usage Documentation is a guidance for installing, operating, and maintaining the product instead of serving as the quality guaranty for the products. SHINING 3D makes all efforts to ensure the applicability of the Product Usage Documentation, but reserves the right of final interpretation. Images and diagrams in the product documentation are presented to provide convenience to user understanding. In the event that any images or diagrams are inconsistent with the physical products, the later shall prevail. In addition to the mandatory provisions of laws and regulations, the contents of the Product Usage Documentation are subject to changes without further notice.
6. SHINING 3D shall not be held responsible for any damages and/or losses caused by human factors, environmental factors, improper storage and use, or any other factors other than due to the quality of the product. SHINING 3D also shall not be held responsible for any indirect anticipated profit loss, loss of reputation and other indirect economic losses. Except as otherwise expressly provided by laws and regulations, the total liability assumed by SHINING 3D (regardless of cause) shall not exceed the purchase price of the products you paid to SHINING 3D.
7. Disputes arising from this Declaration and the Product Usage Documentation thereof shall be governed by the laws of the People's Republic of China, excluding its conflict of law rules. In the event that certain provisions are in conflict with the applicable law, these provisions will be reinterpreted in full accordance with the law, while other valid provisions will remain in force.
8. All disputes between you and SHINING 3D that arise from, shall first be resolved amicably through negotiation. If a dispute cannot be resolved through friendly negotiation, any party may submit the dispute to the Court of Xiaoshan District, Hangzhou City, Zhejiang Province, People's Republic of China for litigation and settlement.
9. In the event of any questions about the contents of this Declaration and application of Product Usage Documentation, please contact us by the contact information provided in the User Manual. Thank you for your cooperation and support! We hope that our products can bring you a great experience of using.

Getting Started

This manual mainly introduces the [hardware](#) and [software functions](#) of AoralScan Elite.

AoralScan Elite is designed to provide powder-free intraoral color scanning with higher speed, bringing greater accuracy and less time-lag for image acquisition. It can be used to scan a single tooth, multiple teeth, and whole dental arches. The captured 3D digital images of teeth and soft-tissue areas are designed to be used in conjunction with the supplied software programs.

Dental Order System module helps manage the patient information and scanned records. Scan module assists you in acquiring digital images. Exporting scanning data (in STL/OBJ format) to CAD/CAM systems for different purposes of dental care is also supported.



Hardware introduction

Hardware introduction describes different parts of the scanner and how to connect the scanner and computer, etc.

Connection

→ [How to connect/disconnect the device?](#)

Scanner description

→ [Check the package list.](#)

→ [Check the description of different parts of the scanner.](#)

→ [How to clean the scanner?](#)

→ [How to store the scanner?](#)

Software functions

The recommended process of scanning by the supplied software IntraoralScan is: Calibrate(Optional) → Create an order → Scan → Pre-design → Send the order.

IntraoralScan (Downloaded from the USB flash driver)

→ [How to register an account or log in?](#)

→ [How to activate the device?](#)

→ [How to upgrade the IntraoralScan?](#)

Calibration

Order system

Calibrate the scanner to ensure its accuracy and the quality of scan data. It contains two functions: brightness adjustment and calibration.

- [Precautions before calibration](#)
- [How to calibrate the scanner?](#)

Create, edit, search and delete orders., as well as add new patients and so on.

- [Functions](#)
- [How to check my orders and my patients?](#)
- [How to create a new order?](#)

Scanning

Scan the object and collect digital images, and export 3D data in STL, OBJ, BEB, PLY format.

- [Precautions before scanning](#)
- [How to set scanning parameters?](#)
- [Functions during scanning](#)
- [How to start scanning? \(Take pro-op scanning as an example\)](#)

Pre-design

Various functions are provided for the users to preview and edit the 3D model.

- [How to edit a 3D model?](#)
- [Other functions](#)
- [How to create my oral report?](#)
- [How to extract the margin line?](#)


Hardware

Unpack the Package

Check the carry box for the following items. If any item is missing or damaged, contact the distributor or service provider immediately.

Note

The following figures in the parts list are for reference only, the actual product shall prevail if there is any inconsistency.



No.	Name	No.	Name
1	Scanner body	2	Calibrator
3	Connection cable	4	Cradle
5	Scanner tips: standard tips, big tip (optional) and mini tip(optional)	6	Dust cap
7	Dust cover	8	USB flash driver
9	Quick guide and certification	10	Disinfection guide A and disinfection guide B
11	Package box		

Notes


- The number of scanner tips is selectable.
- The mini tip is intended to be used for patients who can not fully open their mouth, or for scanning narrow areas in the mouth (such as the back teeth area, etc.). Please do not use the mini tip for other purposes.

Caution

- AC plug types vary by country/region.
- Using accessories, peripherals, or cables not supplied with the product or recommended by SHINING 3D can affect the device in the form of increased emissions or decreased immunity to external EMI/EMC occurrences. Non-specified peripherals, and cables in some cases, can also increase leakage current or compromise the safety of the grounding scheme.
- Using accessories or power supply units other than those specified may cause the warranty to void and result in increased emissions, decreased EMI immunity of the device, or even damages to the device and personal injuries.
- Place the USB flash drive in a safe place for later usage.
- We recommend that you keep all the original packaging components in a safe place in case you need to transport or dispose of the scanner in the future. More details see [Care and Maintenance](#).

Hardware Overview


Scanner Tip and Scanner Body



No.	Name	Description
1	Scanner tip	Use the scanner tip to scan the upper, lower or full jaw. The scanner tip can be used for 100 times. The scanner can identify different types of tips and automatically adjust the format.
2	Heater and laser window	The heating device ensure successful scanning by preventing fogging on the mirror. The laser window.
3	Scan button	Single press to start scanning and pause scanning; long press for about 3 seconds to proceed to the next step; long press for about 7 seconds to turn off the scanner.
4	Scanner body	Rotate the scanner body during scanning to obtain the best scanning angle. During the scanning process, the scanner body may heat up, but the temperature will not cause harm to users and patients.
5	Indicator	Indicates the status of the scanner. <ul style="list-style-type: none"> • Green: The scanner is in scanning, heating or standby status. • Breathing green: The scanner is in standby mode, low battery or unconnected. • Orange: Abnormal status. The scanner is not correctly connected or malfunctioned. The scanner tip is not inserted tightly or the scanner is not correctly connected. • Extinguished: 1. The scanner is in the sleep mode. 2. There is no power supply. 3. The scanner is shut down.

Scanner Cradle





 **Note**

- When the scanner is not in use, place it on the cradle.
- If the scanner is placed on the cradle, it will automatically enter the standby mode. When the scanner is not used and placed on the cradle for more than 10 minutes, it will enter the sleep mode.
- As long as there is power connected to scanner, the scanner tip will be heated even if the scanner is in standby or sleep mode.


Scanner USB cable storage

To prevent the USB cable from getting damaged by excessive bending or twisting, you should loosely coil the cable and avoid making kinks or sharp bends. Up to 30,000 uses of the scanner tail cable.



 **Caution**

- Do not roll the cable over the handle of the scanner or even bend the cable sharply. The illustration below demonstrates improper cable storage.
- Do not drag the tail cable of the device when taking the scanner to avoid damage caused by excessive pulling.



Connect/Disconnect the Scanner

Caution

- Please follow the equipment connection diagram when connecting the devices and aware the sequence of cable interface connections. Do not use non-kit cables for connection.
- To ensure equipment performance and safety, use only the original accessories and software provided with the device.
- Using cables or accessories other than those provided in the package might result in equipment malfunctions, damages, or financial losses.
- Ensure the supplied software programs are installed on your computer before the connection.
- If the accuracy of the equipment decreases or if the equipment does not work properly, please consult technical support promptly.
- Install the scanner in accordance with the instructions stated in the Manual.
- Use the scanner only in dental laboratories, dental clinics, and equivalent environment.
- Do not install, place, and use the scanner in dusty and damp environment or in the areas of temperature extremes or in direct sunlight.
- Prepare a flat surface, e.g. your desk, for the scanner and the cradle. Do not place them on a slanted surface.
- Before the installation is completed, do not plug the power adapter into the wall outlet or turn on the scanner until you are instructed to do so.
- Always hold the scanner firmly when lifting from the stand or when using the scanner. Do not shake the scanner.
- Always return the scanner to the cradle when it is not in use. Do not place the scanner in heated or wet surfaces as this can cause damage to the scanner.
- It is normal that the scanner gets warm when in use. Do not block the ventilation holes on the bottom of the scanner. If the scanner overheats, the scanner will stop working.
- Ensure that you use only the supplied power adapter, power cable, and USB cable.


Attach and Detach the scanner tip

Caution


- Do not place your finger(s) on the mirror of the tip when attaching as this may result in damage to the mirror.
- Do not place your finger (s) on the lens of the tip when detaching the tip to prevent damages to the lens.
- Do not attempt to clean the outer units and inner optical components on the front end of the scanner with any sharp objects or other such tools, which may result in scratches and damage to the scanner.

Follow the step below to attach the scanner tip:

Hold the scanner tip firmly with your thumb and forefinger on both sides, and then gently attach the tip facing downward to the scanner.



Follow the step below to detach the scanner tip:




Connect and disconnect the scanner

Follow the steps below to complete the connection and disconnection:


Connection

Follow the steps to complete wired connection.

1. Push the scanner tip hard to the scanner main body to ensure firm attachment.



2. Connect the scanner cable to the Type-C port.



3. Click the shortcut icon of IntraoralScan (Named "Dental Launcher") on the desktop to launch the software.

Disconnection

Caution

Do not roll the cable over the handle of the scanner or even create any sharp bends in the cable after you disconnect the scanner.

1. Quit the IntraoralScan scanning software.
2. Disconnect the scanner with the computer.

Specifications

Parameter	Description
Type name	Intraoral Scanner
Model name	Aoralscan Elite
Device lifecycle	8 years

Scanner	
Scanner tip types	Big tip, standard tip and mini tip
Scan field (L×W)	<ul style="list-style-type: none"> • Big tip: 19 mm × 14 mm • Standard tip: 16 mm × 12 mm • Mini tip: 12 mm × 9 mm
Scanning tip dimensions (L×W×H) (Tolerance: ±2 mm)	<ul style="list-style-type: none"> • Big tip: 95 mm × 30 mm × 26 mm • Standard tip: 93 mm × 30 mm × 26 mm • Mini tip: 92 mm × 30 mm × 26 mm
Scanning depth	22 mm (-2 mm ~ 20 mm away from the window of the scanner tip)
Scanner size (L×W×H) (Tolerance: ± 4 mm)	245 mm × 30 mm × 26 mm
Scanner weight (Tolerance: ±20 g)	124 g (with a standard tip and without cable)
Scanner tip maintenance	Sterilized and disinfected by users (Maximum: 100 times)
Scanner tip connection	Pluggable
Light source	LED and laser
Wave length	Blue laser: 450 nm White LED: 400 nm - 780 nm
Output	STL、OBJ、PLY
Power	Input: 5V DC/3 A

Cradle	
Cradle size (L×W×H) (Tolerance: ±3 mm)	103 mm × 80 mm × 67 mm
Cradle weight (Tolerance: ±20 g)	242 g

Calibrator	
Calibrator size (L×W×H) (Tolerance: ±3 mm)	236 mm × 50 mm × 50 mm
Calibrator weight (Tolerance: ±30 g)	431 g
Calibrator cable length	1.5 m


Environmental Requirements

Operating and storage requirements	Description
Operating temperature	10°C – 40°C
Scanner storage/transport temperature	-30°C – 60°C
Working relative humidity	30%RH – 80%RH
Scanner storage/transport relative humidity	30%RH – 90%RH
Air pressure	70 kPa – 106 kPa

About IntraoralScan

System Requirements

Before installing and running the supplied software programs, your computer shall meet the following requirements:

 **Note**

Your PC shall meet the safety requirements of IEC 60950-1, IEC 62368-1/GB 4943.1.

Computer

Item	Classification	Description
Intel CPU	/	Intel® Core™ i7-8700F or above Recommend: Intel®Core™ i7-11800H
Memory	RAM	Minimum: 16GB Recommend: 32GB
Hard Disk Drive	SSD	Minimum: 256GB Recommend: 2TB
Graphic Card (GPU)	Discrete graphic card	NVIDIA RTX 2060/2070/2080/3060/3070/3080/4060/4060Ti/4070/4080/4090
Operating System	Windows	Windows 10/11 Professional (64-bit)
Ports	Input/Output ports	At least 1 Type-C port and 1 Type-A USB 3.0 (or higher) port

Note

- When running the software for the first time, a prompt of restarting the software to perform the PC performance test will be popped up. After test, the software will automatically configure the scanning speed to suit your computer's performance. Click Yes to restart the software.
- Your PC shall meet the safety requirements.
- The default virtual memory is 4.8GB, and insufficient virtual memory may cause software crashes.

Register and Log in

Open the software and enter the login interface.

Register


If there are no SHINING 3D accounts, click **New User? Click here to register** to enter the registration interface.

Steps

1. Select a role type and fill in registration information such as contact information and phone number.
2. Read **Privacy Policy** and select the check box.

3. Click **Register**.

4. Click **Back** to return to the login interface. Enter the account and password to log in.




The registration form is displayed on a light blue background. At the top left, there is a back arrow and the text 'Back'. The form includes a country selection dropdown set to 'China 中国' with a globe icon and a downward arrow. Below this is the text 'After registered, this cannot be changed'. The form contains several input fields: 'Email or phone', 'Contact', a 'verification code' field with a 'Send' button to its right, 'Password', and 'Confirm Password'. A 'Click for Smart Verification' button with a green circular icon is positioned above the verification code field. A note below the password field states: 'The password must be between 8 and 32 digits long and contain at least a combination of letters, numbers and special characters'. At the bottom left, there is a checkbox for 'Privacy policy'. A large blue 'Register' button is centered at the bottom of the form.

Note

- You can join in the temporary institution directly for 3 days without approval.
- If you found an institution or are approved to join in a common institution in 3 days, then you will be moved out from the temporary institution.
- One user can only join in the temporary institution once.
- After registration, messages about dental cloud will be sent to the phone or e-mail you entered.

Log in

- The users can log in with a dental cloud account or register a new account.
- The users can log in by **password** or by **verification code**.
- After logging in, registration information and connection can be checked.



Note

The account information can be deleted in the drop-down menu of the account login box.

Activation

Activation for the first time

When using the device for the first time, it needs to be activated.

1. [Connect the device following the instructions.](#)
2. After successfully connecting the device, launch the software while being online.
3. [Register or log in to your account.](#)
4. Once you are in the main interface of the software, a prompt for device activation will appear. Click **Yes**.
5. Fill in the company, industry application and organization type (Dental hospital, Stomatology department of general hospital, Dental clinic, Lab, Reseller and Others). Then, click **Submit**.
 - Industry application: Including Prosthodontics, Orthodontics, Oral implantology, General dentistry, Pediatric dentistry, Hearing-aid and Other. Multiple selection is supported.
 - Organization type: Including Dental hospital, Stomatology department of general hospital, Dental clinic, Lab, Reseller and Others.


The screenshot shows a light blue dialog box titled "Activate" with a close button (X) in the top right corner. Inside the dialog, there are three required fields, each marked with a red asterisk:

- *Company:** A text input field with the placeholder text "Please fill in the real name".
- *Industry application:** A dropdown menu with the placeholder text "Please select your industry application".
- *Organization type:** A dropdown menu with the placeholder text "Please select your organization type".

At the bottom of the dialog, there are two buttons: a blue "Submit" button with a checkmark icon and a grey "Cancel" button with an X icon.


6. Click **Activate** to activate the device. After successful activation, you can start using the device normally.

Authorization renewal


Click  **Settings** → **About** to check the authorization remaining days.

When the authorization of the device expires, the users can click **Request an Authorization Renewal**. On the window of Authorization Tools, select authorized days, and click **Apply**. After the application is approved, the authorization will be renewed.

Upgrade Software

Click " > **About**" on the upper right corner of the homepage.

- Click **Check for Updates** and it will automatically check whether there is a new version. If there is a updated installation package on the network, you can choose to download and install it now or next time.
- If you have chosen **Check updates automatically**, the software will download the updated package when connected to the network. And you can also choose to install it now or next time on the pop-up box.


 **Note**

Please make sure the computer installed with the software has been connected to the network before upgrading.

Calibrate the Scanner

Under these circumstances, we recommend that you shall execute the calibration for the scanner to ensure the accuracy of scanned data:


- The initial setup of the scanner is completed.
- The scanner has been used for a period of time (e.g. 2 weeks).
- Scanner brightness adjustment is recommended once every 3 months.
- The scanner is accidentally dropped or strongly vibrated.


 **Caution**

- Before calibration, please connect the calibrator with the cable.
- Please keep the calibrator in bag to prevent from dust when it is not used.
- Do not drop or strike the calibrator.
- Do not splash water into the calibrator.

Follow the steps below to perform the calibration:


1. Gently take the tip off from the scanner.






 **Caution**

- Do not place your finger(s) on the mirror of the tip when detaching as this may result in damage to the mirror.
- Store the detached tip in a safe place, e.g. a dental instrument tray, for future use.

2. Gently take the calibrator onto the front end of the scanner.



3. Click  on the main interface to enter the calibration interface.
4. Ensure the scanner is plugged into the calibrator firmly. Click **Start** to calibrate.
Normally the calibration takes approximately 7 minutes.
5. After the calibration is finished, click  in right corner to exit.
6. Gently take the Calibrator off the scanner.
7. Reattach the scanner tip to the scanner for later use or put the dust cap onto the scanner to prevent damage and dust.

 **Caution**


Make sure that the Calibrator is removed from the scanner after the calibration is done. Otherwise, the Calibrator temperature may get very high.

Create Order

Interface

After logging in, you can check your orders.

Preview



① Home and dental cloud connection



: Click to return to the homepage.

Dental cloud connection

- : Normal connection.
- : Failed connection.
- : Unstable connection.

② Settings

Check Users information and connection, Calibration, Record and screenshot, Endoscope, Order settings and Help.

Users information and connection

Click to check user profile and the connections.


When adding new connections, click **New** to search a target organization and click  to send the application.

Caution

If you fail to add a new connection and a tip informing that necessary information is not completed, please log in the dental cloud to complete the relevant information, then start to re-establish the relationship.




Screenshot

Record the screen, take a screenshot and check screenshots.

- Click **Screen Record** to record the screen and click  to end it.
- Click **Screenshot** to capture the picture and you can edit the picture by clicking it displayed in the lower right corner.
- Click **Find the screenshots/screen records** to view the storage path.

Endoscope

With Endoscope, users can check the teeth in advance and create the corresponding orders.

Click  to start Endoscope.  to take screenshots and click  to check the path where the screenshots are stored.

Caution

Please connect or wake up the scanner when using Endoscope.

Help



Click **Help** for the User Manual, Shortcut Instructions, User Guide (Enabled by default), Support Center, Remote Assistance, Feedback, Official Account Customer Service and Information.

- Support Center: Check our technical support.
- Feedback: Submit your issues, suggestions or requirements or others.

Please fill in your contact information and the device serial number, and describe what problems you encountered. You can upload attachments for better description.

- Official Account Customer Service: Scan the QR code to follow our official account.
- Information: Check new functions and other information about the software.

③ Create/ Import Orders

-  : Click to create a new order, as well as a scanning order. Details in [Create/Import Order](#).
-  : Click to import an order from ExoCAD and IntraoralScan.

④ Advertisement











See information of other productions.

Click **Dental Cloud** to enter cloud platform.

⑤ Orders


Including [Order list](#) and [Patient list](#).


Order list

Icon	Description
 Search	Enter order number or patient's name for searching.
 Filter	<ul style="list-style-type: none"> Click to choose filtering conditions (e.g. order status, sources, time) to find the targeted order. If check Hide cloud order, the order stored in the cloud platform (No such order data locally) is hidden and not displayed. Click  again to close it.
 Multi-selection mode	Process the duplication or the deletion of multiple orders. Click  on the right top to exit.
 Sync information	Synchronize the cloud order to the software.
 Card Mode & Table Mode	Switch arraying methods between card & table. Check order information (e.g. patient's information, order status and creating time). Open a certain order by double-clicking it.
 Send	Send order information and scanning data to other labs.
 Duplicate	Create a new order identical to the current one and re-scan.
 Delete	Delete the selected order(s) from the list. Delete the file and scanning data permanently by checking Delete File and Scanning Data.

Order details

Click one order to check details, including its order information, the tooth map and the scanned model. Users can change send settings as well.



- Export: Click  to save the order to the local path.

The users can change the folder name, export path, CAD type and the data format. After changing, click **Confirm** to check the local path.



Export Option ✕

Folder Name

Export Path
 📁



CAD Type
 ▼

Export Data Format
 OBJ PLY STL


- Edit: Click  to edit order information again, including changing information of the patient and the doctor, and editing notes.
- Pre-design: Click  to enter the interface of pre-design.
- Convert: When the dentistry type is **Restoration** or **Other**, converting the order is supported. The users can enable pre-op scanning when converting the order.

 **Note**







If **Pre-op** is already enabled, converting the order is not supported.


- Duplicate: Create a new order identical to the current one and re-scan.
- Delete: Delete the selected order(s) from the list. Delete the file and scanning data permanently by checking Delete File and Scanning Data.
- Exocad: Enter Exocad to design the model.
- Explore: Click  to explore the current local path where the order is saved.
- Send: Click  to send the order. More details can be found in [Send Order](#).

Patient list




Patient Name	Patient Sex	Date of birth	Contact	Create Time	Operation
[Redacted]	Male	2007/6		6/8/23 1:56 PM	
[Redacted]	Male	2007/6		6/8/23 9:40 AM	
[Redacted]	Male	2016/6		6/8/23 9:40 AM	
[Redacted]	Male	2017/6		6/8/23 9:40 AM	
[Redacted]	Male	2017/6		6/8/23 9:40 AM	
[Redacted]	Male	2007/6		6/8/23 9:30 AM	
[Redacted]	Male	1991/1		6/5/23 8:09 PM	
[Redacted]	Male	2002/6		6/5/23 8:03 PM	
[Redacted]	Male	2002/6		6/5/23 8:02 PM	
[Redacted]	Male	1998/6		6/5/23 3:39 PM	
[Redacted]	Male	2008/6		6/5/23 3:37 PM	
[Redacted]	Male	2002/6		6/5/23 1:51 PM	

Icon	Description	Icon	Description
 Search	Enter patient's name to search orders.	 Sync information	Synchronize the cloud order to the software.
 Add a patient	Click the button to enter the New Patient window. Enter name, gender, age and other necessary information to add a new patient.	 Create an order	Create a new order based on the certain patient. Relative information will be filled automatically and can be edited. More details can be found in Create/Import Order .
 Edit patient	Click to edit relative information of the certain patient.	 Merge patient	Merge the same patient for better patient management. More details can be found in Patient Merge .

 **Note**

- If the user edits patient information, creates new orders or adds a patient when the computer is not connected to the network, the relative information will be synchronized to the cloud.
- Only institution manager can merge patients.






Order Settings

Click  to set order parameters. To restore parameters on **Order Settings** interface to the default.

General

Parameter	Description
Language	Users can set the software interface language: Chinese (Simplified), English and other languages. The default language is the language selected during software installation.
Default Dentistry Type	Select a default dentistry type from Restoration, Implant, Orthodontics, Removable Denture and Other .
Dental Notation	FDI World Dental Federation notation and Universal numbering system . The default is FDI World Dental Federation notation.
Order Save Path	Set the order saving path. The default is to save the DentalOrder file in the installation package disk. If the installation package is stored in the C disk, the software installation path is C:\DentalOrder.
Exocad DentalCAD Path	To combine EXO software to IntraoralScan, enter EXO software location.

General

Language	English 
Default Dentistry Type	Restoration 
Dental Notation	FDI World Dental Federation notation 
Order Save Path	D: /DentalOrder 
Exocad DentalCAD Path	../ ../DentalCADApp/bin/DentalCADApp.exe 

Naming Rules

Default naming rule is "Create date_Doctor Code_Order ID". Still, you can define your own naming rule. Select wanted naming rules (order ID is a must) one after another and the final naming will be displayed in earlier selected rule > later selected rule sequence.

Note

When the order ID is checked, you can set that the ID is displayed out of order. After enabling it, the order ID number will be displayed randomly.

Order File Naming Rule

<input checked="" type="checkbox"/> Order ID	<input type="checkbox"/> Order Number
<input checked="" type="checkbox"/> Doctor Code	<input type="checkbox"/> Doctor Name
<input type="checkbox"/> Patient Code	<input type="checkbox"/> Patient Name
<input type="checkbox"/> Operator Code	<input type="checkbox"/> Operator Name
<input checked="" type="checkbox"/> Create date	

Naming Rule: "Order ID_Create date_Doctor Code"
Example: "001_2022-09-23_001"

Saved Format

Parameter	Description
Save Exocad order	The exocad order *.dentalProject is saved while saving the order. You can import the order into the EXO design software for design. The default is on.
Save 3 rd party order	Save third-party software order when an order is finished. A third-party software order is saved while saving the order. The folder named "original order name_thirdParty" is generated under the order path, and the <i>.xml file and the scanning result data .stl/*.ply</i> (the coordinate system is consistent with the third-party design software) are stored. You can create an order into a third-party design software and import the data to design. The default is on.

Save Format

- Save Exocad order
- Save 3rd party order

Data Cleaning

By setting the storage time, project files beyond the time can be deleted automatically or manually, and more space is reserved for storing newly scanned project files.

Automatic cleaning:

1. Check **Automatic cleaning of historical project files**.
2. Set interval time to remind and retain days within the order.

Note

After logging in, you will be prompted whether to delete historical data according to the date you set.

Manual cleaning:

1. Click **immediately clean up**.
2. Set the end date.
3. Click **OK**.

Note

- Only temporary data will be deleted.
- Data before a certain date can be cleared manually, while the data on that date can be retained.
- When the disk space for saving project files is lower than 20G, a prompt will pop up to show that clear the data immediately.
- You can choose to save data for 30 days, 90 days or 120 days.

Data Cleaning  (The software will clear the historical project files after the cleanup setting is turned on)

 Clear

Automatically clean up historical project files

interval day reminder to clean up

Keep orders within days

Create Order

Create different kinds of orders according to the scanning results. Orders will show the types of patients, tooth place and scanning process.

Order Information

ID: 031 Create Time: 6/6/24 4:33 PM

Order Number: Type:

Patient Name*:

Doctor:


Operator:

Lab:

Dentistry Type: Restoration Implant Orthodontics
 Removable dentures Other

Please enter the remarks:

Tooth Selection



Full Crown Pontics Inlay Veneer Antagonist

Implant-Based:

Material:


Scan a pre-op model:

Tooth Shade:

Order Information

Click  to enter order information.

1. Choose the types of patients. **FirstVisit** or **FollowUp**.

2. Click  to enter required information of the patient such as name and age.

 **Note**

If it is a followup, you can enter the patient's name directly.

3. Enter doctor's name to search the doctor. Or click  to add a new doctor: enter the doctor's name (and email or phone if needed) in the **Doctor Configuration** interface and click  to add the new doctor.

4. Select a dentistry type:

- **Restoration:** Replace or restore your missing parts of your tooth structure. The treatment types include **Full Crown, Pontic, Inlay, Veneer** and **Antagonist**, and the dentist will choose the proper types according to the situation.
- **Implant:** Replace damaged or missing teeth with artificial teeth. Implant treatment types include **Full Crown, Bridge, Upper, Lower** and **Full**, and the dentist will choose the proper types according to the situation.
- **Orthodontics:** Deal with irregularities of the teeth and their correction. Treatment types include **Fixed Ortho, Mobile Ortho** and **Invisible Ortho**.
- **Removable dentures:** Check the bite relationships of edentulous patients. The types of removable dentures orders can be **Full denture, Partial removable denture** or **Antagonist**. The **Get bite record** can be natural teeth bite, bite rim or old denture.
- **Other:** Receive other dental examinations and treatments including **Oral Exam, Caries Filling, Oral RCT** and **Period DontTm**.

Tooth Selection

1. Based on the teeth to be treated, select the teeth directly on the teeth map.


Keyboard shortcuts	Description
Left-click	Select one tooth or teeth.
Right-click	Cancel settings of a tooth.
Ctrl+ Left-click	Copy the last defined dental restoration type to the currently selected dental position.
Shift + Left-click	Copy the restoration type defined by the previous tooth position to all teeth between the selected tooth position and the previous tooth position.

2. Select implant type.

Implant Type	Description
Custom Abutment	Add the step of scanning ScanBody .
Custom Abutment (manual positioning)	Add the step of scanning Abutment .
Screw Retained	Add the step of scanning ScanBody for dental implant.
Screw Retained (manual positioning)	Add the step of die scanning.

3. Set **Material**.

4. Set to open **Scan a pre-op model**.

 **Note**

When selecting **Orthodontics**, the dentist should choose the stage.

5. Click the buttons for further operations.

Operation	Description
Scan	Click to enter scanning interface.
Save	Click to save the current order for subsequent scanning.
Explore	Check the saved orders on the pop-up window.

Import Orders

Click **Import order** and select INPROJECT or DENTALPROJECT files. Supports importing orders saved by IntraoralScan software and ExoCAD orders.

When importing scanned orders and clicking **Go to scan**.

To load the last scanned data, click **Yes**. To rescan, click **No**.

Patient Merge

In the **Patient List**, merging patients is supported when these patients are the same. Merging patients make it easy to manage patients.

When importing orders, the imported patient can be merged into an existing patient.





Note

Only manager of the institution can merge patients.

Patient Merge in Software

Steps

1. Click  to enter the window of merging patients.
2. The window of **Patient Merge** shows the basic information of the primary patient and similar patients.
3. Click the similar patients and the curtain patient will be added in the **Selected** on the left. Move the cursor on a certain patient and click  to delete the selected patient.

Patient Merge ✕

Primary patient

👤 123 🏠 2007.05 ♂ Male 📞

Similar patient Selected(4)

👤 123 📞 🏠 2007.05 ♂ Male	👤 123 📞 🏠 2007.05 ♂ Male
👤 123 📞 🏠 2007.05 ♂ Male	👤 123 📞 🏠 2007.05 ♂ Male
👤 123 📞 🏠 2007.05 ♂ Male	👤 123 📞 🏠 2007.05 ♂ Male
👤 123	👤 123

The selected patient files are retained, and after merging, merged pat...

✓ Accept
✕ Cancel

Note

When there are no similar patients, manually enter the patient's name or phone number to search.


4. Click **Accept** and a prompt is displayed, indicating that **The combination of patients cannot be withdrawn**. After confirming, click **Accept** again to merge patients.

Tip ✕




Patients are irrevocable after combining, confirm?

✓ Accept
✕ Cancel

5. A prompt showing "**Merge Success**" indicates that the selected patient(s) has been merged with the primary patient.



Functions

Icon	Description
 Select all	Click the button to select all similar patients.
 Clear all	Click the button on the right of Selected to clear all selected patients.
 Search	Enter the patient's name or phone number to manually search patients


Caution


- The combination cannot be withdrawn. And the combination record will be uploaded to the cloud.
- After merging, the merged patients won't be displayed in the similar patients any more.
- Merging patients is accessible only when connected to the Internet.

Merge Patients of Imported Order


When importing orders, imported patients can be merged into existing patients.

Steps


1. Click  to import orders.
2. When the imported order isn't from the current institution, or its patient information is empty, a prompt of "Merge this order to an existing patient?" is popped up.



3. Click **Yes** to enter the window of **Order attribution modification**.
4. The window of Order attribution modification shows the original patient of the imported order and attribution patient.
Enter the patient's name or phone number to search.




5. Click **Accept** and the imported order will be attributed to the selected patient.

 **Caution**


When the imported order comes from the third-party or from Exocad, and information of the imported patient is totally same with that of an existing patient, then the order will be directly merged into the existing patient.

Scan

Scanning Preparations

 **Caution**

- Concerning hand hygiene and personnel safety when performing a scan, you must wear clean surgical gloves through the whole process.
- If you need to use isolation film or protective film (such as blue film) to wrap the scanner, please do not block the air inlet and outlet of the scanner to avoid affecting the heat dissipation.



Intraoral Environment

- Make sure there is no foreign body or blood in the mouth after gargling. Stop the bleeding if necessary.
- If necessary, ask the patient to keep the tongue still and move it to the other side of the mouth.
- Consider using a dental three-way syringe or cotton to dry the tooth surface before starting the scan.
- Turn off the oral light on the dental chair and start scanning.
- Consider using aspirators or cotton to keep the surfaces dry during scanning.
- If necessary, consider using an oral mirror to help create space while working in the narrow area between the teeth.

Scanner Preparation

- Ensure that the scanner tip, scanner body, and cradle are properly pre-cleaned, disinfected, or sterilized. More Details see [Pre- cleaning, disinfection, and sterilization](#).
- Ensure that the scanner tip has no scratches or is not damaged. Additionally, the tip is firmly attached to the front end of the scanner body.

- Ensure that the scanner connection is ready; it is correctly connected to a power source and powered on, and IntraoralScan is launched and ready to work.
- To avoid condensation on the mirror of the tip when scanning, the scanner tip must have been warmed up.
- Calibrate the scanner and verify the accuracy of the acquisition regularly. More details see [Calibration](#).

Scanning Position and Path

- Avoid direct light from any light source, e.g. dentist chair lamp, to shine on the area you are working on.
- Hold the scanner steady by resting it on the tooth surface and keep the scan tip window in the range of -1 mm to 16 mm from the teeth.
- When scanning, slowly move the scanner and simultaneously check the scan results on the screen to ensure that the scanning is of good quality.
- When scanning, the scanner tip should be centered over the teeth, and each movement should align with the cross-hairs, following the lower and upper dental arch shapes.
- A complete scan data of a single area includes the surfaces of bite, lingual, buccal, interproximal contacts of the adjacent teeth, and 2-3 mm buccal gingiva.
- A complete scan data of a single case includes the lower jaw, upper jaw, and bite registration.
- When scanning, change the scanning angle to 35-55 degrees to create overlaps. It is important to achieve an overlap of at least 30% between each acquisition. If the overlap is small, the alignment may fail.
- To scan the occlusal surface of the teeth, hold the scanner at a 90-degree angle; to scan the buccal and lingual surfaces of the teeth, hold the scanner at a 45-degree angle.
- Inspect the scanned image in the 3D scan view window (IntraoralScan) and pay attention to warning messages.


Heat the Scanner Tip

To ensure the quality of optimal images, you should prevent condensation on the scanner mirror before each scan by heating the scanner tip.

Follow the steps below to warm up the scanner before starting an acquisition:

1. Ensure that the scanner tip, scanner body, and cradle are clean. To clean them, please check [Pre-cleaning, Disinfection and Sterilization](#)
2. Gently and carefully attach the scanner tip to the scanner body with the mirror facing downward.
3. Connect the power supply to the scanner. More details see [Connection](#).
4. Place the scanner in the cradle to secure it in place.
5. When the LED ring light on the scanner body lights up green, the heater automatically turns on and detects the temperature.

6. If the temperature of the scanner tip is lower than the set value for anti-fogging, a notification message of pre-heating and current temperature appears.
7. When the message disappears, the scanner tip has been heated. The scanner is now ready for scanning.

 **Note**



- The heater helps keep the scanner tip temperature in a normal range.
- The scanner tip is being heated whenever power is supplied, even if the scanner is in standby or sleep mode.
- If the heater does not reach the necessary temperature for preventing condensation during scanning, a tip will appear.

Scanning Settings

Click  and then click the **Scan Settings** tab to view or set scanning parameters.

After changing, click  to reset the settings.

General

General function	Description
IMU menu	Enabled by default.
Auto-brightness	Adjust the camera brightness to capture clear images. Enabled by default.
Refined scan	Checked by default. Suitable for scanning refined areas, making processed data more precise and improving data details. In the scenario of teeth restoration and implant, the refined areas will be automatically recognized after enabling this function. While in orthodontics and other scanning processes, the software will automatically start refined scanning, but will not automatically recognize refined areas.
Supports alternate day scanning	Enabled by default. Scanned orders support unlimited add-scan.
Filled holes will be displayed in grey	<p>Check to make the holes of teeth displayed in grey, as shown in the following picture. This function is unchecked by default, making the holes of teeth will be displayed in color.</p> 
Auto bite optimization	<p>Optimize the data automatically when scanning the occlusion.</p> <p>When scanning bites,  indicates that auto bite optimization is checked. Click to uncheck it.</p>
Pick up and scan	Enabled by default. When it's enabled, the scanner automatically starts scanning once it's lifted from the cradle. Otherwise, please click the start button on the software, or press space button on the keyboard/start button on the scanner to start scanning.

Preferences Settings

- Scanning order: Support to set the scanning order according to the operation habits. The default setting is working jaw priority.

- Scanning View: Support to set various scanner views.

View	Description
Scanner view 1	Camera Window and scanning data clockwise rotation 90°.
Scanner view 2	The camera window displays the real scanning scenarios same as the scanned data displayed angle.
Operator perspective	<p>The scanning position of the operator perspective is divided into Patients Front and Patients Rear.</p> <p>Patients Front: The scanner perspective of lower jaw scan does not change, and the upper jaw data are mirror images.</p> <p>Patients Rear: The scanner perspective of upper jaw scan does not change, and the lower jaw data are mirror images.</p> <p>When selecting the operator perspective, you can check Mirror Show Upper Jaw, This parameter is mainly designed for the sitting position of the patient. After checking, the camera window is mirrored up and down, the scanning data is clockwise rotated clockwise by 180°, and the mirror image shows the upper jaw.</p>

Camera

Camera top, bottom, left, and right margins' values are acquired from the camera by default; When selecting **Operator Perspective** as scanning view, you can also manually set the cropping parameters for the top, bottom, left, and right margins of the camera window.

Adjustable left or right margins in the range of [8, 300].

Adjustable top or bottom margins in the range of [8, 400].

 : Restores the margin values in the camera.

Music Settings


Click  , upload the music documents.

PC Performance Test

- Start testing PC performance or restore factory settings.
- After the test, the software will automatically configure the scanning speed to suit your computer's performance.

Scan Interface Introduction





Function preview



① Camera Window

Displays the real scanning scenarios and supports to set camera parameters.

Move the mouse within the image window to adjust the camera parameters.


Icon	Description	Icon	Description
	Click to enlarge the image window, and click again to restore it to the default size.		Take a screenshot of the current image in the image window, and save the image in the "Endoscopelmg" folder under the current order directory.
	Scanning indication switch. Disabled by default. When turned on, the camera window displays the current unscanned data area in purple when scanning.		After taking screenshots, click it to open the path of saved screenshot.

② Guide

Contains the current process, operating diagrams and detailed tips.

③ IMU menu/Intraoral scan/Edentulous scan

- **IMU menu:** Double-press the button on the scanner to open the IMU menu. Move the scanner tip and stop on the options for several seconds to go to the next step/former step/view the model/start scanning.

IMU menu is enabled by default. To disable this function, please click  → [Scan Settings](#) → [General](#) → [IMU menu](#).

- **Intraoral scan:** Checked by default, applicable to intraoral scanning. It is recommended to leave this checkbox unchecked when scanning dental models.
- **Edentulous scan:** It can be checked after checking Intraoral scan. If it is checked, scanning process for edentulous patients will be speeded up.

④ Progress Bar

Displays the current position within the whole progress.




To view the other processes, click **⋮**.

To add extra processes, click **+**.

The progress bar includes upper jaw, lower jaw, bite, additional scan and so on.

Additional scan

Additional scanning is used for scanning groups of additional models for extra information.

-  : Click to add additional scan.
-  : Click to rename the additional scanned model.
-  : Click to delete the additional scan.

Note






- When adding another group of additional scan, the user should finish the former one.
- The disk memory must be at least 25 GB when adding additional scan.

⑤ Preview

Supports previewing the scanned data. To gain a comprehensive view to the model, details can be found in [Operation Skill](#).

















⑥ Operations

Including scan operations.

Icon	Description	Icon	Description
 Manual align	Manually align the extra-matched data with the main data.	 Impression scan	Combine the impression scan with the main data. During this scanning process, user can edit the data, remove the isolated data and so on.
 Start/Pause	Click to start or pause the scanning.	 Reset	Reset current scanning data.
 Finish	Finish scanning and save the current scanning data.		

⑦ Scan frames and time

Show the number of scanned frames, scan and data processing time and the scanner status.

Icon	Description	Icon	Description
 	Connected.	 	Not connected.
 	Connecting.	 	Wrong connection.
 	Without scanner tip.	 	Overheated.
 	In standby status.	 	In scanning status.

⑧ Extra match

Extra-match can align the scanned data with third-party accessories for checking multiple sets of data.

More details can be found in [Extra match](#).

⑨ Tools and additional functions


Tools and additional functions can be used for better scanning and previewing the model. More details can be found in [Scanning tools](#) and [Additional functions](#).


Functions


Click  > **Shortcut Instructions.**

Shortcut Instructions ✕


Model Operation



 Rotate Model



 Scale Model



 Translate Model

Tooth Selection



 Select Teeth


shift + 
 Shift Multi-Select



 Drag multiple selections

ctrl + 
 Copy the information of th...

Scan

Motion Control 

Next Space >3S or 

II Start/pause scanning Space or 

D Data Indicator T

I Intelligent Scan A


Turn on instruction

Pre-design

T Texture C

S Screenshot S

Inspection Report

Endoscope Snap 

Note

After checking **Turn on instruction**, a prompt will pop up when you click a function button which is supported by keyboard shortcuts for the first time.

Scan Model

Action	Instruction
Space	Start scan or move forward.

View Model

Action	Instruction
Hold the Left or Right Button and Move	Rotate the model
Hold the Left and the Right Button and Drag	Pan the model
Scroll the Mouse Wheel	Zoom in or out

Edit Model

Action	Instruction
Shift + Left Button	Switch the tool to the eraser
Shift + Scroll the Mouse Wheel	Resize the brush
Del	Delete the data in red color

Scan Workflow



Pre-op Scanning

Pre-op scanning means to collect and retain intraoral data before treatment.

Create an Order

Pre-op scanning takes restoration of full crown as an example.

Enter necessary information.

1. Choose patient types, **FirstVisit** or **FollowUp**.
2. Enter necessary information of the patient such as **name** and **age**.
3. Add a doctor/operator. Click  and enter the doctor's name, then click  to confirm.
4. Choose **restoration** or **implant** as the dentistry type. Add notes such as teeth shade.

Dentistry type	Description
Restoration	Restore missing parts of your tooth structure. Select the restoration type according to the reality such as full crown, pontic, inlay and veneer.
Implant	Replace a lost tooth with a dental implant. Select the implant type according to the reality such as full crown, bridge, upper jaw, lower jaw and full jaw.

5. Select the tooth.

Click on the number of one tooth to select it (click right mouse button on the number to cancel the selection). Multiple selection is available.

- Click **Ctrl+left mouse button** to copy the restoration type of last tooth to the currently selected tooth.
- Click **Shift+left mouse button** to copy the restoration type of last tooth to all the teeth between the last tooth and the currently selected tooth.

6. Select **Yes** in the pop-up window of Scan a pre-op model.

7. Click **Scan** to enter the scanning interface.

The screenshot shows a software interface with two main sections: 'Order Information' and 'Tooth Selection'.

Order Information:

- ID: 009, Create time: 6/28/21 10:48 AM
- Order Number: 009
- Patient Code: 009, Patient Name: Jane
- Doctor: 001, Doctor001
- Technician: 001, Technician001
- Dentistry Type: Restoration (selected), Orthodontics
- Notes: Shade: None, Additional comment: None.

Tooth Selection:

- A diagram of the upper and lower dental arches is shown. Teeth are numbered 1-32. Teeth 27 and 35 are highlighted in blue.
- A 'Clear All' button is located below the arch diagram.
- Restoration type buttons: Full Crown, Pontic, Inlay, Veneer, Antagonist.
- Implant-Based: No implant
- Material: Acrylic/PMMA
- Scan a pre-op model: Yes (highlighted with a red box)
- Buttons at the bottom: Cancel, Explore, Save, Scan.
- Legend: Full Crown (blue), Pontic (green), Inlay (purple), Veneer (orange), Antagonist (pink).

Start Scanning

Steps of pre-op scanning: pre-op of the upper jaw (the treated tooth is on the upper jaw) > pre-op of the lower jaw (the treated tooth is on the lower jaw) > full jaw > upper jaw > lower jaw.




Note

- Default scanning steps are: pre-op of the upper jaw > pre-op of the lower jaw > full jaw (scanning the bite) > upper jaw > lower jaw > full jaw (aligning automatically)
- The scanning steps of upper jaw/ lower jaw (post-op) is similar to implant jaw. After loading pre-op data, the doctor needs to dig a hole on the tooth, then starts scanning post-op teeth.

Scanning steps:


1. Scan pre-op upper jaw with a disinfected scanner.



Caution

Please keep moving the scanner at a relatively same speed.
The distance between the scanner tip and the surface of the teeth should be 3 mm to 5 mm.

2. Dig a hole on the treated tooth to avoid alignment failure resulting from the inconsistent scanning data of pre-op and post-op.
 - (1) Move the cursor to the treated tooth.
 - (2) Scroll the wheel up and down to zoom in and out the tooth model (make the treated tooth as big as the circle). Select the tooth place and make preparation for digging a hole.
 - (3) Double-click to dig a hole.




3. Scan the post-op teeth.

4. Repeat step 1 to 3 if teeth on the lower jaw need to be treated.


After scanning, it will automatically start editing **Refined Area**. The recognized refined areas are colored blue. Manually selecting refined areas in key parts is supported. This function can provide more partial details in pre-design, and is often used as restoration of sharp edges of inlay and veneer.

5. Scanning the bite part.


Scan the left bite part and the right bite part. Click  or press space to start scanning. The software will start alignment automatically when part data is collected.


When the upper jaw and lower jaw is aligned, click  or press space to pause scanning. Check the bite.

[Multiple bites](#) scanning is supported.

 **Note**

- There is no need to scan full jaw for collecting bite data. You can scan left and right (recommended) or left, right and front teeth to align the models of upper jaw and lower jaw.
- The maximum scanning frames of each bite group is 300. The software will stop the scanning once the frames are up to 300.


6. Click  on the right panel to edit the model. Details see [Data Edit](#).

7. Click  to confirm the changes and exit the edit interface.

Alignment

Click  **Static Alignment** or  **Manual Alignment** to improve the performance of alignment.


Static Alignment

Click  to automatically align the bite once again.



Manual Alignment

Click  to enter the interface of manual alignment.


1. Please select 3 sets of corresponding points in occlusal and upper jaw as well as occlusal and lower jaw respectively.
2. Preview the splicing effect after the corresponding points have been selected.



Caution

If dynamic bite is scanned, when click  to start static alignment or click  to finish manual alignment, the data of dynamic bite will be deleted and a tip will pop up. Please confirm whether to delete the data according to the tip.

Scanning Completed



Click  to complete scanning. Once the data is post-processed, the software will enter the interface of pre-design.

Restoration Scan

Create an Order

Pre-op scanning takes restoration of full crown as an example.
Enter necessary information.

1. Choose patient types, **FirstVisit** or **FollowUp**.
2. Enter necessary information of the patient such as **name** and **age**.

3. Add a doctor/operator. Click  and enter the doctor's name, then click  to confirm.


4. Choose **restoration** as the dentistry type. Add notes such as teeth shade.

5. Select the tooth.

Click on the number of one tooth to select it (click right mouse button on the number to cancel the selection). Multiple selection is available.

- Click **Ctrl+left mouse button** to copy the restoration type of last tooth to the currently selected tooth.
- Click **Shift+left mouse button** to copy the restoration type of last tooth to all the teeth between the last tooth and the currently selected tooth.

6. Select treatment type, including full crown, pontic, veneer and antagonist.




7. Click **Scan** to enter the scanning interface.


Start Scanning

Bite scanning can be put at first step. Click  and you will see the bite scanning points.


1. Scanning bites. A exclamation mark on the bite icon means alignment failure with collected bite data.




2. Scan upper jaw. Click  when scanning completed.



3. Scan lower jaw. Click  when scanning completed.




4. Scan left, right and front bite parts.



5. Align automatically. If it succeeds (✓ displays), click . If it fails (! displays), you need to re-scan bite parts.



Scanning Completed

Click  to complete scanning. Once the data is post-processed, the software will enter the interface of pre-design.

Implant Scan


Create an order


Implant scan takes a case of full crown as an example. Enter necessary information.

1. Choose patient types, **FirstVisit** or **FollowUp**.
2. Enter necessary information of the patient such as **name** and **age**.
3. Add a doctor/operator. Click  and enter the doctor's name, then click  to confirm.
4. Choose **implant** as the dentistry type. Add notes such as teeth shade.
5. Select the tooth and select **Full Crown** as the treatment type.
 - Click on the number of one tooth to select it (click right mouse button on the number to cancel the selection). Multiple selection is available.
 - Click **Ctrl+left mouse button** to copy the restoration type of last tooth to the currently selected tooth.
 - Click **Shift+left mouse button** to copy the restoration type of last tooth to all the teeth between the last tooth and the currently selected tooth.
6. Select **Custom Base** or **Screw Retained** to create a implant order.
7. Select material to restore the tooth according to the case.
8. Click **Scan** to enter the scanning interface.


Start Scanning


Steps of implant scanning: lower jaw (implant a tooth on the lower jaw) > implant on the lower jaw > upper jaw > full jaw.


1. Scan lower jaw (without intraoral scan body) and click  to complete lower jaw scanning.
2. Scan lower jaw implant.
 - (1) Enter the interface of lower jaw implant. Scroll the wheel up and down to zoom in and out the tooth model. Select the tooth place and make preparation for digging a hole.




(2) Double-click to dig a hole.



(3) Click  to lock data and start the process of scanning implants (the scan body must be inserted in the implant position before scanning).

 **Caution**


Locked data is displayed in dark green. Only the unlocked areas can be edited.



(4) Scan the scanbody.

Set and Scan the Scanbody


After setting the scanbody, it can be scanned.




Set the Scanby and Automatically Align

The material and shape of the scanbody (reflective metal and length, etc.) make it hard for the scanner to obtain the data of the whole scanbody. By aligning the scanned scanbody with scanbody model, the factory can have better design precision.

Steps of aligning the scanbody (Optional):





1) Click  to enter the Scan Body Matching window.


2) If there is no scanbody data, please download from the cloud database or click  to import local data.

- 3) Select the tooth for implant on the left.
- 4) Choose one type of scanbody and start matching. After matching, the standard scanbody data will be displayed on the right of the scan interface.
- 5) After setting the scanbody, start scanning it.
- 6) After scanning, the scanned data and the standard scanbody will be aligned automatically.

Set the Scanbody and Manually Align


Manual alignment can be adopted if the effect of auto alignment is not good enough.


- 1) Click  to delete all alignment on the scanbody (only scanbody data on the tooth can be retained).
- 2) Click  to enter the Scan Body Matching window.
- 3) Choose one type of scanbody and start matching. After matching, the standard scanbody data will be displayed on the right of the scan interface.
- 4) Click  to enter the manual alignment interface.
- 5) 3 groups of corresponding points were selected respectively on the left and right scanbody models.
- 6) Click  to finish manual alignment.

- (5) Click  and complete scanning lower jaw implants.

3. (Optional) Scan multiple implant teeth.

- (1) Enter the hole digging interface. Dig hole at the location of the implant and delete collected data.
- (2) Follow the guidance on the left panel and scan implants in groups. Try to stagger adjacent scan body. Insert the first group of scan bodies and then scan again. After scanning, click to pause.





- (3) Click  and add another group of implants.




Note


When you start scanning next group of scan bodies, the data of first group will be locked and displayed in dark green. You can not add scanning to the groups that have been scanned. Pause scanning and click next group to continue if there are more groups.

(4) Click  to complete scanning.




4. Scan bite parts. After scanning and saving upper jaw and lower jaw, the software will start bite scanning.

5. Click  on the right panel to edit the model. Details see [Data Edit](#).

6. Click  to confirm the changes and exit the edit interface.

Scanning Completed

Click  to complete scanning. Once the data is post-processed, the software will enter the interface of pre-design.




Dynamic Bite

Caution

- To collect data of dynamic bite, you need to collect data of static bite first.
- The areas of dynamic bite and static bite should be overlapped.
- The scanner tip should be at the center of upper and lower jaw when scanning.
- Make sure there are contacting points when biting.



Steps:

1. Click  to open dynamic bite.
2. Click  to start scanning.
3. Open mouth and scan the bite motion on the left and right.
4. Scanning completed, the video of dynamic bite is played automatically.
5. Click  and save data of dynamic bite

Note

You can move the cursor and rotate the model to view the bite condition from different perspectives.

Removable Denture Scan

Through scanning of removable dentures, the dentists can check the bite relationships of edentulous patients. In addition, the dentists can finish scanning within one order.

Create an order

More details can be found in [Create Order](#).

Note

- Bite scanning can be put at first step during the scanning process of removable dentures.
- [Multiple Bites](#) is supported when scanning bite.

Natural teeth bite

The suggested scanning process of natural teeth bite as the **Get bite record** is scanning upper jaw, lower jaw and bite.




After scanning, click  to go to pre-design.



Bite rim scanning

For upper jaw and lower jaw

The suggested scanning process is scanning upper jaw, lower jaw and bite rim.



After scanning upper jaw and lower jaw, click  to scan bite rim.

The step of scanning bite rim can be divided into  **Scan: scan bite rim** and  **Align: align bite rim with jaws**.

1. Start scanning bite rim.

Note


Please scan the whole bite rim for better alignment.


2. Align bite rim with jaws:

Auto alignment


The scanned upper jaw and lower jaw will be automatically aligned with the bite rim.

Scanned data will be displayed on the right.  indicates success, while  indicates failure.

Click  to hide/show scanned models.



Manual alignment


If something is wrong with the auto alignment, please select the upper jaw or lower jaw and click  to start manual alignment.

Steps:





- Rotate the scanned models to the position shown in the guide.
- Select 3 sets of corresponding points on the left and right models.




- Preview the effect.




Operations

Icons	Descriptions	Icons	DEscriptions
 Delete	Delete all points on the models.	 Undo	Undo the last selected point.
 cancel	Cancel the changes and exit.	 Confirm	Confirm the manual alignment.


3. Click  to complete the scanning.



For a single jaw

Take the upper jaw as example, the suggested scanning process is upper jaw, lower jaw, upper jaw rim, and bite.





1. Start scanning upper jaw rim after scanning upper jaw and lower jaw.

 **Note**

The step of scanning upper jaw rim can be divided into  **Scan: scan upper jaw rim** and  **Align: align rim with upper jaw**, which is similar to the steps mentioned above. References can be found in [For upper jaw and lower jaw](#).

2. Start scanning the bite of the upper jaw rim and lower jaw.




3. Confirm the alignment and click  to finish scanning.



Old denture scanning

For upper jaw and lower jaw


The suggested scanning process is upper jaw, lower jaw, upper jaw old denture, lower jaw old denture, and bite.





1. Scan the upper jaw and lower jaw.


2. Scan upper jaw old denture. The steps can be divided into  **Scan: scan old dentures** and  **Align: align old dentures with jaws.**

Scan upper jaw old denture:





Align upper jaw old denture with upper jaw:


Click  to automatically align the upper jaw old denture with the upper jaw. Click  to start manual alignment. More details can be found in [Bite rim scanning - For upper jaw and lower jaw](#).



3. After scanning upper jaw old denture, please start lower jaw old denture. The steps are the same with steps in scanning upper jaw old denture.

4. Click  to start scanning bite and bite alignment.



5. Click  to finish scanning.



For a single jaw

Take the upper jaw as example, the suggested scanning process is upper jaw, lower jaw, upper jaw old denture, and bite.





1. Scan the upper jaw and lower jaw.
2. Scan the upper jaw old denture.

Note

The steps of scanning the upper jaw old denture can be divided into  **Scan: scan upper jaw old denture** and  **Align: align upper jaw old denture with upper jaw**. More details can be found in [Bite rim scanning - For upper jaw and lower jaw](#).

3. Scan the bite of upper jaw old denture and the lower jaw. After scanning, the models will be automatically aligned.



4. Confirm the alignment and click  to finish scanning.


Extra Match

Extra-match can align the scanned data with third-party accessories for checking multiple sets of data.

Note

Extra match is supported in various scanning process, except **full jaw scanning** (including **scanning of full jaw bite rim**), **old denture scanning**, **bite rim scanning**, **implant scanning** and **extra scanning**.

Create extra align

Click  to choose Import/Scan, Impression match/Abutment match.

Tip ×


Please import or scan

Import Scan

Please choose alignment type

Impression match Abutment match

After choosing **Import** or **Scan**, the imported data or scanned data will be shown on the right of the interface. The data will be divided into 2 columns: impression data and abutment data.



Import

- After choosing **Import**, the user can choose STL, PLY, OBJ or BEB files from the local path.
- The maximum size of imported data is 500MB.

Scan

- After choosing **Scan**, the software will enter the interface of scanning impression or abutment.
- The user can add another group of data after scanning the former one.
- The scanning process of extra match is the same with other scanning processes.
- When scanning, **Edit**, **Remove isolated data** and **Part lock** are supported.

Other functions

Icon	Description	Icon	Description
	Manual alignment.		Rename.
	Display/hide the selected data.		Import or scan a new group of data.
	Auto alignment.		Click to delete the selected data.
	Display/hide all data.		Exit the interface of extra match and save current operations.


Multiple Bites


When scanning bite, users can add and scan multi-positional occlusions. By checking the multiple bites, the doctor can confirm the occlusal relationship.




Caution

Adding multiple bites will prolong the time of processing.

Add multiple bites

After entering the interface of scanning bite, click  to check the multiple bites.



- **Add:** Click  to add a new group of bites. A total of 6 groups of multiple bites can be added.
- **Rename:** Click  to rename the selected group of multiple bite.
- **Delete:** Click  to delete the selected multiple bite.

Note

- The default group of multiple bite is called **Maximum intercusp**. It can not be renamed or deleted.
- Dynamic bite is supported only for the default group.

Pre-design

On the interface of **Preview Edit**, users can check the multiple bites from different perspectives.

Note

If a group of multiple bite is not scanned, then it will be deleted in pre-design automatically.

More details can be found in [Preview Edit](#).

Send




The data of multiple bites can be sent or exported.

More details can be found in [Send Order](#).

Coded Scanbody Scanning

The high-precision coded scanbody, specially designed for edentulous implants, is used to locate and calculate the position and direction of implants in the 3D model. Regular replacement is necessary to maintain accuracy requirements.

Create an order

1. Click  to enter the interface of creating an order or click  to select a patient and create an order. More details can be found in [Create New Order](#).
2. Select **Implant** as the dentistry type.
3. Download the database of coded scanbody.
 - (1) Click **Coded scanbody**  below the tooth map.
 - (2) Download the database through **Enter the product serial number directly**, **Scan the QR code** or **Local import**.


Coded scanbody
✕

Please fill in the serial number of coded scanbody


Scan the QR code

Optional in case of no network [Local import>>](#)

1. What is a coded scanbody?
 The high-precision coded scanbody, specially designed for edentulous implants, is used to locate and calculate the position and direction of implants in the three-dimensional model. Regular replacement is necessary to maintain accuracy requirements.





2. How to fill in the serial number:
 Please fill in the serial number of the coded scanbody as shown in the picture





✓ Download

⚙ Go Settings

- (3) Click **Download** and check the relevant information of the coded scanbody.
- (4) Click **Confirm** to download the database.
4. Select upper jaw/lower jaw as the working jaw.
5. Select the tooth and select **Full Crown** as the implant type. The selected tooth is shown as .
6. Choose the implant systems.
 - (1) Click **Go choose** to enter the interface of choosing the implant systems. The interface shows the history, database and favorites list of the implant systems.
 - (2) After selecting the tooth, selecting the Manufacturer → Implant → Type → Sub Type in the database.

Or search a certain implant system by enter in its name or click  to import local database.

Click  to download database of implant systems.


 **Note**

When choosing an implant system, click ☆ to add it in the favorites list for future using.
- (3) Click **Apply**. If an implant system is suitable for multiple teeth, click **All Applications**.
- (4) Click **Confirm**.
7. Click **Scan**.

Caution

- Please download the database first.
- The number of the coded scanbodies for a single jaw should be no more than 9.
- To ensure product accuracy, the recommended product usage is 300 times or less.

Import an order

Click  to import local orders with scanbodies.

Note

- If the scanbody database of the imported order is not stored in local path, please download the database first.
- In the interface of editing the order or after clicking **Go Scanning** in Pre-design, a prompt pops up. Follow the prompt to download the database.

Scan

Take both upper jaw and lower jaw as working jaws with pre-op process selected, the suggested scanning process is:

① Pre-op scanning (including pre-op upper jaw, pre-op lower jaw and pre-op bite) → ② Upper jaw → ③ Upper jaw with implant → ④ Upper jaw gingiva alignment → ⑤ Lower jaw, lower jaw with implant and lower jaw gingiva alignment.

Note

- To scan the coded scanbody, Aoralscan Elite is required.
- Standard tip or big tip is required.

① Pre-op scanning


Scan the pre-op upper jaw and pre-op lower jaw with the temporary dentures, as well as the pre-op bite.


Pre-op upper jaw with the temporary denture



Pre-op lower jaw with the temporary denture



 **Note**


When scanning jaws with the temporary dentures, click  to enable **All-on-X scanning** for optimizing scanning process of temporary dentures.


② Upper jaw

Scan the edentulous upper jaw and align the temporary denture with the gingiva.

It is divided into  **Scan** and  **Align**.









Scan step

Scan the upper jaw. After that, click .



Align step




Static Alignment	Manual Alignment
<p>The scanned upper jaw will be automatically aligned with the temporary denture. Scanned data will be displayed on the right.</p> <ul style="list-style-type: none">•  indicates success.•  indicates failure.• Click  to hide/show scanned models.	<p>If something is wrong with the static alignment, please select the upper jaw and click  to start manual alignment.</p> <p>Steps:</p> <ol style="list-style-type: none">1. Rotate the scanned models to the position shown in the guide.2. Select 3 sets of corresponding points on the left and right models.3. Preview the effect. <p>Operations</p> <ul style="list-style-type: none">•  Reset: Delete all points on the model.•  Undo: Undo the last selected point.•  Cancel: Cancel the changes and exit.•  Confirm: Confirm the manual alignment.

③ Upper jaw implant


Please install the coded scanbodies first and then scan them.

Step 1: Install the coded scanbody

1. Prepare the coded scanbody and the special screwdriver.




2. Install a coded scanbody with appropriate length at the corresponding implant.



3. Adjust the orientation by gathering the ends of the coded scanbody towards the palatal side/lingual side. Try to contain 2 coded scanbodies within the scanner window.

4. Fix the coded scanbody with a screwdriver (About **10N·cm** is suggested).




Step 2: Scan the coded scanbody

1. Follow the path showing in the software to scan the whole structure of the coded scanbodies. The markers on the scanbody which have been scanned are displayed in green.
2. Follow the path showing in the software to scan each rod again in fine detail.




Note










Try to contain 2 scanbodies within the scanner window.

3. Click  to finish scanning.


④ Upper jaw gingiva alignment

Scan the connection between the coded scanbody and the gingiva to complete the splice of them.




Static Alignment	Manual Alignment
<p>Click  to align the scanbody with the jaw. The scanned data will be shown on the right.</p> <p> indicates success.</p> <ul style="list-style-type: none">  indicates failure. Click  to hide/show scanned models. 	<p>If something is wrong with the static alignment, please click  to start manual alignment.</p> <p>Steps:</p> <ol style="list-style-type: none"> 1. Rotate the scanned models to the position shown in the guide. 2. Select 3 sets of corresponding points on the left and right models. 3. Preview the effect. <p>Operations</p> <ul style="list-style-type: none">  Reset: Delete all points on the model.  Undo: Undo the last selected point.  Cancel: Cancel the changes and exit.  Confirm: Confirm the manual alignment.


After the alignment, conversion and mark of the scanbodies are needed.


 **Note**

Conversion and mark of the scanbodies can be finished in [Dental Cloud](#) as well.


Convert the coded scanbody

To match the coded scanbody with the implant and make sure that the scanned data is useful, the users can click  to converse the coded scanbody when the alignment is paused or finished.

1. Click  to enter the interface of conversion.



2. Select the type. Click Manufacturer → Implant → Type → Sub Type.

Or enter the name of the scan body to search the certain one. Or click  to import a local file.

Click  to download database of implant systems.

3. Confirm the coded scanbody. All coded scanbodies are selected in green by default.



Note


To converse part of coded scanbodies to other types of implants, please click other scanbodies to unselect them and click **Match** to converse remained ones.

4. Click **Match**.


Note


- The users can click ☆ to add the implant into the favorites list for further use.
- By rotating/zooming in and out the coded scanbodies, the users can check the effect after conversion.


5. Click  to save the current conversion and return back to the alignment.

Click  to reset the conversion.


Mark the coded scanbody


Click  to mark the coded scanbody.

1. Click  to enter the interface.
2. Select the tooth on the tooth map.
3. Select the corresponding coded scanbody.
4. Click **Mark**. The tooth label will be displayed on the coded scanbody.



5. Repeat step 3 to step 4 till all coded scanbodies are marked.


Or click  to mark all coded scanbodies.

6. Click  to save and return back to the alignment.

Click  to reset.







Scan the lower jaw, lower jaw implant and lower jaw gingiva alignment.

The steps of scanning the lower jaw, lower jaw implant and lower jaw gingiva alignment are the same with those of the upper jaw.

After scanning, click  to finish the whole procedure.








Scanning Tools and Additional Functions

Scanning tools

Scanning tools include  [Edit](#),  [Remove isolated data](#),  [Refined area](#),  [Part lock](#),  [Undercut](#) and  [Swap jaws](#).

Edit

Edit the model with the following tools.

Icon	Description	Icon	Description
 Rotate	Click Rotate and move the cursor to adjust the view of 3D model.	 Brush	<ol style="list-style-type: none"> 1. Click Brush and the cursor becomes a blue brush. Drag the slider to change its size. 2. Move the brush to the areas you want to delete. 3. Click and the area covered by the blue circle will be deleted.
 Free selection	<ol style="list-style-type: none"> 1. Click Free Selection. 2. Move the cursor and press the left mouse button to draw a closed curve. 3. Release the mouse and the selected areas will be deleted. 	 Undo	Undo the last operation. Multi-clicking for undoing multi-operations.
 Redo	Restore the operation that has been undone. Multi-clicking for redoing multi-operation.		Cancel
 Confirm	Confirm the operation.		

Remove isolated data

Delete smaller and isolated model data that are not connected to the main model. Help delete unrelated model data quickly.


Refined area

The area will display more partial details in a higher resolution grid during the post-op process.

This feature is best used with Tooth Color function to ensure that the quality of the tooth color is true in the refined area.

Select **Refined Scan** in the **Settings** -> **Scanning Settings** to automatically enter the process of refined area selection when the upper/lower jaw scanning is completed.

Part lock

When the scan is paused during the scanning process of **pre-op or upper/lower jaw**, click  **Part Lock** to set the scanning lock area to ensure the data in the locked area will not change in subsequent scanning.

Undercut

It is used to check and calculate the undercut area on the marked teeth (in accordance with the order).







During the scanning process of upper/lower jaw, pause the scanning and click **Undercut**.

The area selected by the brush will be calculated automatically. Details see [Undercut](#).

Swap jaws


Switch data of upper and lower jaw after the scanning is completed or post-processing is completed. This function is used when the upper and lower jaw is reversed during scanning.


Additional functions

Icon	Description
 Texture	Set whether the model shows colored model or not.
 Data indicator	Turn on the tooth color display. If the gray color is displayed when scanning, it means the data quality in these area is not good enough, and the area should be scanned again. When the grey color turns into true color display, it means the data quality has reached the requirements of better post-processing.
 Intelligent Scan	Applied to intraoral scanning . After clicking it, miscellaneous data will be automatically deleted when scanning, such as data of buccal, lingual and soft tissues. When you use the software for the first time and enter the scanning interface, the button is gray and you need to wait for 1-2 minutes for initialization. After the initialization is completed, the button can be clicked.
 View Lock	Check the model from one fixed perspective when scanning, and the automatic perspective tracking function is not enabled. This function is not enabled by default.
 Metal tooth scan	Improve the scanning speed and sweep the whole situation for metal teeth. 


Send Order

Upload the order to the cloud or a third-party platform for the download and use by other authorized institutions.

On the interface of pre-design, click **Send** to send the order. The users can also click  on the order list or double-click a certain order to send it.




Steps

1. On the interface of sending an order, the users can check order information, tooth list, the scanned model and change send settings.
2. Check whether to send back the data. When it is enabled, the users can download and check relevant data which is uploaded to the dental cloud.
3. Click  to check **Send EXO order** and **Send 3rd party order** or not.
4. Set the **Send Settings**.
 - (1) Select the target institution in the drop-down list.

Note

Unconnected institutions are not displayed.

- (2) (Optional) Click  to add attachments.
5. Click **Send** to send the order.




Pre-design

Preview Edit

Observe the scanned model from different angles and then optimize scanning data.

Intraoral data

- Change the opacity and switch the viewing angle to edit the stereoscopic data.
- If multiple bites are added and scanned, the users can switch between different groups of multiple bites to check the models.

Icon	Name	Description
	Upper Jaw	Check the option to view the upper jaw data. Drag the slider to change the upper jaw opacity value.
	Lower Jaw	Check the option to view the lower jaw data. Drag the slider to change the lower jaw opacity value.
	View from Different Perspectives	Check the view of the model.

Application list

The application list shows multiple functions in Pre-design.








- Click **More Applications** to check more on the pop-up window.
- Click **Adjust sequence** to adjust the sequence of applications showing in the interface.
 - (1) Drag the icon of an application to adjust the sequence.
 - (2) Click **Apply** to save the adjustment.
- The last used application will be displayed before **More Applications** by default.



Edit



Click  to edit the model with the following tools.

Icon	Name	Description
	Rotate	Click and rotate the model.
	Brush	<ul style="list-style-type: none"> • Hold left mouse button and move to select areas and delete them. • Drag the slider to change the brush size.
	Free Choice	Hold the left mouse button to select. The area selected will be deleted.
	Undo	Clicking for undoing operations.
	Redo	Clicking for redoing operations.
	Cancel	Discard the changes.
	Confirm	Confirm the changes.






Remove Isolated Data

Remove Isolated Data can quickly delete small data isolated to the main model.

Fill Holes


Steps

1. Click the icon  to enter the **Fill Holes** interface.
2. Drag the slider to adjust the application range.
3. Click  to confirm the operation and exit; click  to discard the process and exit.




Texture

Enabled by default. When enabled, the model will be colored.




Enabled




Disabled

Smooth

Click  to reduce noise and improve the quality of the model.

Fit View


Click  to adjust the size of the model to be adaptive to the interface.

Oral Health Report

Create an oral report

Steps

1. On the pre-design interface, click .



Settings Of Examination Report

Official version [Enter beta version](#)

Hospital

Target Adult Child

Template Type

Note

A beta version is available. Click **Enter beta version** to get new features. More details can be found in [ConsulReport](#).


2. Select a hospital, a target (Adult/ Child) and a template type.
3. Click **Accept** to edit the oral report.

Note

- When enabling assistant mode, the software will start assistant identification automatically before editing the oral report.
- Assistant Mode is only accessible when the target is adult.

Edit the oral report

1. When enabling assistant mode, click **Accept** to enter the main interface to edit the oral report.
2. The main interface displays the bite status of the jaws. Click the picture of a symptom and the model will show the position of the symptom.




Note


- Import an order which was created in the software with lower versions and click oral report, and then a tip is popped up to confirm whether to generate a new report. Click **Yes** to clear operations and enter the interface of editing a new report. Click **No** to preview the report.
- Assistant mode is not enabled for preview-only reports.



Assistant Mode

- A tip of **Abnormal data not identified yet** is popped up when the software can't identify diseases. Otherwise, the identification results are shown in the **Issue** window.




- The results are watermarked with AI. Click the picture to zoom it in.

















Operation	Description
Show description /Show notes and plan	Check to show the description, notes and plan of the symptoms.
Select best (all)	Check to select the first one of the symptoms.
Reset option	Click and a tip of "Whether to clear all operations?" is popped up. Click Yes to delete all selections and notes.
Hide this symptom	Check to hide this kind of symptom.
Accept	Save all the selections and enter the interface of editing oral report.
 Select teeth	Click to select the tooth of the symptom.
 Delete	Click to delete the certain symptom.









- Select no symptom and click **Accept**, then a tip of "Symptom select incompleted, whether to nex step?" is popped up. Click **Yes** and unselected symptoms are deleted.

 **Note**



Selected symptoms are displayed in the main interface of oral report, while the unselected symptoms are completely deleted.

Main Interface of Oral Report

Icon	Description
One-click filling	Click to set descriptions and plans in issue view.
Delete	Click and a tip of "Whether to delete all issues?" is popped up. Click Yes to confirm.
 Preview	Preview the oral health report. Share: Create a QR code which can be scanned by the phone to view the report. Export: The report will be exported in PDF and PNG.
 Assistant Mode	Assistant Mode will automatically detect teeth problems and take screenshots of certain teeth. Click the button when there are already identifications, a prompt of "Whether to perform auxiliary identification again?" is popped up. Click Yes to restart auxiliary identification and current results are deleted.
 Endoscope	Take a screenshot of the model. <ol style="list-style-type: none"> 1. Click and the enlarged image window is displayed on the left. 2. Scan the teeth. 3. Press scanning button on the scanner or click  on the image window to take a screenshot. 4. Select the treated tooth and problems in the new pop-up. Adding notes is also supported. 5. Click  to take more screenshots.
 CT Image	Add an image and edit it. Treatment plan and other information of the treated teeth can be imported. Change image: Click to change the image. <ul style="list-style-type: none">  Mark: Click the treated tooth to mark it with different colors and enter relative information.  Plan: Click plan and a box is displayed. You can enter your treatment plan here.  Brush: Click to paint directly on the image. Drag the slider to adjust its thickness.  Eraser: Clear the marks in the model.  : Multi-click for undoing multi-operations.  : Multi-click for redoing multi-operations.  : Discard the changes.  : Confirm and save the changes.

Icon	Description
 Grab image	<p>Take a screenshot of the scanned model.</p> <ol style="list-style-type: none"> 1. Click the button of screenshot and an enlarged screenshot box is displayed. 2. Press left mouse button within the black box and move the cursor to rotate the model. 3. Press right mouse button outside the box and move the cursor to rotate the model. 4. Press right and left mouse button outside the box and move the cursor to move the model. 5. Scroll the wheel to zoom the model in and out. 6. Double click in the box to confirm and select the treated tooth and problems in the new pop-up. Adding notes is also supported. 7. Click  to take more screenshots.
 IR grab	<p>Take a screenshot of the near-infrared image.</p> <p> Note</p> <p>If the order has activated the near-infrared image, the IR grab function can be displayed.</p> <ol style="list-style-type: none"> 1. Click the button and an enlarged screenshot box is displayed. 2. Click  to adjust brightness and contrast. 3. Press left mouse button within the black box and move the cursor to rotate the model. 4. Press right mouse button outside the box and move the cursor to rotate the model. 5. Press right and left mouse button outside the box and move the cursor to move the model. 6. Scroll the wheel to zoom the model in and out. 7. Double click in the box to confirm and select the treated tooth and problems in the new pop-up.
 Graffiti Mode	<p>Display models in various angles including close jaw, open jaw and CT.</p>
 Recycle Bin	<p>Deleted symptoms are collected in the recycle bin.</p> <p>Click the Restore button on the top right corner of the picture to restore it.</p> <p>Click the recycle bin button again or click Exit on the top right corner to exit.</p>
 Exit	<p>Click to exit the oral report.</p>


Other functions

Name	Description
3D label	<p>Not enabled by default. When it's enabled, symptoms under this perspective will be tagged. Click the tag and the picture of this symptom will be selected on the right. Modifying notes, descriptions and plans, selecting tooth and deleting this symptom are supported.</p> <p> Note</p> <p>When the jaws are open, the 3D tags are hidden if the software can't determine whether the symptoms are on upper jaw or lower jaw.</p> 
Opened jaw	Not enabled by default. When it's enabled, the model is displayed in open jaw.
Direction flag	Not enabled by default. When it's enabled, an indicator avatar will be displayed in the bottom left corner. The indicator avatar will change synchronously as the model is moved, providing users with a more intuitive view of the model.
Graffiti mode	<p>Not enabled by default. When it's enabled, the cursor turns into brush.</p> <p>Press the left mouse button and move the cursor to draw on the model.</p> <p>Press left mouse button and move the cursor to rotate the model.</p>

Coordinate Adjustment


Adjust the position and the angle of model from three views, until the model being approximately horizontal and covering the central shadow area of the coordinate.

Click  under the "pre-design" process to enter the interface.



You can choose to manually or automatically adjust the coordinates.

Instruction	Short-cut Key
Rotate model	Press and hold the left (or right) mouse button, and meanwhile move the cursor.
Move model	Press and hold the left and right mouse button, and meanwhile move the cursor.


Click  and the model will be placed in a proper position.


Click , the system is to restore the model to the initial position.

After adjusting the coordinates, click **Go to send** to upload data or implement other operations.

Mark Teeth

Mark teeth is used to manually mark the corresponding position of the teeth specified in the scan order in the scanned model for subsequent extraction of the margin line and detecting of the undercut.

Click  under the "Pre-design" process to enter interface on which you can mark the teeth, and the tooth map on the left side shows the teeth defined by the order.



Note

IntraoralScan will automatically recognize teeth marks from the imported project.

Steps:


1. Select the tooth needed to be marked on tooth map.
2. Move the cursor to locate the corresponding tooth on the scanned model, double click it to mark.
3. Repeat step 1 and step 2 if you want to mark two or more teeth.
4. Double click the right tooth on the model to change the mark.
5. To delete all teeth markers, please click **Delete** on the left side.

Note

To delete single tooth marker is not supported currently.

Bite Detection

Check if it is a normal occlusion and adjust it. Adjust and view the occlusal gap between the upper and lower jaws.





Click  under the "Pre-design" process to enter the occlusion detection interface.

Color Bar Mode


The occlusion status is shown by the color bars.

- Click the button on the left to change occlusion parameters.
- Double click one touching point to detect the gap.

Color	
Green	The gap
Red	The touching area
Blue	The bite-through area

Description	Operation
Click 	Exchange the upper jaw and the lower jaw.
Click 	Open or clamp jaws for imitating the biting action.
Calculate	Click  to calculate the bite amount. Move the cursor to the calculated area to check the value.
Reset	Click  to reset the adjusted value to the initial value.

Action	
Double Click	Check the specific occlusion value.
Click Up / Down Buttons on the Left Bar	Adjust the numerical range to occlusion.




Closed

Open

Section Mode








Click  to enter the section mode, from which the occlusion status can be easily checked (the positions of teeth in upper jaw and lower jaw). Manually select two points on the section view to measure the distance.

- In this mode, a round section with a diameter of 30mm will be created on the 3D model automatically.
- Click the small sphere in the center of the section and the the arch line will be displayed. Hold the sphere and move the cursor, the section view on the left will show the 2D intersecting part of the section with the model in real time. Click on the other part of the section and hold, then move the cursor, the section will be rotated.
- The section view on the left will show the 2D wireframe of the intersecting part, from which the occlusion status can be easily checked.


Note

Manually select two points on the wireframe to measure the distance. When selecting a third point, the second point is the start point while the third one is the end point.

Icon	Name	Description
	Vertical Arch Line	After enabling it, when moving the center sphere, the section will be vertical to the arch line.
	Restore the Default Scale	Scroll up and down the mouse wheel on the section view, and the wireframe will be zoomed in and out. Click  to restore the default scale.
	Zoom in	Click to zoom in the section view. Click again, the view will be restored.
	Clear Selected Points	Click to clear all selected points.

Extract Margin Line

Extract and save the margin line of the marked (corresponding) teeth.

Click  under the "Pre-design" process to enter the margin line extraction interface.


Caution


Please mark the teeth before extracting its margin line. More details see [Mark Teeth](#).

Note

IntraoralScan will recognize margin lines automatically when importing projects.

Auto Drawing


Click  to extract margin line for teeth with apparent edge line automatically.









Manual Drawing

Steps:

1. On the teeth position marking list, select a tooth.
2. Under tooth map, click to select margin line extraction method.



Approach	Description
 <p data-bbox="194 360 312 427">Sectional view</p>	<p data-bbox="427 282 1374 349">The curvature of the profile along the tooth surface normal to the occlusion plane is displayed, and the effect is as follows.</p> 
 <p data-bbox="194 768 352 835">Drawing line extraction</p>	<p data-bbox="427 692 1445 801">Long press the left mouse key to draw a line, and release the mouse to extract part of the margin line, each drawn line is automatically connected with the previous one, and finally a closed margin line is extracted.</p>
 <p data-bbox="194 978 320 1046">Segment extraction</p>	<p data-bbox="427 902 1453 1093">Left click the mouse to generate a control point, move the mouse along the margin of the tooth, a line segment will be automatically generated between previous control point and current mouse position, left click to generate another control point; repeat the above steps before you set the final control point to overlap with the first one to form a closed margin line by double clicking.</p> <p data-bbox="435 1104 549 1171">  Note </p> <p data-bbox="427 1182 1390 1249">Line segment extraction is a complementary tool for magnetic extraction, suitable for teeth with no clear margins.</p>
 <p data-bbox="194 1388 312 1456">Magnetic switcher</p>	<p data-bbox="427 1312 1366 1379">Used to control the magnetic suction function switch of the line segment extraction method.</p> <p data-bbox="427 1391 1449 1749">After the magnetic suction function is turned on, the line between the current position of the mouse and the previous control point will be automatically attached to the margin (large curvature), the attaching effect is related to the current mouse position; when you left click to add a new control point, the attached line segment between the first two control points will not change again, repeat the above steps before you set the final control point to overlap with the first one to form a closed margin line by double clicking. When the magnetic suction function is off, a straight line will be created between the current mouse position and the previous control point (i.e. line segment extraction). Magnetic suction extraction is suitable for teeth with clear margins.</p>




3. Click to add a new control point and a margin line will be drawn automatically between two control points. Repeat the above steps and you will draw a closed margin line by double clicking near the first control point.
4. The model size can be adjusted by scrolling up and down the wheel. You can view the margin line by pressing right mouse button then rotating the model, or pressing right and left mouse buttons to move the model.

5. Edit margin lines to make them closer to tooth edges.

- **Draw line extraction mode:** The newly drawn line replaces part of originally extracted margin line.
- **Segment extraction mode:** Move the mouse to the extracted margin line, the red point indicates the location of the control point, if you left click the mouse key, a new control point will be added; select a control point, the control point color changes from green to red, you can move or delete the control point.

6. ✓ on the left of the tooth list indicates the line has been extracted.


Buttons

Icon	Name	Description
	Delete	Delete the currently extracted margin line of the tooth.
	Undo	Undo the last operation.
	Redo	Restore the operation that has been undone.




Modified Model

After marking teeth and drawing margin lines, you can set the undercut direction, fill the undercut and sink margin lines on the **Modified Model** interface.

Click  to enter the Modified Model interface.










Steps:

1. In the upper left tooth mark list, click a tooth in the upper or lower jaw to adjust the margin line again, if needed.
2. Check one or more teeth in the tooth mark list in the upper left corner.
3. Click  to enter the **Undercut Direction** interface to set the undercut direction.
4. Click  and set the value at the bottom of the interface to fill the undercut.
5. Click  and set the value at the bottom of the interface to sink the margin line.

Note

After filling the undercut or sinking the margin line, drag the slider in the upper left corner of the interface to view the comparison.


Button


Icon	Name	Description
	Undercut Direction	Click  to enter the Undercut Direction interface to set the undercut direction.
	Delete Undercut Direction	In the tooth mark list, click  to delete the corresponding tooth mark.
	Undercut Direction Group	In the tooth mark list, the icon will appear after checking two or more teeth to set the undercut direction all at once.
	Undercut Filling	Set the value at the bottom of the interface to fill the undercut
	Margin Line Sink	Set the distance and expansion for margin line sinking.




Undercut

Undercut is used to calculate and view the undercut area of the marked (corresponding part in the order) teeth.

After the margin lines are detected, click  to enter the undercut interface.






Steps:

1. Choose one tooth on the upper jaw or lower jaw.
2. Double click left mouse button on the model or use **brush** to select the undercut range.
3. Click  to calculate the undercut.
4. Move cursor can view the undercut value.

Note

Adjust the range of the undercut value through the bar on the left.

Buttons

Icon	Name	Description
	Brush	Click the Brush and set the range of undercut generation. Drag the slider to change its size.
	Calculate the undercut	Display the undercut effect.
	Auto-fill	The model interface will automatically be adjusted to fit the size of window.



Dynamic Bite

Dynamic bite refers to the bite relationship between upper jaw and lower jaw in occlusion. Normally, other dental scanners collect data of static bite, which may pay less attention on occlusal problems with the impact of muscle stretching, joint motion and fissures.

Bite movement is displayed by videos.



After completing the [dynamic bite scanning](#), click  and enter dynamic bite interface to preview the bite movement.

Icon	Name	Description
	Occlusion Detect	Check the occlusion situation.
	Open Jaw	Choose contacting points in upper jaw and lower jaw for malocclusion cases.

ConsulOS

Overview

ConsulOS is a software targeted at orthodontic simulation of teeth data, which can be scanned by IntraoralScan and FScan.

Features

1. Automatically simulate the change of the face and teeth before and after the orthodontic setup;

2. Switch flexibly between the teeth model, providing different perspectives to observe the simulation effect;
3. The enhanced tooth texture makes the model closer to the real teeth of the patients;
4. Upload multiple programs and generate various reports of orthodontic simulation at a time.

Quick guide

Create program

Creating programs after importing teeth models, adjusting the model, confirming teeth and the program are supporting.

→ [Create program](#)

Preview program

Generating orthodontic reports, previewing orthodontic process via animation, and comparing the models before and after orthodontics are supported.

→ [Main interface](#)

Manual setup

Adjust the dental arch and teeth to customize individual teeth.

→ [Manual setup](#)

Add brackets

Simulate the real orthodontic effects with brackets.

→ [Add brackets](#)

Create program

When the user opens ConsulOS from IntraoralScan or FScan, the software will create the program based on the order from IntraoralScan or FScan. The user can add intraoral data in FScan.

Caution

- Select **Orthodontics** as the **Dentistry Type** to enable the orthodontic simulation function.
- Lip line should be extracted in advance in FScan. More details can be found in [Lip line extraction](#).

There are three steps: [Model adjustment](#), [Teeth confirm](#) and [Program confirm](#). [Change program](#) is also supported.

Click **Auto creation** to start.

Note

- Click **X** on the top right corner to exit the window of creating program. And the original model is displayed on the main interface.
- Click **+** on the top of the main interface again to create more programs.
- Up to three new programs can be created.

Model adjustment

In this window, the software will automatically adjust the model and show multi-views. The user can manually adjust the perspectives of models.


Teeth model only


There are four views of the teeth model.


- If the user adjust the model in one view, the other views will be changed correspondingly.


Create program ✕

1 Model adjustment >2 Teeth confirm >3 Program confirm
















↺↻↷

Next step

Operation


Operation	Description	Operation	Description
Rotate the model	Press and hold the left or right mouse button and move the cursor.	 Guidance	Click to pop up the guidance. Click  to set whether to automatically pop up guidance.
Move the model	Press and hold the left and right mouse button at the same time and move the cursor.	 Reset	Click to restore the models.
 Zoom in	Click to zoom in the view.  Note Zoom-in can be used when importing teeth model only .	 Undo	Click to undo operations.
 Redo	Click to redo operations.		

Teeth confirm

In this window, models of upper jaw and lower jaw are displayed separately. The software automatically recognizes each tooth with different colors and displays teeth numbers.

Users can preview the upper jaw and lower jaw from different perspectives, check the teeth numbers, edit teeth areas, remove teeth and add teeth. Click the number on the **Table of contents** to check the certain function.

Create program



① Views

Click to change the views of the model.

② Tooth number

Enabled by default. When enabled, the teeth numbers are displayed.



With tooth number



Without tooth number

③ Texture

When enabled, the model is colored.



With texture



Without texture

④ Area

Enabled by default. When enabled, teeth are differentiated from each other with different colors.








With area




Without area

⑤ Modify Area


Click  to enter the interface of tooth selection. Select a certain tooth to modify its area.


-  : Click to select the tooth area. Drag the slider to adjust the thickness of the brush.
-  : Click to delete the selected area. Drag the slider to adjust the thickness of the eraser.
-  : Undo the last operation.
-  : Redo the last operation.

- **Draft:** Click to save the current operation and return to the tooth selection interface to select other teeth.
- **Cancel:** Cancel all operations and exit Modify Area.
- **Confirm:** Save all operations and return to the teeth confirm interface.


 **Note**

When tooth regions are wholly erased, a tip of "Delete all teeth areas is not allowed" will pop up.

⑥  **Change tooth number**

Click  to change tooth number.



- **Cancel:** Click to cancel all operation and exit.
- **Confirm:** Click to confirm and save all operations.

 **Notw**

The red number means tooth position is repeated. When confirming, a tip will pop up.


⑦  **Remove teeth**

Steps


1. Click  to enter the interface of removing teeth.
2. Select the tooth which should be deleted and the tooth number is displayed in red.
3. Click  to delete it. If more than one tooth needs to be deleted, repeat steps 1-3.
4. Click **Confirm** to save the change, or click **Cancel** to cancel all operations and exit.

⑧  **Add teeth**

Steps

1. Click  to enter the interface of adding teeth.

2. Double-click the position where the tooth should be added.
3. Select the tooth number in the pop-up window.
4. Click **Confirm** to save, or click **Cancel** to cancel all operations and exit.
Repeat steps 1-4 to add more teeth.


 **Note**

The tooth number should be different.

5. If the tooth area is not satisfied, please double-click the tooth area to re-select or modify the area.

Program confirm

In this window, the user can change the program name, choose the shape of dental arch and select tooth types.





1. Click the name box to change the program name.
2. The shapes of dental arch are: **Oval**, **Squirrel** and **Vesica piscis**.
3. Tooth types: normal tooth, extracted tooth, virtual tooth and missing tooth.

Note

Click **Reset** to reset all tooth types.

Change program

To change program, please click  to enter the interface of **Model adjustment**. Click  to enter the interface of **Teeth confirm**.

Caution

After changing the program, the software will clear all programs and a tip will pop up. Click **Continue creating** to confirm the changes.

Orthodontic Simulation

Orthodontic simulation


After creating the program, the user can enter the main interface.


The user can generate the QR code of the simulated orthodontic report, preview the orthodontic effect via animation, and compare the programs with the original states. In addition, the user can [adjust the setup manually](#) and [add brackets](#) as well.


Generate report


On the report, the patient can check the 3D models before and after the orthodontic simulation, and preview the process via a video.

1. Click **Report** on the top right corner of the main interface.





2. Click  to upload the programs and generate the orthodontic report. Saving the QR code or copying the link for future use are supported.

3. Click  to check the QR code again.


4. After changing the program, click  to update the report.

Animation

Animation is supported for previewing the orthodontic process and final effect. Click  to play the animation.

Click  to loop the animation.

Note


Preview the animation of one program is supported. If comparing 2 or more programs, click  and then a tip will pop up to select a program.

Compare the programs

Comparing programs and the original state is supported.

Selecting the programs (the original state can also be selected) on the top of the main interface to start comparison.









Up to 3 new programs and the original state can be compared at the same time.



Other Functions

On the main interface, there are a few functions to better preview the models.

General functions


Name	Description	Effect
 Overlay	Not enabled by default. Click to overlay the model of the original state on the simulated model. The original model will be displayed in blue.	
 Tooth position	Not enabled by default. Click to show tooth number on the model.	
 Brackets	After adding brackets, click to show brackets on the model.  Caution Please add bracket in advance. More details can be found in Add brackets .	
 Occlude	Not enabled by default. Click to check the occlusion status with the color bars.	
 Texture enhancement	Click to change the effect of the simulated model to the effect of the scanned model. The simulated model effect is enabled by default.	
 Interproximal contact	To check the interval between the adjacent teeth. Not enabled by default. To adjust the interval, please start Manual setup .	


Manual setup

Note

- When comparing the original state with programs or comparing programs, **Manual setup** is supported. Please select the certain program on the pop-up window.
- In the original state, if the user adjusts the teeth and confirm the adjustment, a new program will be created. If there is already 3 programs, then manual setup is not supported in the original state.

Through manual setup, the user can adjust the simulated orthodontic effect and customize individual teeth according to specific requirements.

Click  on the top of the main interface when selecting only **one program** to enter **Manual setup** interface.




After adjustment, click **Confirm** to save the changes.

Dental arch

Manually adjusting indexes of the dental arch by dragging the slider is supported.

When checking the teeth model, click  to adjust indexes of the dental arch.




Adjust teeth

Use  Shift,  Rotation and  Torsion to adjust the teeth.

Drag the slider to adjust the sensitivity.


Interproximal contact

Click  to check and adjust the interval between the adjacent teeth.

- There is a point on the contacting part of adjacent teeth. Different colors indicate different intervals.
- Click one point and adjust the interval by **+** and **-**.
- The tags showing the intervals can be pressed and dragged to anywhere to avoid blocking the model.


Note

Only ****Tooth Model**** can check and adjust interproximal contact.






Add brackets

Through adding brackets, the user can simulate the real orthodontic effects with brackets.

Click  on the top of the main interface when selecting only one program to enter **Add brackets** interface.

Steps:





1. After entering the **Add brackets** interface, click  to add.
2. If there are existing brackets, click  to replace the original brackets, or click  to delete the original brackets.
3. Choose the bracket type and brands.

Add brackets ✕

Bracket type

Self-ligating metal bracket ▼

Bracket brands

-  **PROTECT** 普特 PROTECT
-  **Shinye** Shinye
-  **XIHUBIOM** 西湖·尔 XIHUBIOM
-  **3B** 三比正萌 ORTHO Xingchen 3B

Cancel Confirm

4. Click **Confirm** to apply the brackets.

AccuDesign

Introduction

AccuDesign is a model generation software developed independently by Shining 3D, and it is used for editing scanned data to export digital models. Besides, featuring the combination of automation and personalization, the software also supports adding specially designed attachments to digital models living up to 3D printing requirements, mainly applying in dental implant and orthodontics.

Note

- The software can be used independently, but it only supports importing model files in STL, PLY, OBJ, and BEB formats, or project files in ACCUDESIGNPROJECT format.
- Additionally, it also supports being called from scanning software such as IntraoralScan, from which the direct use of the model data is allowed, but it requires the use of Shining 3D intraoral scanners or Shining 3D desktop scanners.

Import the Model

Enter the **Welcome to AccuDesign!** interface, where you can choose to **Create new project** or **Import project**.



Note

If AccuDesign is called from other scanning software (such as IntraoralScan), the software will load the models from the current order. Deleting the models or importing other models are not supported; only 3D preview of the current models is available.


Create a New Project

Click **Create new project**, import model files (in the format of STL, PLY, OBJ or BEB) locally and select **Base**, **Type of honeycomb**, **Articulator** and **Text**.

Import the Model

- Click **Upper Jaw** or **Lower Jaw** import frame to import models.
- The imported model is displayed in the default occlusal view. You can perform operations such as translation, rotation, and scaling on the model (as shown in the table below). In addition, you can click the button to delete the model or click the button to reset the model to its initial state.




Order	Shortcuts
Move	Move the curser while hold down the left and right mouse buttons
Rotate	Move the curser while hold down the right mouse button
Scale	Rotate the wheel up or down

 **Note**

- It supports importing non-segmented models, which means it can import data for a single upper jaw model, a single lower jaw model, or a model including both upper and lower jaws.
- It is recommended to import the correct upper jaw and lower models that match the diagram information within the specified frame, otherwise it may affect the automatic generation of the model.

Base

Support three base types, namely **Model with foundation** (as default), **Model without foundation** and **Orthodontic ABO**; the description of three types is shown in the table below.


Base type	Demonstration	Commons	Differences
Model without foundation		/	/
Model with foundation		The base is a bottom mesh and has drainage holes, and its default thickness is 2mm	Pull the base directly downward from the model edge to extend it. The default height of the base is 3 mm
Orthodontic ABO			The default height of the base is 13 mm. A matching dental arch for the base is added, with a basal angle of 70°

Note


You can hover the cursor over  button to view the guide chart.

Type of honeycomb


Support four types of honeycomb, namely **Baseplate honeycomb** (as default), **Throughout honeycomb**, **Solid** and **Hollow**:



Baseplate honeycomb



Throughout honeycomb



Solid




Hollow

- The usage of print materials: Hollow < Baseplate honeycomb < Throughout honeycomb < Solid.
- The hardness of printed model: Solid > Throughout model > Baseplate model > Hollow.

In a word, when more printing materials are consumed, the hardness of the printed models can be ensured; please select the honeycomb structure according to the actual situation.

Note

- If you select **Baseplate honeycomb**, you should set the printer type: **Shining 3D's printers** or **Other printers**.
- You can hover the cursor over  button to view the guide chart.

Articulator


Support multiple articulator options (the default option is **Empty**): large/medium/small articulator, articulator for quadrant models, etc.


Note

Articulator options can not be applied when only one model is imported.


Text

Support customizing **Text Content** and **Text Position** (Left, Middle, Right).

 **Note**


- The maximum chapter limit is 100.
- You can hover the cursor over  button to preview the position of the text.

Click **Auto generation**, and wait for the software to automatically design and generate the digital model.


 **Note**

- When the interface displays "Processing, please wait...", please wait until the interface finishes loading.
- If a popup window appears with the message "Generation failed, please readjust!", please follow the guidance on the right side of the interface for adjustments:
- If the automatic alignment of the model fails, then you will be automatically navigated to [Coord adjust](#).
- If the recognition of the model area or shell extraction of the base fails, then you will be automatically navigated to [Area adjustment](#).
- If the addition of customized text or articulator fails, then you will be automatically navigated to please proceed to [Attachment adjustment](#).

Import a Project

Click  **Import project**, import project files (in the format of ACCUDESIGNPROJECT) locally.

1. Click **Import project** will pop up a **Import project** dialog box.
2. Double click a project file in the dialog box, or click to choose a project file and click **Open** button to import the project.


 **Note**

If the imported project has not yet proceeded with all steps, then you will be automatically navigated to the unfinished steps; or you will be automatically navigated to [Printing view](#).






Edit the Model

Coord Adjustment

On the interface of **Coord adjust**, adjust the position and the angle of model from three views (front view, top view, left and right side views) until the model is horizontally aligned in all three views, and meanwhile roughly overlaps with most darker colored shadow area on the coordinate system, to ensure the model is positioned at the center of the coordinate system.




Button introduction


-  : Create new project.
-  : Import project.
-  : Export project.
-  : Export model.
-  : **Options** button. Click the button to expand the dropdown list.









Note

If AccuDesign is called from other scanning software (such as IntraoralScan), the functions mentioned above are not supported.

-  : **Setting** button. Click the button to open the **Settings** popup window, where you can **Select Language**, set up **Model thickness** and **Auto pop-up guidance** (switched off by default).

Note

- Support setting language as Chinese or English. If AccuDesign is called from other scanning software (such as IntraoralScan), the language corresponds with that in the scanning software and the language setting is not displayed.
- Support entering number directly or clicking  button to increase or decrease the thickness (mm) of models, and the setting will be effective after restarting the software.
- The instructional animated image will be displayed by default on the right side of the step interface when the **Auto pop-up guidance** is enabled.

- : **Help** button. Click the button to show or hide the instructional animated image.
- : Enlarge button. Click the button to enlarge the instructional animated image; and click  in the top right corner to shrink the image.
- : Enlarge button. Click the button to enlarge three views to help for delicate operations; and click  in the top right corner to shrink the image.
- : **Reset** button. Click the button to reset the model to its initial coordinates.
- : **Undo** button. Click the button to undo the previous action performed on the model; you can click it multiple times to undo multiple actions in succession.
- : **Redo** button. Click the button to redo the previous action performed on the model; you can click it multiple times to redo multiple actions in succession.

Operation introduction

Support moving or rotating the model, as shown in the table below.

Order	Shortcuts
Move	Move the cursor to any area of the three views, then hold down the left and right mouse buttons, or hold down the mouse wheel button, while move the cursor.
2D Rotation	Move the cursor to any position within the light grey area of the three views, then hold down the left or right mouse button, while move the cursor.
3D Rotation	Move the cursor to any position within the dark grey area of the three views, then hold down the left or right mouse button, while move the cursor.

After the model alignment, click **Next step** to proceed with manual operation or click **Auto generation**.











Area Adjustment

On the interface of **Area adjustment**, use **Brush**, **Transmission** and other tools to manually select area to be printed.




Note

When adjusting the area, you can hold down the right mouse button to rotate models or rotate the wheel to zoom in/out, which helps in selecting the area more precisely.

Tool/Button	Instructions
 Brush	<ol style="list-style-type: none"> 1. Click the button, and the button will switch to  ; in the operation interface, the dark gray circle at the cursor represents the coverage area of the brush, and you can use the slider to enlarge or reduce the size of this circle, which helps in selecting the area more precisely. 2. Move the brush to the desired area, then hold down the left mouse button to brush over it to select (highlighted in green) the outer layer data of the scanned model that has been brushed; before brushing, you can click to use the Transmission tool to select both the inner and outer sides of the jaw simultaneously. 3. After selecting the area, you can use the region size adjustment slider to repeatedly enlarge or reduce the selected area's size.
 Eraser	<ol style="list-style-type: none"> 1. Click the button, and the button will switch to  ; in the operation interface, the dark gray circle at the cursor represents the coverage area of the eraser, and you can use the slider to enlarge or reduce the size of this circle, which helps in deselecting the area more precisely. 2. Move the eraser to the desired area, then hold down the left mouse button to brush over it to deselect the outer layer data of the scanned model that has been brushed; before brushing, you can click to use the Transmission tool to deselecting both the inner and outer sides of the jaw simultaneously. 3. After deselecting the area, you can use the region size adjustment slider to repeatedly enlarge or reduce the size of the deselected area.
 Transmission	<ol style="list-style-type: none"> 1. Click the button, and the button will switch to . 2. When using the Brush or Eraser, you can simultaneously select or deselect the corresponding area on the opposite side of the jaw while brushing, which helps in quickly applying or removing the brush strokes in the desired area.
 Automatic filling	<p>Click the button to fill all connected regions.</p>
 Reset	<p>Click the button to reset the model to its initial coordinates.</p>
 Undo	<p>Click the button to undo the previous action; you can click continuously to undo multiple actions in succession.</p>
 Redo	<p>Click the button to redo the previous action; you can click continuously to redo multiple actions in succession.</p>

Base Adjustment

On the interface of **Base adjustment**, drag the base manually to set the height of models by entering values.




Note

When adjusting the base, you can hold down the right mouse button to rotate models or rotate the wheel to zoom in/out, which helps in selecting the area more precisely.

Precise adjustment: Enter values in the input box (range: 2mm ~ 100mm) to accurately adjust the height of the base, and preview in real time.

Note

If the input value is greater than 100, it will be automatically set to 100; if the input value is less than 2, it will prompt "not valid".

Manual dragging: Move the cursor to any position within the base (blue area) to be generated; when the cursor turns to , hold down the left mouse button to drag the base up or down to adjust the height of models.

Attachment Adjustment

On the interface of **Attachment adjustment**, you can add the articulator, text or other attachments to models.


Articulator Adjustment

Support multiple articulator options (with the default articulator selection being the type set during model import): large/medium/small articulator, articulator for quadrant models, etc; the articulator is used to securely hold the upper and lower jaw models, facilitating the examination of the inherent positional relationship between the patient's upper and lower jaw joints.


Note

- Articulator options can not be applied when only one model is imported.
- When adjusting the attachment, you can hold down the right mouse button to rotate models or rotate the wheel to zoom in/out, which helps in selecting the area more precisely.

Move the articulator


1. Click the articulator, a circular transparent adjustment button will appear on it, surrounded by a white circular track.
2. Move the cursor to the transparent area, and when the cursor turns to , hold down the left mouse button to drag the articulator.

Rotate the articulator







1. Click the articulator, a circular transparent adjustment button will appear on it, surrounded by a white circular track.
2. Move the cursor to the white track, and when the cursor turns to , hold down the left mouse button to rotate the articulator along the track.

Text Adjustment

Support adding text labels, such as the owner's name or numerical identifier.

 **Note**

When adjusting the text, you can hold down the right mouse button to rotate models or rotate the wheel to zoom in/out, which helps in selecting the area more precisely.

1. In the **Please enter text** area, enter the text that needs to be engraved on the model; in addition, the **Style** ( **Concave** or  **Convex**), **Height** (0.1mm ~ 2mm) and **Size**.
2. Click  **Add text** button, move the cursor to any position within the area where the text can be added, you can preview the text to be printed.
3. Click on the desired place on the model to add the text, and the corresponding column for text edition will appear under the text box:
 - click  to edit the text;
 - click  to add mirror effect (engrave mirrored text on the other jaw);
 - click  to delete the text.

Text adjustment




Style

Height 1.0mm

Size 3

Note

- The **Mirror** button  will be grayed out when only one model is imported.
- You can hold down the left mouse button to drag the text, so as to adjust its position.
- When you click on the text on the model and it appears in blue color, it means the text is in editing mode, and you can modify the style, height, and size of the text.

Printing View





On the interface of **Printing view**, preview the printing effect of models and set printing parameters.




Note


- If AccuDesign is called from other scanning software (such as IntraoralScan), setting parameters of printing is not displayed.
- When previewing models, you can hold down the right mouse button to rotate models or rotate the wheel to zoom in/out, which helps in selecting the area more precisely.

① Model list


- Click  button on the left side of the model, the corresponding model (shown by default) will be hidden, and the button will switch to .
- Drag the transparency slider to adjust the transparency level of models.
- Click  button to collapse the Model List (expanded by default), and the button will switch to .

② Print settings




Click  button to set **Printer Series, Printer Serial No., Material Brand, Material Name, Layer Thickness** and **One-Click Printing Type**, and click **Save** to save the settings.

Click  **Go to print** button to start up the AccuWare software and the models to be printed will be imported automatically.

Note

- Please install the AccuWare software in the first place.
- If you have not clicked  to adjust printing parameters, a **Print Settings** window will pop up, and you can click **Save** to start up AccuWare.

Export the Model

After entering any step for designing models, in the top right corner you can click  **File** >  **Export project**, or click  **Export model** to save models locally.

Note

If AccuDesign is called from other scanning software (such as IntraoralScan), the designed models will be displayed on the **Model List** when exiting AccuDesign and return back to other software.

On the interface of **Printing view**, you can click  **Go to print** in the bottom right corner to start up AccuWare to print models.

MetronTrack

What is MetronTrack?

MetronTrack is a tool that can compare and analyze the difference of a patient's intraoral data at different times. It can automatically identify and mark tooth positions. The measurement data through the measurement analysis system can serve as a reference for dental diagnosis. After loading multiple sets of model data, it can monitor the changes in tooth positions at different times and intraoral data about tooth wear for orthodontic treatment.

Characteristic

Support Identifying Teeth at Different Times.

Identify teeth at different times, including deciduous dentition, mixed dentition, and permanent dentition. In the mixed dentition phase, Moyers Prediction is accessible. In the permanent dentition phase, Crowding and Bolton Ratio Measurements are accessible.

→ [Reference](#)

Support Moyers Prediction, Crowding Measurement, Bolton Ratio Measurement, Overbite and Overjet Measurement, Molar Relationship Measurement and Curve of Spee.

Measure the degree of crowding at different times and the proportion of the size of the upper teeth to the lower teeth. Based on the measured data, it helps to achieve proper alignment.

→ [Reference](#)

Support Model Comparison.

Monitor tooth position changes and intraoral data about tooth wear for orthodontic treatment.

→ [Reference](#)

Support Measurement Report and Uploading.

Support saving and uploading reports of different measurement data and provide a centralized display for user's reference.


→ [Reference](#)

Functions

Segment



On the interface of segment, models of upper jaw and lower jaw are displayed separately. The software automatically recognizes each tooth with different colors and displays teeth numbers. Users can preview the upper jaw and lower jaw from different perspectives, check the teeth numbers, edit teeth areas, remove teeth and add teeth.

Interface



① Timeline

Via the timeline, users can switch between one set of data and two sets of data or choose the data of mixed dentition or permanent dentition.

- **Mark Teeth** only supports one set of data. Move the  to switch.
- **Measurement Analysis** supports one or two sets of data. Click the node  to select model.
- **Model Comparison** only supports two sets of data.

Note

After entering the functions in the **Mark Teeth** interface, the timeline can't be switched.

② Setting

Report

Save and upload [Report](#) of different measurement data.

Help

Click to open the guide window.

Setting








Select language and choose [Measurement Template](#).

③ Function Bar



Modify Area

Click  to enter the interface of tooth selection. Select a certain tooth to modify its area.


Icon	Name	Description
	Brush	Click Brush to select the tooth area. Drag the slider to adjust the thickness of the brush.
	Erase	Click Erase to delete the selected area. Drag the slider to adjust the thickness of the eraser.
	Undo	Undo the last operation.
	Redo	Redo the last operation.
	Draft	Click Draft to save the current operation and return to the tooth selection interface to select other teeth.
	Cancel	Cancel all operations and exit Modify Area.
	Confirm	Save all operations and return to the segment interface.

Caution

When tooth regions are wholly erased, a tip of "It is not allowed to remove all tooth regions!" pops up.







Change Tooth Number

1. Click  to enter the interface of changing tooth number.
2. Click the tooth to re-number it.
3. Select the tooth number according to the actual situation of the patient.

Note

- The outer tooth number is permanent tooth (permanent dentition for adults), and the inner tooth number is deciduous tooth (period of mixed dentition for children).
- When the mixed dentition is selected, only the Moyers Prediction function is supported; when the permanent dentition is selected, the Crowding Measurement, Bolton Ratio Measurement and model comparison function are supported.




Icons	Description
	Click to switch to primary teeth.
	Click to switch to permanent teeth.
	Click to cancel all operations and exit.
	Click to confirm and save all operations.

Caution

Do not repeat the same tooth number; the red number means tooth position is repeated.


Remove teeth

Steps


1. Click  to enter the interface of removing teeth.
2. Select the tooth which should be deleted and the tooth number is displayed in red.
3. Click  to delete it.
4. If more than one tooth needs to be deleted, repeat steps 2-3.
5. Click  to confirm and exit the operation interface.

Add Teeth

Steps




1. Click  to enter the interface of adding teeth.

2. Double-click the position where the tooth should be added.
3. Select the tooth number in the pop-up window.
4. Click ✓ to confirm.

 **Caution**

Repeat steps 1-4 to add more teeth. And the tooth number should be different.

④ Visible Control

Icon	Name	Description
	Number	Enabled by default. When enabled, the teeth numbers are displayed.
	Texture	Enabled by default. When enabled, the model is colored.
	Area	Enabled by default. When enabled, teeth are differentiated from each other with different colors

⑤ Preview

Support previewing the scanned data. To gain a comprehensive view to the model, please use the following shortcuts for switching perspectives.

- Hold the Right Mouse Button: Rotate the model.
- Hold the Left and Right Mouse Button: Move the model.
- Scroll the Mouse Wheel: Zoom in or zoom out.

Measurement Workflow



Select **measurement items**

Manually adjust dental arch curve/tooth width





Obtain dental arch/tooth width measurement data



Save the [measurement report](#)

Measurement Guidance

Click  to enter the interface of measurement analysis.




Measurement Items

Moyers Forecast



For patients in the mixed dentition stage, the software will predict the crowding of the permanent teeth that are yet to erupt based on the existing mandibular permanent teeth.

 **Caution**

Moyers forecast can only be used in mixed dentition state teeth. When only a single upper jaw is imported, Moyers forecast is not accessible.

The dental arches of the upper jaw and lower jaw will be displayed.

Adjust the dental arch manually:

1. Click the dental arch when it turns into blue and a window pops up.
2. Move the control point on the dental arch to adjust the arch.

Adjust the tooth width manually:

1. Click a tooth and a window showing its labio-buccal direction and mesial direction pops up.
2. Rotate the tooth to adjust the labio-buccal direction and distal direction. The tooth width is adjusted as well.

Dental Crowding Measurement

The dental crowding measurement is used to measure the crowding degree of the teeth. Crowding teeth usually happen when the teeth are not regularly arranged. By analyzing the teeth model, the sum of teeth width and the circumference of the dental arch are measured and the difference between the two values represents the degree of crowding. Crowding measurement is an important indicator in determining whether tooth extraction is necessary.

Moyers forecast can be used for crowding measurement in mixed dental stage. While in permanent dentition stage, moyers forecast is not accurate. Thus, **dental crowding measurement** is used to measure the crowding degree in the permanent dentition state.

The dental arches of the upper jaw and lower jaw will be displayed.

Adjust the dental arch manually:

1. Click the dental arch when it turns into blue and a window pops up.
2. Move the control point on the dental arch to adjust the arch.

Adjust the tooth width manually:

1. Click a tooth and a window pops up.
2. Rotate the tooth to adjust the labio-buccal direction and distal direction. The tooth width is adjusted as well.

Bolton Ratio Measurement

The Bolton Index refers to the proportional relationship between the widths of the upper and lower teeth. It helps preliminarily determine whether the upper and lower teeth are proportionally aligned or if there is an abnormality in their width relationship, serving as a reference for diagnosing dental malocclusions. Besides, the Bolton Ratio Measurement enables to analyze the size of the upper and lower teeth to achieve effective adjustment of teeth alignment. Therefore, this measurement is also valuable as a reference in orthodontic treatment.


Caution

The result of the Bolton ratio can be influenced by factors such as the malocclusion type, gender, ethnicity, extraction patterns, measurement errors, tooth thickness, and torque, which may introduce errors in the Bolton Ratio Measurement. Therefore, it is important to consider these factors and to use this measurement as a reference alongside other clinical assessments for an accurate treatment plan.

Manually adjust tooth width

- 1 Click a specific tooth to display the adjustment panel in lower right corner to adjust tooth width.
- 2 Rotate the tooth to align its mesiodistal and buccolingual directions with the coordinate axes. The measurement value of the tooth width will be updated in real-time.

Overbite and Overjet Measurement


Click  to enter Overbite and Overjet Measurement.

Steps

1. Click a specific point on the model to display the adjustment panel in lower right corner.
2. Manually adjust the point in the panel, bringing a corresponding change in **Overbite** and **Overjet** data.
Rotate the tooth to align its mesiodistal and buccolingual directions with the coordinate axes.
3. For example, click **Middle point of Incisal edge in upper jaw** and manually adjust its position in the panel in lower right corner.
 - The position changes between central incisor in upper jaw and central incisor crown in upper jaw reveal the overbite data.
 - Move the **Middle point of Incisal edge in upper jaw** to adjust the distance between this point and **Middle point of Incisal edge in lower jaw**. The overjet data on measurement table will be updated.

Position	Overbite	Distance(mm)	Overjet
No contact	Open jaw	≤ 0	Crossbite
Within incisal edge 1/3	Normal	0~3(including 3)	Normal
Incisal edge 1/3~1/2	I ° deep overbite	3~5(including 5)	I ° deep overjet
Incisal edge 1/2~2/3	II ° deep overbite	5~7(including 7)	II ° deep overjet
Above incisal edge 2/3	III ° deep overbite	>7	III ° deep overjet


Molar Relationship Measurement

Click  to enter **Molar Relationship Measurement**. The points are displayed differently in primary dentition, permanent dentition and mixed dentition.

Steps

1. Click a specific point on the model to display the adjustment panel in lower right corner.
2. Manually adjust the point in the panel, bringing a position changes between buccal groove and buccal cusp. Rotate the tooth to align its mesiodistal and buccolingual directions with the coordinate axes.
3. Take the permanent dentition as example.
 - Click the point **Mesiobuccal cusp in upper jaw** on the model.
 - Drag this point at desired position, changing the distance between this point and other three points in lower jaw, including Distobuccal cusp, Mesiobuccal cusp and Buccal groove.
 - Accordingly, the measurement value will be updated, defining Molar relationship.

Position	Molar relationship
The Mesiobuccal cusp in upper jaw is located between Distobuccal cusp and Mesiobuccal cusp in lower jaw.	Neutral relationship
The Mesiobuccal cusp in upper jaw is located in front of Mesiobuccal cusp in lower jaw.	Mesial relationship
The Mesiobuccal cusp in upper jaw is located behind Distobuccal cusp and Mesiobuccal cusp in lower jaw.	Distal relationship


 **Note**

When one or two teeth are missing in each pair of upper and lower jaws, the points are not displayed.

Curve of Spee

Click  to enter Curve of Spee measurement.

Click a point from left or right view to manually adjust Curve of Spee, bringing a change in measurement value updated on the table.


 **Caution**


If too many teeth are missing, the Curve of Spee measurement is not accessible.

Measurement Template

Two templates are provided in MetronTrack for your reference.


Steps

1. Click  in the right top corner.
2. In the pop-up setting window, select one template from the dropdown menu: default 1 or default 2.













3. Click  to check, and then click **Select this template**.

Measurement Report

Measurement results and images are visible in the report for use's reference. For different measurement options, reports vary from Moyers Prediction, Crowding Measurement, Bolton Ratio Measurement, Overbite and Overjet Measurement, Molar Relationship Measurement, Curve of Spee and Analysis of tooth movement. Click  in the upper right of the interface, and the measurement report will pop up.


Note


- Users can select the scan data from different dates on the timeline, which will be displayed as model data in the report.
- Users can choose different measurement methods for the report output, as shown in the table below:

Icon	Name	Note
	Edit the Report	Edit patient information.
	Moyers Prediction	Display upper/lower jaw model views and Moyers Prediction data.
	Crowding Measurement	Display upper/lower jaw model views and Crowding Measurement data.
	Bolton Ratio Measurement	Display upper/lower jaw model views and Bolton Ratio Measurement data.
	Overbite and Overjet Measurement	Display upper/lower jaw model views, right view/lower jaw, and Overbite and Overjet data.
	Molar Relationship Measurement	Display upper/lower jaw model views, right view/lower jaw, and Molar Relationship Measurement data.
	Curve of Spee	Display upper/lower jaw model views and Curve of Spee data.
	Analysis of Tooth Movement	Display upper/lower jaw model views and Tooth Movement data.
	Save PDF	Download in PDF format. Export the report and name it as 'Name_Report date'. The exported report will contain the scan data for the date on the current timeline.
	Upload	Upload the report to cloud.

Model comparison

Model Comparison is a tool that can compare and analyze the difference of a patient's intraoral data at different times. It can monitor multiple sets of tooth position changes and intraoral tooth wear data for orthodontic treatment.

Click  to enter model comparison interface.



① Model Alignment



The models from different treatment periods are aligned according to the feature points, both automatically and manually.

Click  and  in the upper right corner of the interface to select jaw-plane view of upper/lower jaw.


Note

In the model comparison interface, the automatically aligned total-jaw model will be displayed by default.

Manual Alignment

1. Left-click to add new points on model 1 and model 2. Press and hold the right mouse to move them. After the added points are set, click  to manual alignment. The overlay effect of two models will be displayed.
2. Click  to reset if needed, leaving only the default three points remained.

② Analysis of teeth movement

1. Click  to enter Analysis of teeth movement. Default is the overlay effect of two models.
2. Click one tooth to highlight the corresponding data on the mobility scale table.

- When you select a specific measurement value on the mobility scale table, the model will display the corresponding tooth and its movement reference line.

The screenshot shows the SHINING3D MetronTrack software interface. At the top, there are three tabs: 'Mark teeth', 'Measurement analysis', and 'Model comparison'. Below the tabs is a 3D model of a dental arch with two models, Model1 and Model2, overlaid. A 'Mobility scale table' is displayed at the bottom, showing measurement data for various teeth (18-28) across different categories: Intrusion/Extrusion, Translation (Buccal/Labial/Lingual), Translation (Mesial/Distal), Rotation (Mesial/Distal), Angulation (Mesial/Distal), and Inclination (Buccal/Labial/Lingual). The table has a grid with columns for tooth numbers and rows for measurement types. The value for tooth 12 in the 'Translation: Mesial (M) / Distal (D) (mm)' row is highlighted in blue.

Measurement data	18	17	16	15	14	13	12	11	21	22	23	24	25	26	27	28
Intrusion (I) / Extrusion (E) (mm)	-	0.0E	1.0I	0.9I	0.3E	0.5E	1.8E	1.4E	1.5E	0.3E	3.0I	1.1I	1.9E	0.7I	0.9I	-
Translation: Buccal (B) / Labial (La) / Ling...	-	1.7B	1.9B	2.1B	1.2B	0.5U	4.2U	0.1U	0.3U	0.1U	2.6La	2.7B	2.4B	2.3B	2.1B	-
Translation: Mesial (M) / Distal (D) (mm)	-	0.7D	0.1D	1.0M	1.3M	1.5M	1.0M	0.7M	1.2D	1.6D	0.7M	0.7D	1.0D	1.9D	1.7D	-
Rotation: Mesial (M) / Distal (D)	-	8.2°D	9.2°M	10.2°M	32.0°M	26.4°M	22.4°D	4.2°D	0.3°M	22.2°M	39.6°M	30.7°M	13.7°D	1.2°M	9.3°D	-
Angulation: Mesial (M) / Distal (D)	-	7.9°M	0.9°M	9.3°D	19.3°D	13.1°M	38.1°D	4.3°D	4.6°M	2.2°D	11.1°M	2.1°D	36.3°M	1.4°D	1.5°D	-
Inclination: Buccal (B) / Labial (La) / Ling...	-	1.5°LI	0.6°LI	1.1°B	0.9°B	4.0°La	4.5°LI	1.4°La	6.8°LI	9.1°LI	15.0°La	0.6°LI	36.5°LI	8.0°B	8.9°B	-

Note

Movement value of new teeth and dental implants cannot be measured.

③ Ribbon Chart


- Click  to view the ribbon chart for teeth comparison. Drag the slider in the right bottom to set the range.

Note

- The movement variables of the teeth in the new data are compared with the scanned data from earlier times.
- Positive values are set to red, which indicate a shift toward cheek; negative values are set to blue, which indicate a shift toward tongue.
- The function of ribbon chart and sectional view cannot be enabled at the same time.


- Move the cursor to one single tooth and its numerical value will be shown in detail.
- Support switching the comparison area: jaw or tooth.


- The photoscope can be dragged to model for viewing. Left-click the photoscope to display the magnification effect of Model 2. Otherwise, the photoscope will display the magnification effect of Model 1.

 **Note**


Colors in the photoscope indicates the differences between Model 1 and Model 2 in this target area. The darker the color, the more different they are.


④ **Sectional View**

- Click  for the sectional view of teeth model.
- Move and rotate the section to the proper position and angle you want.



- Comparison between new and previous scanned data is shown in the bottom-right corner. Users can click two different points to measure the distance.



 **Note**

- The darker lines are the newer scanned data, and the lighter lines are scanned data from earlier times.
- Compare the two lines to detect the changes of teeth movement.

4. Support switching the comparison area: jaw or tooth.


⑤ **Display Effect**

 **Texture**

When this function is enabled, the model is colored.




Show Texture




Hide Texture

 **Overlay**

Enabled by default. When enabled, one model will be overlaid with the other model.



Show Overlay



Hide Overlay

Temp crown

Before making a temp crown, designs and adjustments should be made according to form of the tooth that needs restoration. Temp crown design software is easy-to-use and convenient, and it can build temp crown models in a relatively short time.

With the prevailing popularity of 3D printing technology, temp crown made by digital 3D printing gradually replaces the traditional handmade temp crown by dentists.



Note


- Temp crown is a full crown that temporarily sticks in the patient's oral cavity after the tooth preparation and before the crown restoration.
- The full crown is a restoration that can cover the whole tooth. The full crown can repair the form of the broken tooth and restore its function.



Caution


Only when the order is **Restoration**, can the temp crown design be proceeded.

Click  >**Tempcrown** in Pre Design interface.




Pre Design


1. Select the tooth that needs the temp crown design on the model list.
2. Double click the tooth that needs repair, and extract margin line for teeth automatically.




3. Adjust the margin line by dragging the point.

 **Note**

Press and hold down the left button of the computer mouse, and choose and drag the point (When the cursor is placed on the point, its color will turn into green from red.), and then adjust the margin line.


4. Click  , to adjust the direction of the path of insertion. There will be an arrow on the tooth that needs restoration. Double-click the arrow to adjust the direction.









5. Adjust the position to check the position of the temp crown. The ribbon chart will show on the crown in real time.

Crown Design

Click **Crown Design** and there will be a crown model on the tooth that needs restoration. The crown model is light blue, and the thickness of the crown is dark blue.



Icon	Name	Description
	Sculpture	The sculpting tools enables operations including extrusion, intrusion, and smooth. The data can be denoised to make the model surface smoother and improve the data quality. Add or reduce model data of the model surface according to the actual needs.
	Parameter adjustment	Enables the user to modify the parameter of Inner wall reserve gap and Inner wall edge width .
	Sectional View	It's a two-dimensional measurement function. A screenshot plane will be automatically created after entering Sectional View. The section can be automatically positioned on the crown by double-clicking on the crown model. The two-dimensional window in the lower right corner shows the corresponding section line. You can drag, zoom in and out, and reset the position. Jaw, crown, and minimum thickness are presented by different color lines. The section line is updated in real time as the section moves, and two points can be selected in the section line for measurement.
	Interference detection	The virtual thickness model is displayed in real time on the crown model and cannot be modified. Interference detection shows the contact ribbon chart between the tooth and the adjacent teeth and the opposite jaw.
	Texture	The texture function is turned on to displays the colored model.
	Margin line	With the margin line function on, you can adjust model to its internal side to view the margin line and check if the margin line and the crown are in line.


Once you have completed the crown design, you can click **To Print** to send the project file containing the crown model to the printer for printing.

CrelBT

Introduction


CrelBT (Indirect Bonding Trays) is a design software for digital bracket tray. It can virtually position the bracket on the surface of the tooth. The custom guide correlates the position of the brackets with that of the tooth surface, facilitating the bonding of the brackets onto the desired teeth positions.

After collecting intraoral data from the patient, the doctor can use the IBT design software to arrange teeth, position and place brackets, and generate bracket guides for the patient's dental model. If the model needs to be printed, it can be printed through the AccuWare software developed by Shining 3D as well.

 **Caution**

- Only when the order type is Orthodontics, can the design be performed on IBT.
- To print the guide, please install AccuWare first.

Workflow



Guide


Place brackets according to the FA points. Manually adjusting the positions of brackets are supported. The database containing various brands of brackets allows users to customize their brackets and guides.

Generate guides

Adjust the guides and preview the position.

Functions

Settings

Click  to open the setting window.

Item	Description
Language	You can set the software interface language: Chinese (Simplified) and English. The default language is the language selected during software installation.
Dental Notation	FDI World Dental Federation(ISO) notation and Universal numbering system. The default is FDI World Dental Federation(ISO) notation.
Guide parameter	Click Create a New Template to set default parameters.
Auto pop-up guidance	Enabled by default. Video guidance is available.

Create a New Template

You can change the value of the gap between the guide and the model, diameter of both ends of the connecting rod, diameter of middle part of the connecting rod and gap between the guide and the brackets.

Create a New Template



Template Name

The gap between guide and model ⓘ

mm

Diameter of both ends of the connecting rod ⓘ

mm

Diameter of middle part of the connecting rod ⓘ

mm

Gap between the guide plate and the brackets ⓘ

Overall adjustment of bracket gap in the upper jaw

Tooth position	Gap Value (mm)	Tooth position	Gap Value (mm)
11	-0.05	21	-0.05
12	-0.05	22	-0.05
13	-0.05	23	-0.05
14	-0.05	24	-0.05
15	-0.05	25	-0.05

Overall adjustment of bracket gap in the lower jaw





Tooth position	Gap Value (mm)	Tooth position	Gap Value (mm)
----------------	----------------	----------------	----------------



Save



Cancel

Name	Illustration
The gap between the guide and the model	
Diameter of both ends of the connecting rod	
Diameter of middle part of the connecting rod	
The gap between the guide and the brackets	

Steps


1. Click **Create a New Template** to set the template's parameters.
2. Enter the **Template Name**. The name must be unique and cannot be left blank, with a maximum of 20 characters.
3. Adjust the **The gap between guide and model**, **Diameter of both ends of the connecting rod**, and **Diameter of middle part of the connecting rod**.
4. Modify the **Gap between the guide plate and the brackets**.
You can adjust the overall gap between the guide and the upper and lower brackets and buccal tubes, or

change the gap value for individual teeth positions. The default value is -0.05 mm. Values are rounded to two decimal places and are adjustable from -0.10 mm to +0.10 mm.

5. Click **Save** to return to the settings interface, where it will automatically switch to the newly created template.

Caution


When switching templates, the parameters of the new template will take effect the next time a guide is generated. Existing guides will not be affected.

After changing the items, click  and a confirm dialog is popped up. Click **Yes** to apply the changes.

Segmentation

On the interface of segmentation, models of upper jaw and lower jaw are displayed separately. The software automatically recognizes each tooth with different colors and displays teeth numbers. Users can preview the upper jaw and lower jaw from different perspectives, check the teeth numbers, edit teeth areas, remove teeth and add teeth.

Interface





Functions

Name	Description
------	-------------





Enlarge


Click  to enlarge the model of upper jaw or lower jaw to the full screen. Click  to restore.





Modify Area


Click  to enter the interface of tooth selection. Select a certain tooth to modify its area.


 : Click Brush to select the tooth area. Drag the slider to adjust the thickness of the brush.


 : Click Erase to delete the selected area. Drag the slider to adjust the thickness of the eraser.

 : Undo the last operation.

 : Redo the last operation.

 : Click Draft to save the current operation and return to the tooth selection interface to select other teeth.


 : Cancel all operations and exit Modify Area.

 : Save all operations and return to the segment interface.





Caution


When tooth regions are wholly erased, a tip of "It is not allowed to remove all tooth regions!" is popped up.



Change Tooth Number

Click  to enter the interface of changing tooth number. Click the tooth to re-number it.

 : Click to cancel all operation and exit.

 : Click to confirm and save all operations.



Caution

The red number means tooth position is repeated. When confirming, a tip is popped up.







47 Tooth position is repeated, please check!



Remove Teeth




Steps

1. Click  to enter the interface of removing teeth.
2. Select the tooth which should be deleted and the tooth number is displayed in red.
3. Click  to delete it.
If more than one tooth needs to be deleted, repeat steps 1-3.
4. Click  to confirm. Click  to cancel all operations and exit.



Add Teeth

Steps




1. Click  to enter the interface of adding teeth.
2. Double-click the position where the tooth should be added.
3. Select the tooth number in the pop-up window.
4. Click  to confirm. Click  to cancel all operations and exit.



Caution

Repeat steps 1-4 to add more teeth. And the tooth number should be different.

Options

Name	Description
 Number	Enabled by default. When enabled, the teeth numbers are displayed.
 Texture	Enabled by default. When enabled, the model is colored.
 Area	Enabled by default. When enabled, teeth are differentiated from each other with different colors.

Arrangement


The tooth arrangement ensures precise alignment of teeth and prevent collisions and excessive gaps between them, which can assist users in [marking FA points](#).

In the **Tooth Arrangement** interface, the software automatically arranges the teeth based on the tooth segmentation results, and users can manually edit the automatic arrangement results.

Interface





Model:



- Click  to show/hide the upper and lower jaws to more intuitively observe the arrangement of the teeth.
- Drag the slider to adjust the opacity of the upper and lower jaws.



Operations


Item	Description
------	-------------

 <p>Shift</p>	<p>Left-click and hold on a single tooth, then move the cursor to move the tooth.</p> 
------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------

 <p>Rotation</p>	<p>Left-click and hold on a single tooth, then move the cursor to rotate the tooth to any angle.</p> 
---------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

 <p>Torsion</p>	<p>Left-click and hold on a single tooth, then move the cursor to twist the tooth left or right.</p> 
----------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------


 <p>Interproximal contact</p>	<p>Adjust the interproximal distance between teeth. When the distance is zero, only a small dot is displayed. Otherwise both a small dot and a label value are displayed. Click the dot or label to adjust the distance from -1 to 2.</p> 
------------------------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

- 
- : Contact
 - : Collision
 - : Gap



Occlusion

Display the occlusion status by the color bar. It can quickly and intuitively help determine which teeth have occlusion problems. Click the ▲ or ▼ beside the bar to change occlusion parameters.



- **Green:** The gap.
- **Red:** The touching area
- **Blue:** The bite-through area



Reset



Undo all operations and reset the model to the automatic tooth arrangement state.



Undo/Redo

Undo the previous step / Revert the undone operation.


Options

Item	Description
 Texture	Enabled by default. When enabled, the model is colored.
 Grid	When enabled, the model background shows grid for reference when adjusting teeth vertically or horizontally.



FA points

On the interface of FA, the software will automatically generate the FA points and the user can preview the three views of the FA points on the 3D model as well as the reference lines and adjust them as needed. The teeth in the three views will correspondingly change with the adjustments of the FA points.

Interface



Model:

- Click  to show/hide the upper jaw or lower jaw.
- Drag the slider to adjust the opacity of the upper jaw or lower jaw.
- Click  to show/hide the reference lines.










Operation


Move/Rotate the FA point

- Movement: Click and hold the FA point that needs to be adjusted, and move the cursor to adjust the FA point position.
- Rotation: Click the FA point you want to rotate, hold the green arrows on either side of the FA point and move the cursor to rotate the FA point.




Fine movement/rotation

- Fine Movement: Press and hold the keyboard arrow keys for fine movement ( /  /  / ).
- Fine Rotation: Hold down control + arrow keys at the same time to perform a fine rotation ( +  /  /  / )










Move/Rotate the FACC line




With the FA point of a tooth position selected, move the cursor over the FACC line (horizontal axis) of the tooth position, and when the line is blue, drag the line to rotate the FACC line.



Toolbar

Operation	Description
 Import Panoramic X-ray	<p>Click to import an image (PNG, BMP, JPG, JPEG formats supported). The viewer will appear as a floating window, allowing you to zoom and move the image. On first use, a tutorial will guide you through the operations. After the tutorial closes, hover over  to view it again.</p>  <p>The screenshot shows a floating window titled "Panoramic Radiograph Viewer" with a close button. It displays a panoramic radiograph with a dashed vertical line. On the left, under "Shift", it lists "Mouse Operation" (Hold down the left mouse button and drag) and "Gesture Operation" (Press and hold with one finger and drag). On the right, under "Scale", it lists "Mouse Operation" (Middle Scroll Button) and "Gesture Operation" (Double finger expansion to zoom in, Double finger contraction to zoom out).</p>
 Reset	<p>Rest the FA points of all dental position to the default position. Click the button and a tip is popped up. Click Yes to confirm.</p>  <p>The tip dialog box has a title bar "Tip" with a close button. It contains an exclamation mark icon and the text: "Reset the FA points of all dental position to the default position, Are you sure?". At the bottom, there are two buttons: "Yes" (with a checkmark) and "No" (with an X).</p>
 Undo	<p>Undo the last operation.</p>
 Redo	<p>Redo the last operation.</p>


Options

Item	Description
 Texture	Enabled by default. When enabled, the model is colored.
 Area	Enabled by default. When enabled, teeth are differentiated from each other with different colors.
 Grid	When enabled, the model background shows grid for reference when adjusting teeth vertically or horizontally.




Bracket

You can select and edit brackets in this interface, where bracket information and model information are displayed. You can also choose to show or hide the upper and lower jaw and brackets, or adjust their opacity.

Interface




- **Bracket Information :**
Display the brand and model of the bracket or buccal tube.
- **Model :**

- Click  to show/hide the upper jaw or lower jaw.
- Drag the slider to adjust the opacity of the upper jaw or lower jaw.
- Click  to show/hide the reference lines.
- Click  to show/hide the collision parts in blue between the brackets and the teeth.

Bracket selection

Steps

1. When you first enter the main bracket interface, the bracket selection interface will automatically appear.

Alternatively, you can click  to select bracket brand and model.

X
Select Bracket

! The bracket database has been updated Update Now

Bracket
Clear all selections

Brands

Please select the bracket brand

Type

Please select the bracket type

Add method

Add automatically ?

Buccal Tube
Clear all selections

Brands

Please select buccal tube brand

Type

Please select buccal tube type

Add method

Add automatically ?

Upper tooth
 Lower tooth

Upper tooth	17	16	15	14	13	12	11	21	22	23	24	25	26	27
Angulation														
Torque angle														
Anti-rotation angle														

✓ OK

✕ Cancel

2. Select the brand and model of bracket or buccal canal (or both) as needed, and the software will display the axial inclination, rotational angle and anti-rotation angle of each tooth in the upper and lower positions.

info

You can select the brand and model of brackets and buccal tubes separately, such as brand A brackets paired with brand B buccal tubes. Besides, you can also choose only one of the brackets and buccal tubes.


3. Click **Confirm**, brackets are automatically added to the model.

info

- If the imported teeth model has brackets, which are not included in the database, you have to choose whether to continue to use the current bracket data. If not, you need to choose a new bracket model.
- If there is an update to the bracket database, an update prompt will appear at the top of the Bracket Selection interface.


Bracket editing

The main interface allows the user to preview and edit three views of the brackets on the 3D model.













Move/Rotate the bracket

- To move a bracket: Click and hold the brackets you want to move, then move the cursor to adjust the bracket position.
- To rotate a bracket: Click the brackets you want to rotate, hold the green arrows on either side of the brackets and move the cursor to rotate the brackets.




Fine movement/rotation










- Fine Movement: Press and hold the keyboard arrow keys for fine movement ( /  /  / ).
- Fine Rotation: Hold down control + arrow keys at the same time to perform a fine rotation ( +  /  /  / ).





Delete the bracket

- To delete a bracket: Select a bracket to be deleted and click  to delete the bracket.

Toolbar

Operation	Description
 Import Panoramic X-ray	<p>Click to import an image (PNG, BMP, JPG, JPEG formats supported). The viewer will appear as a floating window, allowing you to zoom and move the image. On first use, a tutorial will guide you through the operations. After the tutorial closes, hover over  to view it again.</p> 
 Replace	<p>Select the bracket and if there is another bracket which can be replaced, click  to confirm the tooth bit.</p>
 Reset	<p>Undo all operations.</p>
 Undo	<p>Undo: Undo the last operation.</p>
 Redo	<p>Redo: Redo the last operation.</p>
 Delete	<p>Click to delete the selected bracket.</p>

Options

Item	Description
 Texture	Enabled by default. When enabled, the model is colored.
 Grid	When enabled, the model background shows grid for reference when adjusting teeth vertically or horizontally.

Guide

You can create and edit guides on this interface.


Interface



List to be generated:

A list includes the selected tooth positions for generating a guide.

Model:

- Click  to show/hide the upper or lower jaw and brackets.
- Drag the slider to adjust the opacity of the upper or lower jaw and brackets.

Guide create



Steps

1. Selecting Brackets:








- **Single Bracket:** Click directly on a bracket to select it.
- **Multiple Brackets:** Click and drag the cursor over multiple brackets to select them all at once.
- **Preset Teeth Regions:** The toolbar provides pre-defined regions of teeth for quick selection. Simply click on the desired region in the toolbar to select all brackets within that area.
- **Deselect Brackets:** Click a selected bracket again to deselect it.

2. Click **+** to add the selected brackets into the **List to be generated**.


3. Click **✓** to confirm the selected brackets and create guides.


4. Select a group of guides in the **Guide List** to preview the guides on the model. The selected guides are displayed in blue. Others are displayed in grey.


Toolbar

Icon	Description
 <p>Multiple teeth selection</p>	Select different regions of teeth.
 <p>Reset</p>	Clear all selected brackets.
 <p>Add to list</p>	Add the selected brackets into the List to be generated .
 <p>Cancel</p>	Cancel all operations and exit.
 <p>Confirm</p>	Save all operations.
 <p>Delete</p>	Delete the selected guides. Click  and a confirmation dialog is popped up. Click Confirm to delete.





Guide edit




1. Select a group of guides in the **Guide List** and click  to enter the interface to edit guide.
2. Click to adjust the position and size of the guides.

3. Click  to confirm the edit.
4. Click **To Print** to print the guides on AccuWare.

Toolbar


Operation	Description
 Undo	Undo: Undo the last operation.
 Redo	Redo: Redo the last operation.
 Cancel	Cancel all operations and exit.
 Confirm	Save all operations.

Options

Item	Description
 Texture	Enabled by default. When enabled, the model is colored.

Export model with brackets

Steps

1. Click  in the top right corner.
2. Select the **Export model with brackets** option from the menu.
3. The file browser will launch, allowing you to choose a folder to save the model file.
4. If an export success message appears, click the **View** button to open the folder's directory.

Note

If an export failure message appears, you can click the **Retry** button to attempt the export again.

Report




After creating the guide, you can click  in the top right corner to generate a report.

This report includes the patient's basic information, attending dentist's name, guide details and different types of model view.

Caution

This report is for reference only. Please use it in conjunction with clinical conditions.

Functions

Name	Description
 Edit report	Edit the basic information of the report.
 Export PDF report	Save the report as a PDF file to a local path.
 Upload report	Click to upload the report to the cloud. After uploading, you can scan the pop-up QR code to view and share the report, save the QR code image locally, or copy the link.

CreSplint

Splint Introduction

CreSplint is an intelligent design software for splints developed by Shining 3D. After collecting intraoral data from the patient, the doctor can invoke CreSplint design software.

With this software, when importing a teeth model, the software will automatically design an expected splint quickly and conveniently according to the teeth shape by changing the path insertion, the edge line and relevant parameters.


The doctor can manually edit the shape of molar pad or retainer as well.

After designing, the doctor can save the file to a local disk or export it into the printer and print the splint.

CreSplint supports **Molar pad** and **Retainer**.

The **Molar pad** is specifically designed to prevent tooth wear for people who grind their teeth when sleeping.

The **Retainer** is an appliance used after orthodontic treatment to protect the results and restrain the teeth in the expected positions, ensuring both aesthetic appeal and functionality of the teeth.

 **Caution**

- Only when the order type is **Orthodontics**, can the design be performed on Splint.
- To print the splint, you need to install AccuWare first.

Characteristic



Edit Path Insertion

Change the direction in which the molar pad or the retainer is positioned to change the shape of the molar pad/retainer.

→ [Edit path insertion](#)



Edit Edge Line

Edit the edge line to change the range of the molar pad/retainer.

→ [Edit edge line](#)



Edit Parameters

Edit **Undercut Ratio**, **Thickness** and **Interval between the molar pad/retainer and teeth**.

→ [Edit parameters](#)



Visualization Tool

View the molar pad/retainer clearly via **Sectional View** and **Texture**.

→ [Visualization tool](#)

Settings

Click to open the setting window.

Setting



Select Language

English



Dental Notation

FDI World Dental Federation notation




Guide



Caution

If the user opens CreSplint from IntraoralScan, the settings of **Language** and **Dental Notation** in CreIBT follow those in IntraoralScan and are not displayed.

Item	Description
Language	Users can set the software interface language: Chinese (Simplified) and English. The default language is the language selected during software installation.
Dental Notation	Select FDI World Dental Federation notation or Universal numbering system . The default is FDI World Dental Federation notation.
Guide	Show the video guidance. Checked by default.

After changing the items, click  to save.

Data edit

After importing the model, the software will automatically create molar pad or the retainer.

In the main interface, users can edit the path insertion, edge line and parameters.


Note




The operations of editing molar pad and retainer are similar. This manual takes editing molar pad as an example.



Model list

By moving the sliders on the left upper corner of the interface, the user can respectively change the opacity of the models and the pads.



Click  to hide the model or the molar pad/retainer.


Click  /  /  to check the model from different perspectives.


Path insertion



The path insertion refers to the direction in which the molar pad or the retainer is positioned. Different path insertions have different undercut ranges, affecting the shape of the molar pad or the retainer.

Steps

1. Click  to enter the window to edit the path insertion.
2. Click the blue arrow pointing to the upper/lower jaw and the software will automatically adjust the path insertion.
3. Click  to confirm the path insertion after changing.

 **Note**

- Click  to exit editing path insertion. If there are already some changes, then a tip is popped up.


 **Tip** 





Confirm exit?The current operation is not applied after exiting

- Click **Exit** to cancel the changes.

Edge line


The edge line wraps the range of the molar pad or the retainer. Editing the edge line changes the range.

Click  to enter the window to edit the edge line. The software will automatically generate edge line.


Object	Operation
Edit edge line	Select a control point and move the cursor to edit the edge line.
Delete a control point	Double-click right mouse button to delete the certain control point.
Add a control point	Press the line and move the cursor to add a new control point.
Restore the edge line	If it is a straight line, click  to restore to the initial state. If it is a curve, click  .
Confirm	Click  to confirm the changes.
Cancel	<ul style="list-style-type: none"> After changing, click  and a tip is popped up. Click Exit to cancel all operations and return back to the main interface. The tip isn't popped up if there is no change.

Note

- When **Molar pad** is selected, the edge line is a straight line by default. While when **Retainer** is selected, the edge line is a curve by default.
- When switching the edge line shape, a tip is popped up. Click **Change** to confirm.




Parameter

Click  to enter the window of editing parameters.


In this window, users can edit **Undercut Ratio**, **Thickness** and **Interval** between the molar pad/retainer and teeth.

- **Undercut Ratio** refers to the ratio of filling undercut in the current path insertion. 100% indicates that all undercut is filled. 0% indicates that no undercut is filled.
- **Thickness** refers to the thickness of the molar pad or the retainer.
- **Interval** refers to the interval between the molar pad/retainer and the teeth.







Visualization tool

Sectional view


Click  to enter the interface of sectional view. A screenshot plane will be automatically created on the model.


The 2D window in the lower right corner shows the corresponding intersecting parts between the plane and the model.



Object	Description
Move the plane	Press the ball on the middle of the plane and move the cursor.
Rotate the plane	Press the other part except the ball of the plane and move the cursor.
Zoom in/out the 2D window	Click  and  .
Reset the position of the plane on the model	Click  .

Texture

Click  to display the colors.



AccuWare

Click  in Pre Design interface.

After importing the data into AccuWare, the user can slice the model and send it to the printer. For more, see [AccuWare User Manual](#).

Tools

Introduction to Dental Cloud

Dental Cloud is an integrated platform which provides service support and service invoke for users. Users in different identities can customize and choose the services they need and search, acquire, transmit and share data.

Dental Cloud is a platform for connecting clinical and technical workers. On this platform, users can establish connections between labs, doctors and hospitals/clinics, manage orders, manage organizations and maintain account information.

Visit the website at: [Dental Cloud](#) [🔗].

Use ScanBinder

1. On the start menu, select DentalLauncher > ScanBinder to open the scan binder.

 **Note**


Support binding of ExoCAD and Dentalwings design software.

2. Select the path of DentalDB (supporting Chairside 2.1, ExoCAD 2018, 2019 and 3.0).
3. Click **Bind**. When **Bind Scanner Succeed** pops up, it is indicated that the binding is successful.
4. Start DentalDB program.
5. Create an order. Set the scanning type as **Digital Model Scan**, define restoration types, save and scan.
6. Start the scanning software for scanning and data processing.
7. After the data processing is complete, click **Go to Send** to open the send interface, and select the target lab.


Use DentalViewer


Through the DentalViewer, you can preview the edited scanning data directly, that is more convenient and fast.

On the start menu select **DentalLauncher** > **DentalViewer** to open DentalViewer as follows.











Operation Guide


1. Click  to add the scanning data generated by IntraoralScan.

 **Note**




support dragging files to the DentalViewer interface to edit the data.

1. Edit or optimize the scanning model data through the right toolbar.

Icon	Name	Description
	Edit	Details in Data Edit .
	Remove Isolated Data	Quickly delete small data isolated to the main model.
	Fill Holes	1 Click the icon to enter the fill holes interface. 2 Drag the slider to adjust the application range. 3 Click  to confirm the operation and exit; Click  to discard the process and exit.
	Texture	Click the icon to view the texture on the scanned model.
	Smooth	Click the icon to clean and improve the quality of the model.
	Fit View	Click the icon to adjust the size of the model to be adaptive to the interface.

3 Click  to save the edited model data.

Related operation

Icon	Name	Description
	Transparency	Adjustable transparency of model data. Drag the slider to the far left, and the model data becomes transparent.
	View	View the data from six different directions.
	Delete Selection	In the model list, delete the selected data file.

Care and Maintenance

Pre-cleaning, Disinfection, and Sterilization

The whole set of scanner, including scanner tip, scanner body, scanner cradle and the calibrator, requires proper care, cleaning, and handling.

As individual part may be processed differently, please read and follow the given Instructions to help you effectively.

1. Scanner body, cradle and calibrator maintenance
2. Scanner tip maintenance

Note

- All parts are provided non-sterilized. Please follow the Instructions before the first use.
- Follow the Instructions to pre-clean, disinfect, and sterilize each part of the scanner. Using other methods not approved by the Instructions may damage your scanner and void your warranty.
- Only disinfect or sterilize the specified part(s). Do not attempt to disinfect or sterilize all parts of the product. The Company is not liable for any damages due to improper disinfection and sterilization.
- To ensure safety and effectiveness, it is recommended to use equipment, materials, and disinfectants that have been approved by local regulatory authorities for sterilization and disinfection.

⚠ Caution

- Before pre-cleaning, disinfection and sterilization, please wear a pair of clean medical gloves.
- Ensure that you have completely cut off the power supply and all connections from the scanner.
- After sterilization, cool the scanner tip for a period of time to prevent possible heat injuries to the user and the patient.
- To prevent cross-contamination, proper pre-cleaning, disinfection and sterilization of the scanner after each use is necessary.
- When the scanner tip is detached from the scanner, always protect the subtle units and the inner optical components on the front end of the scanner body by putting on the supplied dust cap.

Scanner Body, Cradle and Calibrator Maintenance


The scanner body, cradle and calibrator require an intermediate-level disinfection.

⚠ Caution


Before disinfecting the scanner body, please ensure that the scanner tip is detached from the scanner body, and the scanner body is covered by the dust cap.

Follow the steps below to complete the disinfection:

1. Disconnect the power of the scanner (see more details in [Connection and disconnection](#)).
2. Hold the scanner tip firmly and then gently take the tip off from the scanner.




3. Store the detached tip in a safe place, e.g. a dental instrument tray.
4. Hold the supplied dust cap and blocks to the matching slots on the front end of the scanner body.
5. Use new cotton gauze moistened with 70%-75% solution of ethanol to wipe the surface of the scanner body.
6. When done, store the scanner body in a clean and safe place.
7. Use new cotton gauze moistened with 70%-75% solution of ethanol to wipe the surface of the cradle and the calibrator.
8. When done, store the cradle and the calibrator in a clean and safe place.

 **Caution**

- When detaching the scanner tip, do not put your fingers on the lens of the scanner tip. Otherwise, the lens may be damaged.
- Put the scanner into a dust-proof bag when it's not in use to avoid collision or accidental drop.
- Avoid using any kind of detergent as some detergents or surfactants might penetrate the surface of the scanner body and then damage the device.
- Please do not clean the outer units and inner optical components on the front end of the scanner with any sharp objects, which may cause scratches and damage to the scanner.
- Do not clean the intake and exhaust vents with any sharp objects or other such tools.

Scanner Tip Maintenance

The scanner tip is the most essential part of the scanner as it is inserted into the patient's mouth during scanning. Therefore, in order to prevent cross-contamination, the tip must be thoroughly cleaned and sterilized before and after it touches a patient.


 **Caution**

- When immersing the scanner tip into the disinfectant solution, please follow the Instructions of the disinfectant solution.
- The scanner tip can be sterilized under high temperature up to 100 times and must be disposed of afterwards.
- Before cleaning and disinfection/sterilization, please wear clean surgical gloves and goggles.

There are two effective and approved methods:

Method 1: Cleaning and high-level disinfection.


Method 2: Cleaning and high temperature autoclave sterilization.

 **Caution**


High-level disinfection and steam sterilization must NOT be combined.

Cleaning steps

1. Disconnect the power of the scanner, and pull the scanner tip off the scanner body.




2. Pre-clean the scanner tip with distilled water for 3 minutes to remove stains like the saliva or blood.
3. Brush the inside and the outside of the tip with enzyme surfactant for 3 minutes. Repeat the step for at least 2 times.
 - When cleaning the inside of the tip, insert the sponge brush into the tip from both the front and rear ends, and move the brush lightly in tiny circles.
 - When cleaning the outside of the tip, brush the surface back and forth lightly.
4. Rinse the tip thoroughly with distilled water for at least 3 minutes.
5. If there is any stain or fingerprint on the lens, repeat step 3 and 4.

 **Caution**

Rinse the tip with distilled water every time and discard the used water. Do not reuse the water for rinsing or any other purpose.


6. Dry the scanner tip with a soft lint-free cloth, and put it in a well-ventilated place to ensure it's totally dry, or put it in a dryer for 10 minutes.
7. Inspect the lens of the tip to make sure it is clean and free of damage.

 **Caution**

If the lens of the tip has cracks or scratches on it, contact your local distributor or service provider.

High-level disinfection steps

1. Carefully fill the container with phthalaldehyde at a concentration of 5.5g/L, and immerse the scanner tip in the disinfectant for 12 minutes.
2. Take out the scanner tip from the disinfectant, and rinse it with distilled water 3 times and 1 minute for each to remove disinfectant residues.

 **Caution**

Discard the used distilled water. Do not reuse the water for rinsing or any other purpose.

3. Flush the tip with distilled water for at least 3 minutes.
4. Dry the tip with a soft, lint-free cloth and put it in a well-ventilated place to make sure it is totally dry, or put it in a dryer for 30 minutes.
5. Inspect the lens of the tip to make sure it is clean and free of damage.
6. If immediate use of the scanner tip is required, reconnect it; if not, store it with other dental instruments, and make sure it is totally dry.

High temperature autoclave sterilization steps

1. Fill the scanner tip with medical gauze and seal it in the autoclave bag.
2. Put the wrapped scanner tip into a sterilizer for 30 minutes at 121°C under a relative pressure of 102.9kPa (or 4 minutes at 134°C under a relative pressure of 205.8kPa), and then dry it for 30 minutes.

Caution

After drying, cool the scanner tip to room temperature to avoid a scald.

3. Inspect the lens of the tip to make sure it is clean and free of damage.
4. If immediate use of the scanner tip is required, reconnect it; if not, store it with other dental instruments, and make sure it is totally dry.

Scanner Storage

In case you need to transport the device, we strongly recommend that you keep the original packaging after unpacking the Scanner. Shipping the device without its original packaging material may cause possible product damage and result in additional service fees.

If the original packaging is no longer available or damaged, carefully package each part of the scanner with bubble wrap to protect against any possible damage during transportation.

Storage for Transport

- Make sure that the scanner is clean before placing it in the original carry box/package to avoid any possible contamination.
- Place each part of the product, e.g. the tip, scanner body, cradle, power adapter, in the original package carefully and prevent kinks of the cable.
- Make sure that each cable is rolled up and tangle-free before placing it in the original carry box.
- Before closing the lid, make sure no part of the product is protruding from the package.

Daily and Long-term Storage

- Always place the scanner in the cradle when it is not in use.
- When the scanner tip is detached from the scanner body, always protect the subtle units and the inner optical components on the front end of the scanner by putting on the supplied dust cap.
- Ensure the scanner is clean before long-term storage.
- Avoid storing the scanner and accessories in areas of extreme temperatures or under direct sunlight.
- Before storing the scanner, make sure the scanner tip, scanner body and cradle are thoroughly dry.

Contact Us

Contact Us

Technical Support:

<https://support.shining3ddental.com/en/support/home>

Facebook:

<https://www.facebook.com/shining3d.dental/>

Instagram:

<https://www.instagram.com/shining3ddental/>

YouTube:

<https://www.youtube.com/@SHINING3DDental>

Address and phone number:

APAC Headquarters

Shining 3D Tech Co.,Ltd.

Hangzhou, China

P: +86-571-83821925

Email: sales@shining3d.com

No. 1398, Xiangbin Road, Wenyan, Xiaoshan, Hangzhou, Zhejiang, China, 311258

EMEA Region

SHINING 3D Technology GmbH.

Stuttgart, Germany

P: +49-711-28444089

Email: sales@shining3d.com

Americas Region

SHINING 3D Technology Inc.

San Leandro, United States

P: +1415-259-4787

Email: sales@shining3d.com

2450 Alvarado St #7,

San Leandro, CA 94577