

SUPPLIES

SOLUTIONS FOR HEALTH CARE PROFESSIONALS

**Welcome to our new 2020/2021
Henry Schein Brand Catalog, your resource
for the latest products from Henry Schein.**

- Henry Schein Brand products are designed to bring exceptional value to your practice
- We offer a comprehensive range that meet your needs; from gloves, to curing lights, restoratives to endodontic items, you will find many product options to enhance your productivity and increase patient satisfaction.
- Our products have met high standards of quality to earn the Henry Schein Seal of Excellence — your 100% guarantee of satisfaction

**Henry Schein is committed to helping you be more
successful, so that you can focus on delivering the best
quality patient care.**

You can Rely on Us.



Anesthetics & Needles	1
Impression Materials	3
Consumables	9
Infection Control	17
Surgical	29
Cements & Liners	33
Restoratives	37
Pins & Posts	45
Crown & Bridge	47
Finishing & Polishing	51
Preventives & Whitening	55
Hand Instruments	61
Rotary Instruments	71
Endodontics	95
Small Equipment	103
X-Ray & Imaging	107
Laboratory Products	111

We would like to introduce you to some of our new product ranges that you will find throughout the catalog

Nitrile Examination Gloves New Generation

Latex-free textured gloves.

Lower risk of allergic reactions.

Complies to EN 455-1, 2, 2, EN 374-1,2 and EN 4204.



Updated Glass Ionomer ranges in restoratives, baseliners and luting cements with the addition of capsules and light cured versions.



Natural Elegance Bulk Fill

Ceramic-based light-curing nanohybrid composite which enables curing at 4mm layer thickness.



MTA

Bioceramic root canal sealer with high radiopacity.



Fluoride Varnish

Desensitising varnish with 5% sodium fluoride. Available in 3 flavours.



Anesthetics & Needles



Anesthetics/Needle Containers	2
Syringes	2
Hypodermic Needles	2
Dental Needles	2

Anesthetics/Needle Containers



HS Sharps Container

- Single-use container for the safe disposal of clinical waste such as needles and scalpels.
- Universal needle disconnector on the lid so the container can be permanently locked for final disposal.

0.25L 9881898	1.5L 9881901	3L 9881908	6L 9881905
0.5L 9881899	2L 9881907	3.5L 9881903	10L 9881906
1L 9881900	2.2L 9881902	4L 9881904	

Syringes



Cartridge Syringe Chrome-Plated

- Chrome-plated brass syringe

Aspirating 'A' type Syringe 1.8ml
9006402



Cartridge Syringe Chrome-Plated

- Chrome-plated brass syringe

Aspirating 'C' type Syringe 1.8ml
9006403

Hypodermic Needles



HS Hypodermic Needles - Long Bevel

- Disposable, long bevel luer lock hypodermic needles.
- Precision lancet point according to ISO 9626 with exact, smooth bevel cut.
- Innovative anti-coring treatment of needle tube and consistent low penetration load of needle and skin.
- Needles are regularly biomechanically and clinically analysed to ensure product quality.
- Sterile - shelf life 5 years after sterilisation.

Needles x 100

21G x 40mm - Green 9003334	27G x 40mm - Grey 9006444
23G x 30mm - Blue 9003336	30G x 12mm - Yellow 9003630

Dental Needles



HS Disposable Dental Needles

- Metric Disposable Plastic hub needles.
- Painless penetration aided by long lancet point.
- Sterile.

Needles x 100

27G short x 25mm - Purple 9000250	30G X-short x 12mm - Blue 9000253
27G long x 42mm - Green 9000251	30G short x 16mm - Yellow 9000752
30G short x 21mm - Pink 9000252	

Impression Materials



Alginate	4
Addition Cured Silicone (VPS)	4
Condensation Silicone	5
Bite Registration	6
Mixing Guns & Tips	6
Impression Trays	7
Tray Adhesives	8
Solvents & Cleaners	8

Alginate

HS Maxima Krom Alginate Colour Change

- Fast-setting, dust-free alginate (Class A, Type 1) with colour change.
- Colour changes to violet to load tray, pink to indicate when ready to insert in patient's mouth, and white when ready to remove from mouth.
- Oral setting time of 60 seconds.
- Total setting time of 2 minutes 30 seconds (including working time).
- Stable for up to 48 hours when kept in a closed plastic bag.
- Vegetarian and vegan friendly. Gluten-free.




Intro Pack

Alginate Plus 2 x 500g, Scoop and Water Measure Packed in a Plastic Container
5720877

Refill 500g
5720878

HS Maxima Alginate Plus

- Fast-setting, dust-free alginate (Class B, Type 1)
- Ideal for study models and orthodontics.
- Stable for up to 48 hours when kept in a closed plastic bag.
- Mint flavour.
- Homogeneous.
- Vegetarian and vegan friendly. Gluten-free.



Intro Pack

Alginate Plus 2 x 500g, Scoop and Water Measure Packed in a Plastic Container
5720875

Refill 500g
5720876

Addition Cured Silicone (VPS)

HS VPS Hydro Putty

- Hydrophilic vinyl polysiloxane impression material.
- High tear resistance.
- Superb elastic recovery from deformation.
- Superior adhesion between putty and wash.
- Mixing time: 30 secs.
- Working time: fast 1 min 30 secs; regular 2 mins.
- Setting time: fast 2 mins 30 secs, regular 3 mins 30 secs.



Soft Putty 500ml/800g

Regular
9001448

Fast
9001450

HS VPS Hydro Wash Light Body

- Hydrophilic vinyl polysiloxane wash material in cartridge format.
- Thixotropic - flows easily under pressure.
- Superior adhesion between putty and wash.
- Superb elastic recovery from deformation.
- High tear resistance.
- Total working time (incl. mixing time): light body regular 2 min; light body fast 1:30 min; regular body 2 min.
- Time in mouth (minimum): light body regular 3:30 min; light body fast 2:30 min; regular body 3:30 min.
- Setting time (includes mixing time): light body regular 5:30 min; light body fast 4 min; regular body 5:30 min.



Cartridges 50ml x 2 plus Tips

Light Body Wash Regular Set
9001451

Light Body Wash Fast Set
9001453

Regular Body Wash
9001455

Cartridges 50ml x 4

Light Body Wash Regular Set
9001452

Light Body Wash Fast Set
9001454

Regular Body Wash
9001456

Addition Cured Silicone (VPS)



HS VPS Hydro 380

- For use in Pentamix®, MixStar®, Renfert® or other automatic mixing machines.
- 5:1 ratio addition-cured vinyl polysiloxane material.
- Good readability; precision detail.
- Hydrophilic.
- Smooth flowing.
- Excellent mechanical properties.
- Setting times (incl. working time): Soft Regular total 5.30 mins; Fast total 4.00 mins.

Cartridge 380ml x 2, Bayonet Rings x 2 and Dynamic Mixers x 20

Regular Set
9002133

Fast Set
9002134

Condensation Silicone



HS Eurosil Max 2 Putty

- Eurosil Max 2 is a C-silicone for either two-stage or single-stage technique.
- Very high initial flowability.
- Elastic recovery - >98%.
- Yellow - 75 Shore A.
- Soft White - 65 Shore A.
- Mixing time: 30 sec.
- Total working time: 1:15 min.
- Time in mouth: 3:15 min.
- Setting time: 4:30 min.

Putty

Yellow 900ml
9008532

Yellow 10kg
9008533

Soft - White 900ml
9008534

Soft - White 10kg
9008535



HS Eurosil Max 2 Correction

- Low viscosity C-silicone wash material suitable for two-stage technique and single-stage technique.
- Excellent precision detail and hydrophilic.
- Hardness - 35 Shore A.
- Mixing time: 30 sec.
- Total working time: 1:30 min.
- Time in mouth: 3:30 min.
- Setting time: 5 min.

Correction Paste

Blue 140ml
9008536

HS Eurosil Max 2 Catalyst

- Catalyst for use with all Henry Schein Eurosil Max, Max 2 and Eurosil Lab, Lab 2 C-silicones.
- Light pink gel.

Catalyst Paste

Pink 60ml
9008538

Bite Registration



HS VPS Premium Hydro Bite

- Outstanding hardness and rigidity
- Easy-to-use
- Thixotropic properties and a non-slumping consistency.
- Fast setting (60 secs).

Hydro Bite Cartridge 50ml x 2
9001268



HS VPS-CADCAM

- Scannable bite registration.
- VPS CADCAM bite registration for the CAD/CAM (CEREC®) machines.
- Reflection active and can be used for powder-free optical registration of antagonists for CAD/CAM and CIM (computer-integrated manufacturing) impression technology.
- Green in colour. Can be trimmed by stones or acrylic burs.
- Fast set working time: 30 secs. Setting time: 90 secs (incl. working time).

50ml x 2
9002623

Mixing Guns & Tips



HS T-Mixer Dispensing Gun

- Mixing gun 1:1/2:1 with comfortable design.
- Less pressure required.
- Robust construction.

T-Mixer Dispensing Gun
Dispensing Gun 1:1/2:1 Cartridges
9009653



HS T-Mixer 1:1 Tips

- Mixing tips and intraoral tips for use on A-silicone or bite registration cartridges.
- Colour of cap on cartridge indicates colour of mixing required.
- Up to 30% less waste material due to innovative and compact construction.

T-Mixer Tips Yellow x 50
9009620

T-Mixer Tips Pink x 50
9009621

T-Mixer Tips Green x 50
9009622

Intraoral Tips Yellow S50 x 100
9000768

HS VPS Hydro 380


- Accessories for HS VPS Hydro 380.



Accessories

Dynamic Mixers Set
9002657

Bayonet Rings Set
9002656



HS Dynamic Mixing Tips for Pentamix

- Mixing tips compatible with Pentamix 1 and 2 machines.

Mixing Tips Clear x 50
5720912

Impression Trays

Multi-Trays




HS Multi-Trays

- Master, counter impressions and bite registration in one easy step.
- Manufactured from thin plastic, with retention rims.

Anterior x 35 1007896	Posterior x 50 1009981
Quadrant x 35 1003919	Posterior Sideless x 50 1014060

Disposable Plastic Impression Trays



HS Thermoplastic Impression Trays with Integral Handle

- Perforated, light-blue plastic dentate impression trays with fixed handle.
- Strong rigid design.
- Thermoplastic can be custom adjusted by heating or trimming.

Upper x 12	Lower x 12
Large #1 9004563	Large #2 9004564
Medium #3 9004565	Medium #4 9004566
Small #5 9004567	Small #6 9004568
Anterior #9 9004571	Anterior #10 9004572
Partial x 12	
Upper Left/Lower Right #7 9004569	
Upper Right/Lower Left #8 9004570	



HS Disposable Impression Trays

- Rigid, plastic disposable perforated impression trays in edentulous and dentate formats.
- Removable, autoclavable metal handle.

Trays x 25 plus Handle		Edentulous Lower	No 14 Average 9006191
Edentulous Upper	No 13 Average 9006190	No 4 Large 9006183	No 18 Small 9006193
No 3 Large 9006182	No 17 Small 9006192	No 6 Average 9006185	Orthodontic Lower
No 5 Average 9006184	Orthodontic Upper	Dentate Lower	No 16 9003993
Dentate Upper	No 15 9003992	No 10 X-Large 9006187	
No 9 X-Large 9006186		No 12 Large 9006189	
No 11 Large 9006188			

Metal Trays



HS Stainless Steel Impression Trays

- Stainless steel impression trays with rimlock.
- Perforated or solid format.

Solid Upper	Solid Lower
No 2 Large 9001635	No 20 Large 1016043
No 3 Medium 1018054	No 21 Medium 1012373
No 4 Small 1014471	No 22 Small 1015255
No 5 X-Small 9001634	No 23 X-Small 9001636
Perforated Upper	Perforated Lower
No 1 Large 1010221	No 20 Large 1013664
No 3 Medium 1012094	No 21 Medium 1010433
No 4 Small 1013021	No 22 Small 1012003

Tray Adhesives



HS Alginate Adhesive

- Alginate adhesive with a fresh apple fragrance.
- Forms strong adhesive film for alginate.

Liquid 15ml
9002480

Spray 300ml
9002481

Solvents & Cleaners



EuroSept® Plus Orange Oil Cleaner/ Solvent

- Cleaning solution containing natural orange oils.
- Removes residues of zinc oxide, alginate and plaster, as well as impression materials, traces of blood and other contaminants.
- Cleans impression trays, instruments, glass plates, ceramic, stainless steel and other surfaces.
- Solvent contains less orange oil than the cleaner.

Cleaner
250ml
9882716
250ml with Spray
9882717

Solvent
250ml
9882719
500ml with Sprayhead
9882718




EuroSept® Plus Tray Cleaner Concentrate

- Mild alkaline liquid concentrate for cleaning impression trays, mixing accessories and dental instruments.
- Suitable for cleaning corrosion-resistant metals and sensitive materials such as anodised aluminium.
- Suitable for use in ultrasonic baths or as an overnight soak.
- Use at 5-10% concentration.

1L
9882707

5L
9882708

Consumables



Air/Water Syringes	10
Evacuation	10
Cotton Products	11
Paper Products	12
Beakers	13
Headrest Covers	13
Bibs & Aprons	13
Dispensers	14
Plasticware and Containers	15
General Disposables	15

Air/Water Syringes



HS Syringe Tips with Metal Insert - 3 in 1

- Disposable 3 in 1 syringes with metal tip suitable for most DCI push ring and Adec style air/water syringes which have not already been adapted.
- Excellent through flow.


Syringe Tips White x 250
9003774



HS Syringe Tips without Insert - 3 in 1

- Disposable air/water syringe, sturdy and resistant to breakage.
- Spray pattern is comparable to regular tip.
- Gentle on the patient's soft tissue.
- Syringe tip is easy on the o-ring in syringe head, preventing water cross-over.

Syringe Tip without Insert x 250
Transparent 76mm
9006798



HS Syringe Tips with Plastic Insert - 3 in 1

- Disposable air water syringe tip; sturdy and resistant to breakage.
- Spray pattern is comparable to regular tip.
- Gentle on the patient's soft tissue.
- Syringe tip is gentle on the o-ring in syringe head, preventing water cross-over.
- No adaptor/converter required on original equipment.
- Compatible with Adec style or DCI.

Syringe Tip with Plastic Insert x 250

Green 9006794	Purple 9006796
Pink 9006795	Yellow 9006797

Evacuation



HS Part Vented Evacuators

- Disposable S-Curve, vented 1 side.
- Use either way round.

White 11mm x 100
9000739

HS High Volume Evacuators

- Autoclavable.

Blue - 13cm, ø16mm x 10
9001370
Grey - 13cm, ø16mm x 10
9001371




HS Aspirator Tubes

- Plain aspirator tube.
- 135mm length.
- 11mm diameter.

11mm x 100

Black 9883270	Lilac 9883268
Light Blue 9883269	



HS Saliva Ejectors

- Flexible saliva ejectors in 4 colours with removable tips.
- 13cm.

Removable Tips x 100

Transparent 9883108	Green 9883112
Blue 9883110	Red 9883113
Yellow 9883111	



HS Surgical Aspirator Tips

- Stainless steel surgical aspirator tips with 3mm and 5mm diameter openings.
- Autoclavable.

3mm Diameter Tip
9002014

5mm Diameter Tip
9002015



HS Surgical Aspirator Tip

- Sterile surgical aspirators for use in general oral surgery.
- Two section removable blue tip gives either 1.2mm or 2.5mm opening.
- White tube is 4.8mm diameter.
- Ready-to-peel sterile packs.

2.5mm Tip x 20 plus Adaptor
9006996

Cotton Products

Cotton Rolls



HS Maxima® Cotton Rolls

- Cotton rolls with a firm composition to keep soft tissue away from teeth.
- Excellent absorption.

Size 1 x 870
9008526

Size 3 x 435
9008528

Size 2 x 600
9008527

Cotton/Gauze Balls



HS Cotton Wool Balls

- General use cotton wool balls.

Large x 250
9009536

Egg Size x 500
9008742

Gauze



HS EuroGauze Swabs - Non-Sterile

- Non-sterile cotton swabs - 17 threads per square inch.

12-Ply Gauze Swabs x 100

5 x 5cm
9003177

7.5 x 7.5cm
9003178

10 x 10cm
9003175

8-Ply Gauze Swabs x 100

5 x 5cm
9003180

10 x 10cm
9003179

Cotton Pellets



HS Cotton Pellets

- Cotton pellets for applying medicaments or to dry cavities.
- Approx. 600 pellets per pack.

Size 00 4g
9001879

Size 0 4g
9001880

Size 1 10g
9001881

Size 2 10g
9001882

Cotton/Gauze Balls



HS EuroGauze Balls

- Loosely rolled cotton gauze strips.
- Non-sterile.
- 100% fibre-free cotton.

Hazelnut Size x 500

6 x 9cm (length before rolling)
9008739

Walnut Size x 500

15 x 15cm (length before rolling)
9008740

Plum Size x 500

20 x 20cm (length before rolling)
9008741


Large x 250

30 x 30cm (length before rolling)
9008743

X-Large x 250

30 x 40cm (length before rolling)
9008744

Paper Products



HS C-Fold Towels

- Fits both C- and Z-fold dispensers.

White - 2-ply - 250 x 310mm -
16 sleeves of 152 towels
9003323



HS Interfold Towels

- Z-fold towels.
- Fits Z-fold dispensers only.

White - 3 ply - 320 x 220mm -
20 sleeves of 125 towels
9880905

White - 2-ply - 240 x 210mm -
20 sleeves of 160 towels
9003325



HS Facial Tissues

- Soft 2-ply white facial tissues.

210 x 205mm -
30 boxes of 100 tissues
9002377



HS Dry Cleaning Cloth

- Dry cleaning cloth, non-linting.
- Strong and absorbent, use with water or detergents.
- HTM 01-05 compliant.

White x 50
9002938



HS Toilet Rolls

- 100mm wide rolls.

White - 2-ply - 10cm x 200 Sheets - 4pk x 12
9003328

White - 3-ply - 10cm x 250 Sheets x 56 Rolls
9006481



HS Coreless Centrefeed Towel Rolls - Mini

- Coreless centrefeed towel roll.
- Perforated.
- Rolls out from the middle.

White - 2-ply - 230mm x 60m - Pack x 12
9003318



HS Coreless Centrefeed Towel Rolls

- Rolls out from the middle.

White - 2-ply - 200mm x 160m x 6
9003316



HS Disposable Tray Liners

- Filter paper lining for instrument trays.

Tray Liners 180 x 280mm x 250

White
9001318
Blue
9001017
Green
9001018

Yellow
9002521
Pink
9883610



HS Double Tray Paper

- Paper for lining trays and drawers in dental units.

Double Tray Paper 280 x 360mm x 250

White
9001401

Drinking Cups



HS Drinking Cups

- Disposable plastic drinking cups - tall 180ml.

Beakers x 3000

Blue	Green
9006532	9006534
Yellow	White
9006533	9006535

Headrest Covers



HS Headrest Covers

- Disposable headrest covers made from plastic.
- Easy-to-use dispensing box.

24 x 28mm x 250
1124873

Bibs & Aprons



HS Patient Apron

- Disposable 2-ply aprons on a roll.
- 1-ply tissue and 1-ply plastic.

Roll 2-ply x 80

White 54 x 58cm	Light Green 54 x 58cm
9001380	9001379
Light Blue 54 x 58cm	
9001378	



HS Patient Apron

- Disposable, polyethylene aprons on a roll.
- Size: 80 x 54cm.

Dark Blue x 200
9002949



HS Disposable Patient Bib Holder

- Disposable, soft paper bib connector with adhesive ends.

Disposable Bib Holders x 250
9009590



HS Dri-Gards

- Disposable patient bibs.
- 2-ply paper and 1-ply poly film.
- Absorbent textured finish.
- Size 33 x 45cm.

Bibs x 500

Blue	Lavender
9001031	9001038
Green	Yellow
9001032	9001039
White	Dark Blue
9001030	9001402
Rose	Aqua 40 x 45cm
9001033	876330
Aqua	White 40 x 45cm
9001034	876327



HS Contour Bibs

- Disposable bibs with contoured neck.

Bibs x 400

3-Ply paper, 1-Ply plastic - 46 x 56cm
9002613

Dispensers



HS Hand Towel Dispensers

- Dispensers for C-fold or Z-fold towels.
- C-fold dispenser for 9003323, 9792221, 9792218.
- C-fold dispenser dimensions - (W)27.5 x (H)37.4 x (Wall depth)10.5cm.
- Z-fold dispenser for 9003325, 9880905, 9792222, 9792224.
- Z-fold dispenser dimensions - (W)27.5 x (H)37.4 x (Wall depth)14cm.

Z-Fold
9008446

C-Fold
9008447



HS Towel Roll Dispensers

- Centrefeed dispenser.
- Plastic.
- Large dispenser uses paper towels 9003316, 9792215.
- Small dispenser uses paper towels 9003318, 9792217.

White Dispenser

Large (W)250 x (H)365 x (D)230mm
9008448

Small (W)180 x (H)345 x (D)172mm
9008449



HS Maxima® Cotton Roll Dispenser

- Dispenser suitable for sizes 1, 2 and 3 rolls.
- Simply press the top to dispense roll - protecting the other from dust or damp.
- Made from autoclavable plastic.
- Contents not included.

Dispenser with Pusher

White
9006049

Blue
9006050



Glove Dispenser

- Clear acrylic glove holders.
- Contents not included.

Single
9795054

Triple
9795056

Double
9795055




HS RoundTips Dispenser

- Contents not included.

Dispenser for Round Tips

White
9006468



HS Cotton Roll Dispenser

- Dispenser with drawer, suitable for sizes 1, 2, 3.
- Contents not included.

Dispenser with Drawer

Blue
9006048

Plasticware & Containers

Trays & Inserts



HS Plastic Tray Inserts

- Disposable tray inserts.
- Preformed plastic.

Standard x 50	Mini x 50
28 x 18cm	18 x 14cm
9000937	9001160

Denture & Retainer Boxes

HS Retainer Boxes Super Tuff

- Crush-proof retainer box with ventilation holes.
- 7.5 x 6.5 x 2.5cm deep.



Retainer Boxes x 10
Blue 1027932
White 1026785
Red 1027464



HS Retainer Boxes

- Shock resistance, stable resinretainer boxes for storage and transport of orthodontic and prosthetic works.
- Inner dimensions: (width & depth & height): 79 x 69 x 36 mm.

Retainer Boxes x 10	
White 9002210	Pink 9002213
Yellow 9002211	Blue 9002214
Orange 9002212	Green 9002215

General Disposables

HS Disposable Suction Filter

- Disposable filter for suction unit.
- Made from yellow plastic.
- Replace twice a week.
- Filter outer diameter 37mm ± 0.5mm, filter length 73-144mm, filter hole 6mm ± 1mm, inner diameter approx. 25mm ± 0.2mm.

Filters x 12 9006472



HS Disposable Pulp Kidney Dish

- Made of 100% pulp.
- Autoclavable.

Disposable Pulp Kidney Dish x 320 9009416



Endodontic Products From A Name You Can Trust.

Our Maxima K-files and H-Files are made in Switzerland with the finest standards of quality to provide maximum control without finger stress.

You can trust that Henry Schein Brand Maxima Endodontic Products provide unsurpassed quality—backed by the Henry Schein Seal of Excellence, your 100% guarantee of satisfaction.



Maxima K-Files		
	21mm	31mm
Assortment 10-35	100-9209	100-6972
Assortment 40-70	100-4460	100-8692
Maxima H-Files		
	21mm	31mm
Assortment 15-40	900-0678	101-3475
Assortment 45-80	101-3981	101-9055



Made in Switzerland

See pages 95-102 for our selection of Endodontic Products.



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!

Infection Control



Barrier Products	18
Masks	18
Gloves	19
Spittoon Cleaners	21
Instrument & Bur Cleaning & Disinfection	21
Impression Cleaning & Instrument Disinfection	22
Surface Disinfection	23
System Wipes	25
Washer Disinfection	25
Tray Cleaners	25
Orange Solvents/Cleaners	25
Hand Cleaning & Disinfection	26
Testing	27
Pouches & Cassettes	28
Miscellaneous	28

Barrier Products

Protective Sleeves



HS Barrier Sleeves

- Disposable barrier products for handpieces, 3-in-1 syringes, curing lights and operatory trays.

Air/Water Syringe Cover

254 x 64mm x 500
9009372

Slow Handpiece

38 x 203mm x 500
9009043

Curing Light Cover

Curing Light Tip for 7mm and 8mm x 400
9009042

Tray Barrier Sleeve x 500

356 x 267mm
9009376

Halogen Curing Light x 250
9009044

LED Curing Light x 250
9009045

Pen LED Curing Light x 250
9880090



HS Surface Barriers

- Disposable barrier products to aid in the prevention of cross-contamination.

Cover Film Clear with Adhesive Back x 1200

101 x 152mm
9009370

Light Handle Cover (T-Bar) x 500

155 x 107mm
9009371

Tube Cover

101mm x 365m Roll
9009046

Eye Protection



HS Faceshield Visor Refills

- Faceshield visor refills - use with any clip-style frames.
- Protects full face area.
- Wipe clean, antiglare transparent material.
- Excellent antifog performance due to thinner material (125µm)

Faceshields x 25

9009630



HS Faceshield Visor Frames

- Faceshield visor frame made from polycarbonate.
- Side clips to hold the visor.
- Autoclavable to 121°C.
- Do not clean with alcohol-based disinfectant solutions.

Frame

Green Transparent
9009596

Blue Transparent
9009597

White
9009598

Pink Transparent
9009599

Face Masks - Earloop



HS Maxima® Earloop Masks

- Earloop procedure masks.
- Triple layer, fibreglass-free material surrounded by soft, non-irritating, inner and outer facing.
- Latex-free.
- PFE > 99.9% at 0.01 microns.
- BFE ≥ 98% at 0.1 microns.
- Differential pressure < 29.4.
- EN 14683 Type IIR.

Masks x 50

Green
9002168

Blue
9009061

White
9009062

Pink
9009063

Yellow
9009064



HS Maxima® Earloop Masks - Sensitive

- Earloop mask with medical grade cellulose inner for sensitive skin.
- Free of dyes so less irritating to sensitive skin.
- Fluid repellent.
- Latex-free.
- BFE ≥ 98% at 3 microns.
- PFE > 97% at 1 and 0.1 microns.
- Differential pressure < 49.0.
- Splash resistance mm/Hg ≥ 120.
- EN 14683 Type IIR.

Masks x 50

White
9882179

Barrier Products

Face Masks - Earloop



HS Earloop Face Masks

- 3-ply latex-free masks.
- Good breathability.
- Adjustable nosepiece: soft and non-irritating for maximum comfort.
- BFE greater than 98%.
- EN 14683 Type II.

Masks x 50

Blue
9006146

White
9006148

Green
9006147



HS Maxima® Earloop Masks with Eyeshield

- Extra protection pleated mask provides superior fluid protection.
- Full coverage, optically clear shield is treated on both sides with antifog material.
- Light and comfortable.
- Latex-free.
- 3-ply fibreglass-free filter media surrounded by soft, non-irritating gum.
- BFE ≥ 98%.
- Differential pressure < 49.0.
- Splash resistance mm/Hg ≥ 120.
- EN 14683 Type IIR.

Teal Masks x 25
9009065



HS Maxima® Earloop Masks with Extra Protection

- Extra Protection Teal Mask.
- Provide the greatest filtration and fluid penetration resistance BFE ≥ 98%.
- Resistance to synthetic blood penetration at 160mm Hg.
- Light and comfortable for extended wear when necessary.
- Filter media is surrounded by soft, non-irritating inner and outer facing.
- Full-width flexible nose piece provides fit and security.
- Conform EN 14683 type IIR.
- Latex-free.

Teal Masks x 50
9009066

Face Masks - Tie-On



HS Tie-On Face Mask

- 3-ply latex-free mask.
- Bacterial filtration efficiency (BFE): 98%.
- Micro-spunbond polypropylene non-woven filter.
- EN 14683 Type II.

Masks x 50

Blue
9006143

Green
9006144

White
9006145



HS Maxima® Tie-On Surgical Masks

- Tie-on masks free of fibreglass and latex.
- Triple layer.
- Cloth-like ties.
- BFE ≥ 98% at 0.1 microns
- PFE > 99.9% at 0.01 microns
- Differential pressure < 29.4.
- EN 14683 Type IIR.

Masks x 50

Blue
9009068

White
9002165

Green
9002164

Pink
9002167

Gloves - Vinyl Powder-Free



HS Criterion Vinyl Powder-Free Gloves

- Vinyl gloves free from powder, latex and accelerators.
- Non-sterile.
- Skin friendly.
- Conforms to EN 455-1, 2, 3, EN 374-1, 2, and EN-420.

Gloves x 100

X-Small
9008167

Small
9008166

Medium
9008165

Large
9008164

X-Large
9008163

Barrier Products

Gloves - Latex Powder-Free



HS Aloe Vera Non-Sterile Latex Gloves Powder-Free

- Powder-free latex gloves containing aloe vera gel which moisturises and hydrates the skin.
- Excellent tactile sensitivity.
- Green in colour.
- Textured palms and fingertips.
- Manufactured to EN 455-1, 2, 3.

Gloves x 100

X-Small	Large
9002871	9002874
Small	X-Large
9002872	9002875
Medium	
9002873	



HS Economy Non-Sterile Latex Gloves Powder-Free

- Non-sterile gloves with textured surface for better grip.
- Guaranteed maximum protein level of 50mg per gram.
- Powder-free.
- Beaded cuff.
- Conforms to EN 455-1, 2, 3 and AQL 1.5

Gloves x 100

X-Small	Medium
1024184	1025289
Small	Large
1028304	1029168



HS Premium Non-Sterile Latex Gloves Powder-Free

- Chlorinated, polymer-coated, powder-free gloves.
- Guaranteed maximum protein level 50mg per gram or less.
- Textured and beaded cuff.
- Conforms to EN 455-1,2,3.

Gloves x 100

X-Small	Large
1019924	1010445
Small	X-Large
1018787	9001715
Medium	
1013806	

Gloves - Latex Surgical



HS Criterion Surgeons' Latex Powder-Free Gloves - Sterile

- Sterile surgeons' gloves - powder-free with a polymer coating for easy donning.
- Supersoft latex and anatomical design provides comfort and reduces fatigue.
- 3 year expiry date.
- Conforms to EN 455-1, 2, 3.

Gloves x 50 pairs

Size 5.5	Size 7.5
9883672	9883676
Size 6	Size 8
9883673	9883677
Size 6.5	Size 8.5
9883674	9883678
Size 7	
9883675	

Gloves - Nitrile



HS Nitrile Gloves Powder-Free

- Non-sterile nitrile powder-free gloves.
- Textured fingers.
- Latex-free.
- Conform to EN 455-1, 2, 3, EN 374-1, 2, and EN 420.

Gloves x 100

X-Small	Large
9880926	9880929
Small	X-Large
9880927	9880930
Medium	
9880928	



HS Extra Protection Long Nitrile Gloves

- Non-sterile, powder-free nitrile gloves with longer 16" cuff.
- Textured fingers.
- Latex-free.
- Single use.
- Colour: blue.
- Conforms to EN 455-1, 2, 3, and 374.3.

Gloves x 50

X-Small	Large
9880158	9880161
Small	X-Large
9880159	9880162
Medium	
9880160	

EuroSept® Plus Spittoon Disinfection

- Ready-to-use special alkaline cleaner for disinfection and cleaning maintenance of cuspidors.
- Reliably removes residues from limescale, blood and colourants from plaque tablets or mouthwash.
- Effective against bacteria and candida albicans and tuberculosis.

Maxima EuroSept®
Spittoon Foam Cleaner 400ml
9882691

Spittoon Disinfection 750ml
9882690




Limescale

EuroSept® Plus Surface Limescale

- Ready-to-use special cleaner for the removal of limescale deposits.
- Only for use on acid-resistant surfaces.

Surface Limescale 1L
9882677



Instrument & Bur Cleaning & Disinfection



HS EuroSept® Max Instrument Lemon

- Aldehyde-free concentrated instrument solution containing ethanol.
- pH value of 9.
- Used as a cold steriliser or in an ultrasonic bath.
- Renew solution once a week or more frequently if necessary.
- Dilution: 1.5% - 15ml makes 1 litre.
- 3% - 30ml makes 1 litre.
- HBV/HCV/HIV, bactericidal, fungicidal and tuberculocidal; effective within 60 minutes for 1.5% solution and 15 minutes for 3% solution.
- Material containing ABS should be placed in a 1.5% solution for 1 hour only.

5L
9003332



EuroSept® Plus Bur Disinfection

- Ready-to-use disinfection solution for burs and rotary instruments.
- For use in an ultrasonic bath or rotary container.
- Contains corrosion inhibitors.

Bur Disinfection 2L
9882709

Impression Cleaning & Instrument Disinfection



EuroSept® Plus Impression Disinfectant

- Ready-to-use impression disinfectant solution
- Do not mix with other cleaning agents or disinfectants
- Depending on the degree of contamination, the working solution should be replaced at least once per working day
- Bacteria (incl. MRSA) as per VAH, EN 13727 and EN 14561 - 1min; Candida albicans as per VAH, EN 13624 and EN 14562 - 1min; Enveloped viruses (e.g. HBV, HCV and HIV) as per RKI/DVV (with biological load) - 1min; (with contamination) TbB (M. terrae in a suspension test) - 15min; Noroviruses as per EN 14476 - 30min; Rotaviruses as per BGA/DVV - 1min

Impression Disinfection 5L
9883992

EuroSept® Plus Instrument Disinfection FV

- A fast acting, virucidal, liquid concentrate for disinfecting all types of medical instruments.
- Melon fragrance.
- Bactericidal (including MRSA), fungicidal: 2% - 15 min.; 1.5% - 30 min.; 1% - 60 min.; tuberculocidal (M. terrae): 2% - 15 min.; 1.5% - 30 min.; 1% - 60 min. virucidal including HBV/HIV, polio virus, adenovirus, papovavirus, vaccinia virus: 2% - 30 min.
- Conforms to the EU guideline 93/42/EC for medical devices.

Instrument Disinfection FV 5L
9882715

EuroSept® Plus Instrument Disinfection

- Material friendly, active-cleaning instrument disinfectant for routine clinical use.
- Excellent antimicrobial efficacy also gives a high degree of personal protection during instrument re-processing.
- Pleasant lemon fragrance.
- Aldehyde-free formulation.
- Suitable for use on all types of medical instruments.
- Bactericidal (incl. MRSA), fungicidal, virucidal (HBV/HIV/HCV viruses). 1% – 60 mins. 2% – 30 mins. 3% – 15 mins.
- Conforms to the EU guideline 93/42/EC for medical devices.

Instrument Disinfection 2L
9882712

Instrument Disinfection 5L
9882713

EuroSept® Plus Instrument Cleaner Enzymatic Concentrate

- Highly efficient, enzymatic cleaning concentrate for manual instrument reprocessing.
- Suitable for ultrasonic baths.
- Dilute with warm water at 1% dilution ratio (10ml per litre for instruments) and 0.5% (5ml per litre in ultrasonic baths).
- Exposure time is 5 minutes minimum.

Instrument Cleaner Enzymatic Concentrate 1L - dosing bottle
9882711

EuroSept® Plus Surface Alcohol-Free

- Alcohol-free, ready-to-use solution for disinfecting and cleaning of medical inventory, medical devices and other kinds of surfaces.
- Particularly recommended for alcohol sensitive medical devices and surfaces. Based on quaternary ammonium compounds.
- Contains no alcohols, aldehydes or phenols. Application and exposure: bactericidal: 1 min., effective against all enveloped viruses incl. HBV/HIV/HCV/ Vaccinia/BVDV/Influenza/flu viruses: 30 secs.
- Conforms to the EU Biocidal Product Regulation (BPR, 528/2012).
- Conforms to the EU guideline 93/42/EC for medical devices.

**Surface Alcohol-Free
1L W/Trigger**
9882672

**Surface Disinfection
Alcohol-Free 5L**
9882673



EuroSept® Plus Surface Disinfection Alcohol-Free Concentrate

- Alcohol, aldehyde and phenol free surface disinfection and cleaning in a concentrate.
- Perfume free.
- For use on medical devices and surfaces.
- Suitable for cleaning and disinfection of acrylic glass.
- Dilute at either 10, 20 or 25ml per litre.
- Bactericidal including MRSA, fungicidal, virucidal (HBV, HIV, HCV, BVDV, Vaccinia, Papova viruses): 2.5% 15mins, 2% 30mins and 1% 60mins.

**Surface Disinfection Alcohol-Free
Concentrate 1L - Dosing Bottle**
9882674

**Surface Disinfection Alcohol-Free
Concentrate 5L**
9882675

EuroSept® Plus Surface Disinfection FV

- Ready-to-use, alcohol based (ethanol) rapid disinfectant.
- Designed for applications requiring a targeted virucidal disinfection of alcohol resistant surfaces.
- Innovative SLT booster technology ensures an efficient disinfection in cases of high organic soiling (dirty conditions).
- Melon Fragrance.
- Application and exposure: Bactericidal incl. MRSA, tuberculocidal and virucidal incl. Polio-, Adeno-, Rota- and Noro viruses. Surface disinfection (wipe disinfection) (bactericidal, dirty conditions): 1 min.; Tuberculocidal (M. terrae) acc. to EN 14348: 30 sec.; Virucidal acc. EN 14476 incl. HBV/HIV/HCV: 1 min.; Polio- and Noro virus (MNV): 1 min.; Adeno- and Rota virus: 30 sec.
- Conforms to the EU Biocidal Product Regulation (BPR, 528/2012).
- Conforms to the EU guideline 93/42/EC for medical devices.

**Surface Disinfection
FV 1L Spray**
9882670

**Surface Disinfection
FV 5L**
9882671



EuroSept® Plus Surface

- Ready-to-use, aldehyde-free rapid disinfectant based on ethanol for an effective disinfection of medical inventory, medical devices as well as other kinds of surfaces.
- Especially suitable for difficult to reach surfaces.
- The aldehyde-free formulation avoids any discolorations of surfaces.
- Do not use on alcohol-sensitive surfaces (e.g. acrylic glass).
- Application and exposure time: bactericidal: 1 min., effective against all enveloped viruses incl. BV/HIV/HCV/ Vaccinia/BVDV/Influenza/flu viruses: 30 secs.
- Conforms to the EU Biocidal Product Regulation (BPR, 528/2012).
- Conforms to the EU guideline 93/42/EC for medical devices.

**Surface Spray 1L
With Trigger**
9882666

**Surface
Disinfection 5L**
9882667

Surface Disinfection



EuroSept® Plus Surface Lemon

- Ready-to-use, aldehyde-free rapid disinfectant based on ethanol for an effective disinfection of medical inventory, medical devices as well as other kinds of surfaces.
- Especially suitable for difficult to reach surfaces.
- The aldehyde-free formulation avoids any discolorations of surfaces.
- Do not use on alcohol-sensitive surfaces (e.g. acrylic glass).
- Lemon fragrance.
- Application and exposure time: bactericidal: 1 min., effective against all enveloped viruses incl. BV/HIV/HCV/Vaccinia/BVDV/Influenza/flu viruses: 30 secs.
- Conforms to the EU Biocidal Product Regulation (BPR, 528/2012). Conforms to the EU guideline 93/42/EC for medical devices.

Surface Lemon 1L With Trigger	Surface Disinfection Lemon 5L
9882668	9882669

Wipes

EuroSept® Plus Wipes/XL Wipes

- Ready-to-use alcohol-based wipes for disinfecting and cleaning near patient and hygienically critical areas.
- Does not leave any residue on surfaces.
- Lemon fragrance.
- Application and exposure time: bactericidal: 1 min.; tuberculocidal (M. terrae): 1 min; limited virucidal (including HBV, HIV, HCV), rotaviruses: 30 sec.; noroviruses (MNV): 1 min.
- Conforms to the EU Biocidal Product Regulation (BPR, 528/2012).
- Conforms to the EU guideline 93/42/EC for medical devices.

Wipe Refill 150pk	XL Wipe Refill 200pk
9882646	9882649
Wipe Tub 150pk	XL Wipe Tub 200pk
9882648	9882653



EuroSept® Plus XL Alcohol-Free Wipes

- Alcohol-free disinfection wipes for rapid disinfection and cleaning of medical inventory and all kinds of surfaces in medical areas.
- Based on quaternary ammonium compounds, free from alcohol, aldehydes or phenols.
- High-quality lint-free wipes suitable for all professional applications.
- Application and exposure time: bactericidal: 1 min., effective against all enveloped viruses incl. BV/HIV/HCV/Vaccinia/BVDV/Influenza/flu viruses: 30 secs.
- Conforms to the EU Biocidal Product Regulation (BPR, 528/2012).
- Conforms to the EU guideline 93/42/EC for medical devices.

Wipe Refill 200pk	Wipe Tub 200pk
9882654	9882658

System Wipes

EuroSept® Plus System Alcohol Free Wet Wipes

- Alcohol-free disinfection wipes for rapid disinfection and cleaning of medical inventory and all kinds of surfaces in medical areas.
- Based on quaternary ammonium compounds, free from alcohol, aldehydes or phenols.
- High-quality lint-free wipes suitable for all professional applications. Application and exposure time: bactericidal: 1 min., effective against all enveloped viruses incl. BV/HIV/HCV/Vaccinia/BVDV/Influenza/flu viruses: 30 secs.
- Conforms to the EU Biocidal Product Regulation (BPR, 528/2012).
- Conforms to the EU guideline 93/42/EC for medical devices.



EuroSept® Plus System Wet Wipes Rfl Alcohol-Free 32x30cm

9882665

EuroSept® Plus System Dispenser large (empty)for size 32x30cm

9882660

Washer Disinfection



HS EuroSept® Plus Thermotabs

- Enzymatic cleaning tabs for cleaning instruments and devices in a washer disinfectant or ultrasonic bath.

Tablets x 50

9882720

Tray Cleaners



EuroSept® Plus Tray Cleaner Concentrate

- Liquid concentrate for cleaning impression trays, mixing accessories and dental instruments.
- Dilute with water to the level required - this can be in 1% increments between 5 and 10% (50ml - 100ml per litre of water)

Tray Cleaner Concentrate 1L - Dosing Bottle

9882707

Tray Cleaner concentrate 5L

9882708

Orange Solvents/Cleaners



EuroSept® Plus Orange Oil Cleaner/ Solvent

- Cleaning solution containing natural orange oils.
- Removes residues of zinc oxide, alginate and plaster, as well as impression materials, traces of blood and other contaminants.
- Cleans impression trays, instruments, glass plates, ceramic, stainless steel and other surfaces.
- Solvent contains less orange oil than the cleaner.

Cleaner

250ml
9882716

250ml with Spray
9882717

Solvent

250ml
9882719

500ml with Sprayhead
9882718

Hand Cleaning & Disinfection

Hand Disinfection



EuroSept® Plus Hand Disinfection Liquid

- A virucidal alcohol based hand disinfectant to be rubbed into the skin, exceptionally gentle to the skin.
- Its user friendly formulation offers an extensive spectrum of efficacy, combined with excellent skin-care properties.
- Application and exposure: Bactericidal, tuberculocidal, virucidal (HBV, HIV, HCV) virucidal acc. to EN 14476.
- Hygienic hand disinfection (EN 1500): 3 ml - 30 secs.
- Effective against noroviruses (MNV) in 15 sec.
- Conforms to the EU Biocidal Product Regulation (BPR, 528/2012).

**Hand Disinfection Liquid 300ml
with Pump Dispenser**
9882679

Hand Disinfection Liquid 1L Flexo
9882681

Hand Disinfection Liquid 5 Litre
9882682



EuroSept® Plus Hand Disinfection Gel

- An alcohol based hand disinfectant gel free from cumulative longterm acting ingredients that might have an irritating or sensitising effect on the skin.
- Precise dosing and easier, non-drip application through a user friendly gel formulation.
- Compatible with all washing lotions. Free from aldehydes, phenols and quaternary ammonium compounds.
- Application and exposure: Bactericidal, tuberculocidal, virucidal (HBV, HIV, HCV). Hygienic hand disinfection (EN 1500): 3 ml - 30 secs.
- Conforms to the EU Biocidal Product Regulation (BPR, 528/2012).

Hand Disinfection Gel 100ml
9882683

Hand Disinfection Gel 500ml
9882685

Hand Disinfection Gel 1L Flexo
9882686

Hand Disinfection Gel 5 Litre
9882687

Wash Lotions

EuroSept® Plus Washlotion

- pH neutral washing lotion for hands and skin.
- Premium quality skin friendly ingredients and emollients to prevent the skin drying out.

Washlotion 500ml
9882699

Washlotion 5L
9882701



EuroSept® Plus Washlotion Sensitive

- pH neutral, soft and creamy washing lotion for hands and skin.
- Suitable for applying to healthy skin with friendly ingredients and perfume free to prevent allergies.
- Moisturising components and glycerine reduce skin dryness caused by washing.

Washlotion sensitive 500ml
9882702



Hand Cleaning & Disinfection

Hand Washes & Lotions

HS EuroSept® Plus Hand Cream

- High-quality, rich skincare cream for all skin types.
- Soothing and moisturising.
- Easily absorbed.
- Perfume 0.5%.

100ml
9008925




HS Hand Cream Lotion

- Lotion for all types of skin.
- Absorbs quickly without leaving any residue.

500ml
9008402



Testing




HS Bowie Dick Test

- Test set for steam penetration in the autoclave.
- Colour changes to green to confirm correct function.
- Conforms to ISO 18472.

Kit
Cassette plus Tests x 20
9008646

Refill
Tests x 20
9008647



HS Class 6 Indicator Tests Strips

- Steam sterilisation test strips.
- Monitors sterilisation process to ensure proper completion.
- 100 strips performs 200 tests.
- Colour changes to green after successful completion: 3.5 mins at 134°C and 15 mins at 121°C.
- Conforms to ISO 18472.

Strips x 100 (200 tests)
9880811



HS Helix Test

- Helix steam test strips for autoclaves.
- Latex-free self-adhesive backing for storage on files.
- Conforms to ISO 18472.

Kit
Holder plus Sheets x 250
9008649

Refill
Sheets x 250
9008650

Pouches & Cassettes



HS Self Seal Sterilisation Pouches

- Self-seal pouches made from polyethylene and paper.
- Colour indicator panel changes when sterilisation has taken place.

Sterilisation Pouches x 200

83 x 158mm 9002984	134 x 280mm 9002985	254 x 381mm 9002986
89 x 254mm 9002983	190 x 356mm 9882127	



HS Sterireels

- Sterilisation reels for use in a heat sealer.
- Made from plastic with paper backing.
- Selection of widths for various applications.

50mm x 200m 9002991	100mm x 200m 9002993	250mm x 200m 9880487
65 x 200mm x 100m 9002990	150mm x 200m 9880485	300mm x 200m 9880488
75mm x 200m 9002992	200mm x 200m 9880486	

Pouches & Cassettes



HS Instrument Cassettes

- Stainless steel locking instrument cassette tray with blue silicone holding rails.
- Fully vented design for optimal sterilisation.
- Upper and lower rails are interchangeable.
- Suitable for washer disinfectors and autoclaves.

5-Instrument Cassette 182.5 x 72.5 x 32mm 9008543	16-Instrument Cassette 182.5 x 279 x 32mm 9008545
10-Instrument Cassette 182.5 x 145 x 32mm 9008544	

Miscellaneous

HS Scrub Brush

- Autoclavable.



Scrub Brush
9008916

Surgical



Scalpels	30
Sutures	30
Syringes	31
Drapes, Gowns & Accessories	31
Patient Monitoring	32

Scalpels

Handles

HS Air Collection - Scalpel Handles

- Henry Schein Air Collection is a range of single and double ended instruments for operative and periodontal procedures.
- Lightweight, with an average weight of 12g.
- Stainless steel grade 304 handles and 420 tips.
- Suitable for washer disinfectors and autoclaves.

No 3
9008098

No 4
9008099



Sterile Blades



HS Scalpel Blades - Sterile Carbon Steel

- Scalpel blades made from carbon steel - sterile.
- Individually wrapped and sterilised by gamma radiation.
- Exceptionally strong and uniformly sharp for superior cutting.

Blades x 100

No 10	No 12
5720852	5720854
No 11	No 15
5720853	5720856

Blades x 10

No 10	No 12D
1126183	1126180
No 11	No 15
1126175	1126181
No 12	No 15C
1126176	1126186

Sterile Scalpels



HS Disposable Scalpels - Sterile

- Disposable scalpels with plastic handles and carbon steel blades.
- Sterile.
- Individually wrapped and sterilised by gamma radiation.
- Exceptionally strong and uniformly sharp for superior cutting.
- Lightweight scalpel handle.

Scalpels x 10

No 10	No 12	No 15
1126196	1126194	1126187
No 11	No 12B	No 20
1126189	1126188	1126193

Sutures

Absorbable

HS PGA Rapid

- PGA Rapid sutures are designed for short-term closures.
- They are processed to absorb quickly without tissue reaction.
- Ideal for skin and mucosa closings that may benefit from fast absorption.
- Maintains 50% strength at 5 days, 0% at 14 days and is essentially absorbed at 42 days.



Reverse Cutting 45cm x 12

3/0 - 3/8 Circle - 18mm
9003025
4/0 - 3/8 Circle - 18mm
9003027

Taper Point 45cm x 12

3/0 - 1/2 Circle - 19mm
9003029
3/0 - 1/2 Circle - 22mm
9003031
4/0 - 1/2 Circle - 19mm
9003028

HS PGA Violet

- PGA Violet sutures have a smooth, low friction proprietary coating which reduces tissue drag.
- Synthetic sutures ensure closure through the critical wound healing period.
- They retain 50% of their high initial strength at 14 days and 20% of their strength at 21 days.



Reverse Cutting 70cm x 12

3/0 - 3/8 Circle - 24mm
9003036
4/0 - 3/8 Circle - 18mm
9003034

Tapered Point 45cm x 12

4/0 - 1/2 Circle - 18mm
9003035

Tapered Point 70cm x 12

3/0 - 1/2 Circle - 22mm
9003037

Sutures

Non-Absorbable

HS Silk Black Sutures

- Braided black silk suture made from non-capillary silk with excellent handling characteristics.
- The braiding provides a uniform, smooth surface and greater tensile strength.
- For long-term closures.

Reverse Cutting 45cm x 12

3/0 - 1/2 Circle - 15mm
9003246

4/0 - 1/2 Circle - 20mm
9003038

4/0 - 3/8 Circle - 18mm
9880696

Taper Point 45cm x 12

3/0 - 1/2 Circle - 22mm
9003039

Tapered Cutting 45cm x 12

3/0 - 1/2 Circle - 18mm
9003046



Syringes



HS Sterile Luer Syringe - 2 Part

- Disposable 2-part luer syringe made from the highest quality plastic with excellent tolerances for the perfect extrusion of liquids.
- Extended graduation for accurate measurement.
- Single-use, individually wrapped.

Syringes - Luer x 100

1ml	2ml	5ml	10ml	20ml
9003016	9003017	9003018	9003019	9003020

Drapes, Gowns & Accessories



HS Sterile Drape Kit

- Sterile set of standard surgical drapes for oral surgery and or implantology procedures.

Drapes x 4: 45cm x 75cm x 2, 75cm x 90cm, 150cm x 155cm, Drapes for Tubing 8cm x 120cm x 2 and Tape 10cm x 50cm x 2
9002443



HS Surgical Drape

- Sterile blue surgical drape.
- Non-adhesive.

75 x 90cm
9881776



HS Sterile Surgical Gowns

- Sterile, single-use OP gowns.
- Non-woven, hydrophobic material SMS 50gms.
- Four ties system and tie-round type enabling sterility of the users back.
- Knitted elastic cuffs for greater comfort.
- Blue.

Medium	Large	X-Large
9008381	9008382	9008383

Drapes, Gowns & Accessories



HS Shoe Covers - Latex-Free

- Non-sterile, unisex shoe covers with non-skid design.
- Latex-free.

Blue x 50
1049452



HS Bouffant Cap

- Hair and head protection with breathable, low-linting design.
- Latex-free, elastic headband for reduced allergic reactions.
- Fluid repellent material for optimal protection.

Blue x 100
1046058

White x 100
1026614



HS Surgical Drape with Adhesive Strip

- Sterile blue drape with 5cm adhesive strip.

50 x 50cm
9881772

50 x 70cm
9881770

75 x 90cm
9881771

Patient Monitoring

Blood Pressure Monitors



HS Maxima® Blood Pressure Monitor

- Maxima aneroid blood pressure unit for professional use.
- Latex-free and comes with a double tubed cuff.
- Measuring unit is corrosion-free.
- Easy transportable case.

Blood Pressure Monitor
9003510

HS Digital Blood Pressure Unit

- Automatic digital blood pressure monitor with extra large LCD display.
- 2 user settings with 60 memory sets for each user.
- WHO blood pressure classification indicator.
- Date/time display.
- IHB (irregular heart beat) detector function.
- Alarm clock reminding function.
- Unit: 36 month warranty, cuff: 12 month warranty.



Unit

Blood Pressure Unit with Cuff 22-32cm,
AA Batteries x 4, IFU, Storage Bag and Guarantee Card
9008181

Accessories

Large Cuff
9008183

Medium Cuff
9008182

Pulse Oximeters

HS Finger Pulse Oximeters

- Finger Pulse oximeter with OLED display - a compact, portable and convenient method to measure the oxygen saturation in the blood and pulse rate.
- Not for continuous monitoring.
- Suitable for both adults and children.
- Contacts are made from medical rubber only which has no toxins and no harm to the skin.
- Uses 2 AAA batteries - gives up to 30 hours of use.



Finger Pulse Oximeter
9006368

Digital Ear Thermometers

HS Digital Ear Thermometer

- Easy-to-use digital thermometer.
- Results in 1 second, complete with fever alarm and memory option.
- Low battery indicator.



Thermometer
9006438


Covers x 150
9006439

Cements & Liners



Calcium Hydroxide Baseliners	34
Luting Cements	34
Glass Ionomer Baseliners	34
Resin Cements	35
Temporary Cements	35
Implant Cements	36
Temporary Filling Materials	36
Accessories	36

Calcium Hydroxide Baseliners



HS Calcium Hydroxide Temporary Root Canal Filling Material

- Temporary root-canal filling material that can be used for exudation, partial pulpotomy, and apexification.
- Ideal for internal and external root resorption.

Syringe 2g plus 20 Tips
9003663



HS Preline Plus

- Radiopaque calcium hydroxide cement liner.
- Rigid baseline and pulp capping agent to reduce chemical and thermal sensitivity by sealing dentine tubuli.

Light Cured Syringe 4g x 2 plus
Tips x 20
9001013

Base 13g, Catalyst 11g
9001012

Luting Cements



Maxima® Resin Modified Glass Ionomer Luting Cement

- Continuous fluoride-releasing glass ionomer with high adhesion.
- In addition to common areas of usage, also suitable for high-strength zirconium dioxide ceramics.
- Self-activating capsules.
- Low film thickness with good wetting behaviour and precise marginal fit.
- Resistant to moisture and acid to allow secure and quick bonding.
- Excess material is easy to remove.

Capsules x 50
Light Yellow
9004392



LGI Glass Ionomer Luting Cement Capsules

- Luting cement in self-activating capsules.
- Radiopaque fluoride-releasing materials.
- Good wetting behaviour with a low film thickness.
- High translucency.
- Suitable for luting of metal-based crowns, bridges, inlays and onlays.
- Also used for luting of high-strength all-ceramic crowns and bridges.

Capsules x 50
9004393

Glass Ionomer Baseliners

RGI Glass Ionomer Baseline

- Radiopaque, easy-to-mix, HEMA-free baseliner with excellent biocompatibility.
- Used under composite, amalgam, ceramic and metal restorations.
- May also be used for deciduous cavity restorations and glass ionomer/ composite 'sandwich' restorations.
- Light-cure but will also dual-cure.



Intro Kit
Powder 15g and Liquid 6ml
9883066

Refill
Powder 15g
9883067

Liquid 6ml
9883068

LGI Plus Glass Ionomer Luting Cement

- Glass ionomer luting cement for cementing crowns, bridges, inlays, onlays, pins, posts and orthodontic bands.
- High erosion-resistance and excellent biocompatibility.
- Easy to mix (powder and liquid).



Powder 35g and Liquid 20ml
9883056

Resin Cements

HS Resin Cement Self-Etching

- Dual-cured, self-etching resin cement for permanent adhesive.
- No need for etching, priming or bonding.
- For crowns and bridges, inlays and onlays, posts and cores, ceramic crowns and Maryland bridges, veneering of alloys.

Syringe Set

Smartmix Syringe 5ml x 2, Mixing Tips x 5
and Intraoral Tips x 5
9008244



HS Cement Mixing Tips - Brown

- Small mixing tips for use on automix cement cartridges.
- Brown in colour.

Blunt End Mixing Tips x 20
1048918

Pointed Mixing Tips x 20
9008131

Intraoral Tips for Cement x 20
9003061

Temporary Cements



HS Temporary Cement - Non-Eugenol

- Temporary cement without eugenol in paste-to-paste formulation.
- For temporary crowns or trial cementation of permanent restorations.
- Suitable for patients who are sensitive to eugenol.

Base 40g, Catalyst 14g
9003664

Temporary Cement Automix

- Zinc oxide temporary cement for provisional cementation of crowns, bridges, inlays and onlays, and endodontic access.
- Easy to apply due to smartmix syringe, with less waste due to short mixing tips.

Syringe 5ml x 2 plus 3 Tips

With Eugenol
9883490

Without Eugenol
9883491



Implant Cements



HS Premium Implant Cement

- Excellently adapted handling for Implantology
- Non-eugenol cement for implants and long-term temporary crowns.
- Strong retention, easy removal when required.

Smartmix 5ml plus Mixing Tips x 5
9008245

Temporary Filling Materials



HS Max-Filltemp

- Temporary filling material made with eugenol free, zinc oxide in a synthetic resin base.
- Good packable consistency, material sets on contact with saliva.
- Impermeable to root canal medicaments.
- Neutral pH.
- Fast setting.
- Releases fluoride.
- White in colour.

28g x 2
9001614

Accessories



HS Mixing Pads

- Poly-coated mixing pads for various applications.

Non-Skid Poly Pads

3.5 x 3.5cm x 100 Sheets x 12
1008309

7.5 x 7.5cm x 100 Sheets
1002880

7.5 x 15cm x 100 Sheets
1009180

15 x 15cm x 100 Sheets
1009265



HS Mixing Pads

Parchment Pads

10 x 12.5cm x 100 Sheets x 3
1009536

18 x 20.5cm x 100 Sheets x 3
1003951

5 x 5cm x 100 Sheets x 3
9883865




HS Transparent Mixing Pads

- Ideal for mixing all types of cements: transparent background provides a clearer view of how well the powder and liquid have mixed.
- Antislip base.
- All the benefits of a glass slab but in a disposable format.

7 x 8cm x 100 Sheets
9001314

14 x 8cm x 100
Sheets
9000938

Restoratives



Etch Gels	38
Bonds	39
Articulating Paper	39
Alloys & Accessories	40
Light-Cured Composites	40
Compomers	41
Glass Ionomer Filling Materials	42
Matrix Retainers & Bands	43
Pressure Indicators	43
Mixing Accessories	44
Restorative Mixing Guns & Tips	44
Interdental Wedges	44

Etch Gel



HS Maxima® Self-Limiting Etch Gel

- 35% phosphoric acid etch specifically designed to provide an optimum level of viscosity.
- Self-limiting in its depth of etch (average depth of 1.9 microns with 20 second etch).
- Viscous enough to prevent migration caused by gravity and also maintain a thick enough layer to prevent premature drying.
- Contains no glycerine and designed for rapid and complete removal upon rinsing.
- Suitable for etching before bonding composites, sealants or bonding adhesives.
- Blue in colour.

Syringe 60g
9880492



HS Jumbo Etch

- 40% phosphoric acid etch gel.
- 3cc syringes fill from the 60g syringe and are ideal for practices with multiple surgeries.
- Cost-effective.
- Blue in colour.

Jumbo Kit

Syringe 60g plus Syringes 3cc x 5 plus Tips
1025823

Refill

Refillable Syringes
3cc x 5
9883122

Accessories

Jumbo Syringe 60g
1026893



HS Etch Gel

- 40% water-based phosphoric acid gel with excellent viscosity for precise application.
- Syringe delivery system with applicator tips.
- Blue in colour.

Syringes 1.2 ml x 4 plus Tips x 20
1004649

Syringe 5g plus Tips x 3
9001200

Syringe 12g plus Tips x 3
1010711

Accessories

Syringe Tips x 100
9883123



HS Etch Gel for Ceramics (Silane)

- Etch Gel with hydrofluoric acid for use with ceramics, inlays, onlays and veneers.
- Easy, precise application giving a retentive etched surface.
- Composition: 1ml Etchgel for Ceramics contains 0.05g hydrofluoric acid (HF).
- Red in colour.

Bottle 6ml
9002441

Bond

HS Maxima® Natural Elegance 7th Generation Single Bond

- Light-cured, self-etching adhesive.
- Reduced evaporation.
- Excellent shear bond strength.



Bond 5ml
9009415

HS Natural Elegance Single Coat Bond 5th Generation

- Light-cured, multipurpose, one-component adhesive.
- Bond strength between 22-36MPa depending on use.
- Indications include: adhesive bonding of composite materials to natural enamel and dentine; adhesive bonding of compomer materials and dentine sealing.
- Solvent-free.



Bond 5ml
9003552



HS Natural Elegance Self-Etch Bond 6th Generation

- Light-cured, self-etching adhesive system consisting of primer and bond for adhesive restorations.
- Self-etch primer demineralises the dentine and enamel surfaces, solubilising the smear layer without removing it.
- Light-cured bond is solvent-free and filters through the collagen fibres giving a highly retentive and homogeneous bond layer.
- Can be used in dry, moist, wet areas of application.

Self-Etch Bond Kit	Primer	Bond
Primer 5ml, Bond 5ml 9003553	5ml 9003554	5ml 900355

Articulating Paper



HS Articulating Paper

- For marking contact areas.
- Smooth, cushioned surface.

14 Books x 10 Sheets

X-Thin 40 microns Blue
9000733

12 Books x 12 Sheets

Thick 195 microns Blue
9000734

Thin 45 microns Blue
9000737

Thin 45 microns Blue/Red
9000736

6 Books x 12 Sheets

Horseshoe 45 microns Red/Blue
9000735

Alloys & Accessories

Alloy Tablets/Capsules

HS Amalgam Non Gamma 2 Caps

- High silver (69.3%) amalgam capsules (19.4% tin, 10.9% copper, 0.4% zinc).
- Mixture of lathe cut and spherical particles.
- Highly polishable, corrosion resistant.
- Mixing time 5 seconds at 4000/4500 cpm.
- Condensation time < 4.3 minutes, carving time < 5.3 minutes.



Capsules 1 Spill x 50
9002128

Capsules 2 Spill x 50
9002129

Capsules 3 Spill x 50
9002130



HS Stratosphere 2

- Dispersed phase, high silver amalgam alloy with a mixture of lathe cut and spherical particles.
- Free of gamma 2.
- Mixing time over 4000rpm is 5-7 seconds for all spills.
- 69.2% silver, 18.6% tin, 11.9% copper and 0.3% zinc.

Capsules x 50

1 Spill -
Yellow/White
9009036

2 Spill -
Yellow/Grey
9009037

3 Spill -
Yellow/Black
9009038

Light-Cured Composites

Bulk Fill Composites



Natural Elegance Bulk Fill Composite

- Ceramic-based, light-curing, radiopaque nanohybrid composite.
- Used as a base in Class I and II restorations - can also be used in Class V restorations and for locking and splinting of loose anteriors.
- Universal shade allows curing of 4mm layer thickness (base fill) with low shrinkage - 1.25% by volume.
- Homogeneous material offering outstanding handling.
- Compatible with all conventional bonding agents.

Syringe 5g
Universal
9007691

Tips 0.25g x 15
Universal
9007692

Nanohybrid Composites



HS Maxima Natural Elegance NH Composite

- Highly polishable HCM (high cohesive molecular) resin-based nanohybrid composite.
- Indications include direct filling of all teeth, and the bonding of composite and ceramic restorations.
- Easy to condense and sculpt.
- Radiopaque.
- Excellent margin adaptation and wear resistance.
- Duo shades reduce the number of shades needed.

Syringes 4g

A1/D2
9003538

A2/B2
9003539

A3/D3
9003540

A3.5/B3
9003541

A4/M5
9003542

Bleach
9003543

Tips 0.25g x 20

A1/D2
9003545

A2/B2
9003546

A3/D3
9003547

A3.5/B3
9003548

A4/M5
9003549

C2/C3
9003811

Hybrid Composites



HS 20/20 Hybrid Composite

- Light-cured hybrid composite that can be polished to a high lustre.
- Excellent abrasion resistance combined with superior colour matching enables the composite to be used in both anterior and posterior restorations.

Syringe Refill 3g

A2
9006380

A3
9006381

A3.5
9006382

B1
9006383

B3
9006384

C3
9006385

Light-Cured Composites

Nanohybrid Flowable Composites



HS Maxima Natural Elegance NH Flow Composite

- Nanohybrid flowable composite with thixotropic formula for use in even the smallest cavities.
- Indications include Class V fillings, anterior fillings (Class III and IV), small posterior fillings and minimally invasive preparations.
- Can be used for composite and ceramic veneer repairs and the block-out of undercuts.
- Radiopaque.

Syringe 2.3g

A3/D3	A2/B2
9003855	9003854

Microhybrid Flowable Composites



HS Natural Elegance Flowable Composite

- Fluoride-releasing, light-cured, low viscosity composite designed for Class III, V and selected Class IV restorations.
- Flows easily in difficult areas.
- Excellent aesthetics.
- VITA shades.
- Radiopaque.

Intro Kit

Syringes 2g each of A2, A3, B3, C3 plus Tips x 20
1018074

Refill 2g x 4 plus Tips x 20

A1	A3	B3
1019306	1015337	1017388
A2	A3.5	C3
1016773	1012772	1013784

Microhybrid Composites

HS Natural Elegance Microhybrid Composite

- Resin-based microhybrid composite with excellent physical properties.
- For anterior and posterior restorations.
- High viscosity.
- Highly condensable and radiopaque.



Syringe Kit

A2, A3, A3.5, B2, C2, C3 and Universal Opaque Syringes, Flowable Composite, Sealant, Etch and Accessories
9000994

Syringe Refill 4.5g

A1	B2	A1	B2
9004001	9000998	9004004	9000990
A2	B3	A2	B3
9000995	9000999	9000987	9000991
A3	C1	A3	C1
9000996	9004003	9000988	9004006
A3.5	C2	A3.5	C2
9000997	90001014	9000989	9000992
B1	Universal/Opaque	B1	Universal/Opaque
9004002	9001001	9004005	9000993

Tip Refill 0.25g

A1	B2
9004004	9000990
A2	B3
9000987	9000991
A3	C1
9000988	9004006
A3.5	C2
9000989	9000992
B1	Universal/Opaque
9004005	9000993

Compomers



HS Compomer

- Light-cured filling material which combines the advantages of glass ionomer with the ease and aesthetics of composite.
- Suitable for cervical lesions, Class III and V restorations.
- For intermediate Class I and II restorations, and restorations in deciduous teeth.
- Fluoride release.

A1 0.25g x 20	A3 0.25g x 20
9008261	9008363
A2 0.25g x 20	A3.5 0.25g x 20
9008262	9008264

Glass Ionomer Filling Materials



Maxima® Reinforced Glass Ionomer LC

- Biocompatible fluoride-releasing glass ionomer.
- Suitable for larger cavities and the CBF (Composite Bonded to Flow) technique in deep cavities.
- Light-cured with good working time.
- Packable after placement, and easy to model and sculpt without sticking to the instrument - no varnish needed.

Capsules x 50

A2
9004394

A3
9004395

RGI Resin Reinforced Glass Ionomer

- Radiopaque resin-reinforced glass ionomer material.
- Used for deciduous cavity restorations, Class V cervical lesions, and single-surface Class I or II cavities in adult teeth.
- May be used for glass ionomer/composite 'sandwich' restorations and ART or MI procedures.
- Light-cure but will also dual-cure.
- Easy to mix (powder and liquid).
- High compressive and flexural strength with excellent shear bond strength.



Intro Kit

Powder A3 15g and Liquid 6ml
9883062

Powder 15g

A2
9883063

A3
9883064

Refill

Liquid 6ml
9883065

RGI Glass Ionomer XFast 2.0 Capsules

- Self-activating capsules with adjustable consistency due to variable mixing time.
- Easy to place in difficult-to-access or smaller, precise areas.
- Can be shaped and sculpted easily without sticking to the instrument.
- Excellent margin adaptation and compressive strength, which gives superior abrasion resistance.
- Setting time of 2 minutes from placement of the filling.



Capsules x 50

A2
9007998

A3
9007999

RGI Glass Ionomer Filling Material

- Radiopaque glass ionomer filling material with high compressive strength.
- Used for Class III and V cavities, posterior cavities in deciduous teeth and selected permanent teeth, 'sticky' pits and fissures, core build-up, and MI treatments.
- Easy to mix with water.
- Very good translucency and excellent biocompatibility.
- Good chemical adhesion to tooth structure.



Starter Kit

Shades A1, A2 and A3
9883057

Refill 15g


A1
9883058

A2
9883059

A3
9883060

RGI Glass Ionomer Fast 2.5 Capsules

- Biocompatible, radiopaque, high fluoride-releasing glass ionomer with aesthetic tooth-like translucency.
- Capsules require activating.
- Non-sticky consistency and very high compressive and transverse strength.
- Excellent abrasion resistance.
- Stable fillings with no marginal gaps.
- Setting time of 2.5 minutes: ideal for the treatment of children and anxious patients.



Capsules x 48

A3
9007997

Glass Ionomer Filling Materials



RGI Glass Ionomer Regular 3.5 Capsules

- Self-activating capsules with adjustable consistency due to variable mixing time.
- Easy to place in difficult-to-access or smaller, precise areas.
- Can be shaped and sculpted easily without sticking to the instrument.
- Excellent margin adaption and compressive strength for superior abrasion resistance.
- Setting time is 3.5 minutes from placement of the filling.

Capsules x 50

A2	A3
9004396	9004464



RGI Silver Reinforced Glass Ionomer

- Radiopaque silver-reinforced glass ionomer cement that is easy to mix with water.
- Can be used as a temporary or long-term treatment on Class I and III cavities in deciduous teeth, for the repair of amalgam restorations, and as a radiopaque posterior restorative.
- May be used for Class I cavities in selected permanent teeth.
- Suitable as a base under amalgam or posterior composite restorations, core build-ups, and for all cases of cavities where radiopacity rather than aesthetics is a prime requirement.
- Exceptional strength.
- Packable consistency and high wear and erosion resistance.

Powder 15g
9883061

Matrix Retainers & Bands



HS Matrix Strips

- Clear straight matrix, made from polyester fibre (mylar) resists solvent stretch.
- 0.002 gauge.
- Size: 10cm x 1cm.

Strips x 500
1008525



HS Tofflemire Retainer Bands

- Made from high-quality stainless steel.

Contra-Angle	Junior	Universal
9003695	9003696	9003694



HS Tofflemire Bands

- High-quality, medical grade stainless steel tofflemire bands.

Bands 0.002 x 12

No 1 Adult Universal 1003050	No 13 Paedo Universal 1000443
No 2 MOD Wide 1009722	No 14 Paedo MOD Narrow 1001412
No 3 Adult Narrow 1000139	No 15 Paedo MOD Wide 1004501

Pressure Indicators



HS Occlusion Spray

- Micropowdered indicator spray.
- Coats surface with an ultra-thin layer of pigment, which on contact will mark opposing model with colour.

Green 75ml 9001392	Red 75ml 9002390
------------------------------	----------------------------

Mixing Accessories



HS RoundTips

- Disposable and bendable applicator tips.
- 86mm long with round foam head.

RoundTips x 100

Regular Blue
9001393

Regular Green
9001394

Regular Silver
9001395

Dispenser for Round Tips
(Contents Not Included)

Dispenser
9006468

Regular Violet
9001396

Fine White
9003226

Superfine Yellow
9003227



HS Mixing Bowl

- Flexible white silicone bowls - autoclavable.

X-Small 5ml x 3
9001747

Small 15ml
9001743

Medium 25ml
9001744

Large 65ml
9001742

Restorative Mixing Guns & Tips



HS Compute Dispenser

- Autoclavable economy compule dispenser in compact size.

Dispenser
9882046



HS Natural Elegance Flowable Composite Tips

- Accessories for the Natural Elegance Flowable Composite range.

Flow Compotips 18G x 100
9883125

Interdental Wedges



HS Interdental Wedges

- Anatomical wooden interdental wedges.
- The square end makes it easy to manoeuvre into position with a plier.
- Made of high-quality Ahorn wood.

Interdental Wedges x 100



Peach-40
9002628

White-30
9002627

Pink-20
9002626

Green-10
9002625

Pins & Posts



Dentine Pins	46
Fibre Posts	46

Dentine Pins



HS Titanium Posts Dentipins

- Titanium dentine retention pins 0.6mm with automatic breakage.
- Made from non-corrosive titanium.
- Self-centering.
- Smooth cutting, no stress, soft shearing, great retention, superior stability.
- Machined from a solid rod of premium pure titanium.
- Part of the original rod of titanium is precision moulded into a shank for greater stability.

Titanium Starter Kit	Refill	
Pins x 25 plus Drill 9001407	Titanium Refill Pins x 100 9001406	Dentipin Drill x 6 9001408

Fibre Posts



HS Maxima® Fibre Posts


- Highly radiopaque, light-transmitting reinforced glass fibre posts.
- Colour coded for ease of use.
- Suitable for use with self, dual or light-cured cements.
- Cylindro-conical structure fits to the natural anatomy of the canal.
- Biocompatible due to epoxy-free composition.
- Superior stability.

Starter Kit
Posts No 1 - 4 (5 of each), Reamers No 1-4 (one of each)
9006456

Post Refill x 10	Reamers x 3
No 1 White Diameter 1.25 - 0.80mm 9006457	No 1 9006461
No 2 Yellow Diameter 1.50 - 0.90mm 9006458	No 2 9006462
No 3 Red Diameter 1.80 - 1.00mm 9006459	No 3 9006463
No 4 Blue Diameter 2.20 - 1.2mm 9006460	No 4 9006464



Crown & Bridge



Temporary Crown & Bridge Materials	48
Haemostatic Agents	48
Retraction	48
Temporary Crowns	49

Temporary Crown & Bridge Materials



HS MaxiTemp 1:1

- Crown and bridge acrylic in a 1:1 ratio.
- Self-cured material based on multifunctional methacrylic esters.
- Excellent marginal adaption.
- Highly polishable.
- Temporary restorations can be produced in the mouth or on a model.
- Working time: 30 secs.
- Curing time: 4 mins from start of mixing.
- Removal time from mouth: 1-1.
- Fully trimmable at 4 mins (including mixing time).

Intro Pack

Cartridge 50ml, Mixing Tips x 6

A2	A3
9009069	9009070

Accessories

Mixing Tips Blue 1:1 x 50
9003022



HS MaxiTemp HP 4:1

- Crown and bridge acrylic in 4:1 ratio automix cartridges.
- Self-cured material based on multifunctional methacrylic esters.
- Excellent marginal adaption.
- Highly polishable.
- Temporary restorations can be produced in the mouth or on a model.
- Working time: 45 secs.
- Curing time: 6 mins from start of mixing.
- Removal time from mouth: 2-3 mins (including mixing time).
- Fully trimmable at 6 mins (including mixing time).

Intro Pack

Cartridge 50ml, Mixing Tips x 6

A1	A2	A3	A3.5	Bleach
9001284	9001281	9002126	9001282	9004091

Refill Cartridge 50ml x 2

A1	A2	A3	A3.5
9001283	9001279	9002100	9001280

Accessories

Dispenser Gun 4:1
9001261

T-Mixer Tips Light Blue/Orange x 50
9009623

Haemostatic Agents



HS Haemostatic Gel

- Viscous 20% ferric sulphate haemostatic gel.
- Kind to soft and hard tissues.
- Eliminates fluid contamination for optimal impressions and bonding.

Syringe Kit

Syringe 30ml, Empty Syringes x 5 plus Bent Infusor Tips x 10
9003975

Retraction

HS Retraction Cord

- Knitted, non-impregnated retraction cord.
- Made from 100% knitted cotton fibres.
- Packed in a dispenser with cutting blade.

Size 00 Yellow 244cm
9009801

Size 0 Violet 244cm
9009802

Size 1 Blue 244cm
9009803

Size 2 Green 244cm
9009804





HS Polycarbonate Crowns

- Easily scissor-trimmed, contoured and crimped.
- Thin design prevents gingival irritation.

Assorted Kit
Crowns x 180 (3 of each size)
1000058

								
Lateral			Central			Cuspid		
Upper Right		Upper Left	Upper Right		Upper Left	Upper Right		Upper Left
No 2		No 25	No 100		No 15	No 300		No 35
1000184		1007601	1000059		1007505	1009398		1006476
No 20		No 26	No 101		No 16	No 301		No 36
1007747		1004232	1000188		1003926	1009960		1000645
No 21		No 27	No 10		No 17	No 30		No 37
1004697		1003491	1000586		1002516	1001470		1002056
No 22		No 28	No 11		No 18	No 31		No 38
1004973		1001129	1001868		1008020	1004686		1006554
No 23		No 29	No 12		No 19	No 32		No 39
1005853		1003598	1005954		1005567	1004938		1000894
No 24		No 200	No 13		No 102	No 33		No 302
1008658		1006085	1008980		1000599	1005742		1000061
			No 14		No 103	No 34		No 303
			1008628		1001909	1008305		1000194

								
Bicuspid			Incisor					
Upper Long		Lower Long	Lower Long		Lower Short			
No 40		No 50	No 60		No 65			
1002849		1002950	1008330		1009144			
No 41		No 51	No 61		No 66			
1009081		1009403	1006556		1009150			
No 42		No 52	No 62		No 67			
1008950		1009979	1000900		1009169			
No 43		No 53	No 63		No 68			
1008532		1001804	1002868		1009230			
No 44		No 54	No 64		No 69			
1007199		1005750	1009142		1009425			

Infection Control Products

From A Name You Can Trust.

Our trusted range of infection control products gives dentists and hygienists the confidence that they are performing every procedure with the right choice. You can trust that Henry Schein Brand infection control products provide unsurpassed quality—backed by the Henry Schein Seal of Excellence, your 100% guarantee of satisfaction.

See pages 17-28 for our selection of Infection Products.



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!

Finishing & Polishing



Denture Adjustment	52
Polishing Discs	52
Polishing Points	53
Polishing Pastes	53
Polishing Strips	53

Denture Adjustment



HS Denture Adjustment Kit

- For chairside denture adjustment.
- Green/grey 10,000 to 12,000rpm.
- Beige 5,000 to 7,000rpm.
- HP format.
- Bur in kit can be reordered under code 9001930.

Denture Adjustment Kit
9002104

Refill x 6

Grey HP 666
9002122

Beige HP 667
9002123

Green HP 668
9002124

Polishing Discs



HS Finishing and Polishing Discs

- Flexible polishing and finishing discs for composite restorations.
- Aluminium oxide coated.

Finishing and Polishing Disc Kit

9.5mm Discs x 120 (Coarse x 30, Medium x 30, Fine x 30 and X-Fine x 30),
12.7mm Discs x 120 (Coarse x 30, Medium x 30, Fine x 30 and X-Fine x 30) and Mandrel
9008850

Finishing and Polishing Disc Refill x 85

Coarse 9.5mm
9008851

Fine 9.5mm
9008853

Coarse 12.7mm
9008855

Fine 12.7mm
9008857

Medium 9.5mm
9008852

X-Fine 9.5mm
9008854

Medium 12.7mm
9008856

X-Fine 12.7mm
9008858

Accessories

Pop-On Mandrel
9008859



HS Finishing and Polishing Discs - Extra Thin

- Extra thin finishing and polishing discs for scaling, buffing and polishing.
- Ideal for incisal, cervical and approximal use.
- Aluminium oxide coated.

Finishing and Polishing Disc Kit

14mm Discs x 80 (Coarse x 20, Medium x 20,
Fine x 20, X-Fine x 20) and Mandrel
9008431

Finishing and Polishing Disc Kit

10mm Discs x 80 (Coarse x 20, Medium x 20,
Fine x 20, X-Fine x 20) and Mandrel
9008436

Finishing & Polishing Disc Refill XT x 80

Coarse 14mm
9008432

Fine 14mm
9008434

Coarse 10mm
9008437

Fine 10mm
9008439

Medium 14mm
9008433

X-Fine 14mm
9008435

Medium 10mm
9008438

X-Fine 10mm
9008440

Accessories

Pop On Mandrel XT
9008441

Polishing Points

HS Amalgam Polishing Kit

- A selection of polishers in RA format for polishing amalgam and precious metals.
- Black pre-polish.
- Brown polish.
- Green for final polish.
- Speed 5000 to 7000rpm.



Polishing Kit

Polishers x 12
9002102

Refill x 12

Black 0648 Torpedo
9002105

Black 0649 Cup
9002106

Brown 0650 Torpedo
9002107

Brown 0651 Point
9002108

Brown 0652 Cup
9002109

Green 0653 Torpedo
9002110

Green 0654 Point
9002111

Green 0655 Cup
9002112

HS Universal Silicone Polishing Kit

- Silicone polishing kit for resin-based material like composite, compomer, tooth enamel and porcelain.
- Medium grit.
- Right angle (RA) format.
- Speed 5000-12000rpm.
- HP format.



Polishing Kit

Polishers x 12
9002103

Refill x 12

0656 Mini Point
9002113

0657 Point
9002114

0658 Flame
9002115

0659 Torpedo
9002116

0660 Small Cup
9002117

0661 Wide Cup
9002118

0662 Eclipse
9002119

0663 Bullet
9002120

Polishing Pastes



HS Diamond Polishing Paste

- Polishing paste containing diamond particles.
- Produces high gloss on porcelain, composites, gold and other materials.

Paste 2g
9883522

Polishing Strips



HS Finishing and Polishing Strips

- Strips for finishing and polishing composite resins.
- Coated with medium and fine grits; grey and white colours respectively.
- Medium and fine aluminium oxide separated by non-coated area.
- 4 x 170mm.

Finishing and Polishing Strips x 100
9009646

Preventives Products From A Name You Can Trust.

Our trusted range of preventive solutions gives dentists and hygienists the confidence that they are performing every procedure with the right preventive choice. You can trust that Henry Schein Brand preventive products provide unsurpassed quality—backed by the Henry Schein Seal of Excellence, your 100% guarantee of satisfaction.

See pages 55-60 for our selection of Preventive Products.



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!

Preventives & Whitening



Caries Indicators	56
Pit & Fissure Sealants	56
Desensitisers	56
Mouthwash & Tablets	56
Fluoride Treatment	56
Prophylaxis	57
Interdental Cleaning	59
Toothbrushes	60
Whitening Treatment & Accessories	60

Caries Indicators

HS Caries Indicators

- Dye solution used to detect caries in dentine.
- Economical, safe, and effective formula in 15ml bottle with dropper cap.



Red 15ml
1028155

Green 15ml
1024453

Desensitisers



HS Natural Elegance Sealant Kit

- Light-cured opaque fissure sealant for an easy application on primary or permanent posterior teeth.
- 3ml etchant ensures a good flow into the pits and fissures for a long-lasting sealant.
- Low solubility.
- Natural Elegance etching liquid is recommended for use with this product.



Pit and Fissure Sealant Kit/Opaque

Syringe 1g x 4, Etchant Liquid 3ml plus Tips x 6
9001004

HS Desensitising Mousse

- Desensitising mousse contains sodium fluoride, potassium nitrate and 30% xylitol.
- Can be used in whitening trays or directly on the teeth.
- Provides care following tooth whitening procedure by remineralising teeth and relieving sensitivity.
- For use either by a dental professional in the surgery or provided to patient for home treatment of dentine sensitivity.

Kit

Syringe 1.2ml x 2 plus Tips x 2
9008817

Mouthwash & Tablets

HS Mouthwash Tablets

- Mouthwash tablets to be dissolved in warm water.

Tablets x 1000

Mint - Green
9006466

Thymol - Pink
9006467



HS Acclean® Alcohol-Free Mouthwash

- Ready-to-use mouthwash solution for professional cleansing and care of the mouth and throat area.
- Freshens breath while the antibacterial cleansing effect helps maintain oral hygiene and health.
- Rinse mouth with approximately 15ml of solution for 20-30 seconds - do not rinse mouth with water afterwards.
- For external use only, do not swallow.
- Chlorhexidine-free.
- Contains 0.1% benzalkonium chloride and 0.05% sodium fluoride.



250ml with dosing pump
9009663

Fluoride Treatment


HS Fluoride Varnish

- Fluoride-containing desensitising varnish with 5% sodium fluoride.
- 1 ml of product contains 50ml of Sodium Fluoride equivalent to 22.6 mg of Fluoride.
- For treatment of hypersensitive teeth and cervical areas after professional cleaning and calculus removal.
- Immediately releases fluoride and has good flow characteristics.
- Adheres well to moist surfaces and sets within seconds after contact with saliva.
- Does not contain any aspartame, saccharin or gluten.
- Naturally sweetened with xylitol.

Tube 10ml


Cherry
9883139

Mint
9883140



Prophylaxis

Prophy Brushes



HS Prophy Brushes

- Prophy brushes in a screw-fitting format.

Nylon Brush Screw
Type x 12
9008911

Bristle Brush Screw
Type x 12
9008912



HS Prophy Brushes - Latex-Free


- Latex-free prophy brushes in screw type format.
- Colour coded to indicate stiffness.

Brushes Screw-In x 100

Firm - Blue
9008904

Medium - Purple
9008905

Soft - Pink
9008906



HS Prophy Brushes - Latex-Free

- Latex-free prophy brushes in RA format.
- Colour coded to indicate stiffness.

Brushes RA x 100

Firm - Blue
9006524

Medium - Purple
9006525

Soft - Pink
9006526



HS Prophy Brushes

- Prophy brushes in RA format.

Bristle Brushes RA x 50
9000703

Nylon Brushes RA x 50
9000705

Prophy Angles

HS Acclean Prophy Angles

- Consistently smooth running and vibration-free prophy angles for optimum control.
- Disposable.
- Latex-free.
- Individually wrapped.



Soft Cup - Green x 100
9008926

Firm Cup - White x 100
9008927

HS Acclean Oscillating Disposable Prophy Angles

- Oscillating disposable prophy angles.
- Ergonomically designed angle providing a 90° reciprocating 'back and forth' rotating motion for continuous contact and pressure on tooth surface.
- Latex-free.



Soft Cup x 100
9009617

Firm Cup x 100
9009618

Prophylaxis

Prophy Cups

HS Prophy Cups - Latex-Free

- Latex-free, colour coded prophy cups.
- Made from synthetic polymer.



Cups RA x 100

Firm - Blue
9006521

Medium - Purple
9006522

Soft - Pink
9006523



HS Prophy Cups - Latex-Free

- Latex-free, colour coded prophy cups.
- Made from synthetic polymer.

Cups Screw Type x 100

Firm - Blue
9008907

Medium - Purple
9008908

Soft - Pink
9008909

Prophy Paste



HS Acclean Propy Paste

- Propy paste in easy dispensing tube.
- Less risk of cross-contamination.
- For use on natural teeth as well as composite and amalgam fillings.
- Mint flavour.

Tube 150g
9003196

Prophy Powder

HS Acclean Air Perio Polishing Powder Gentle, Extra Fine



- A prophy powder developed for gentle subgingival treatments
- For getting less paradontal issues and to strengthen the gums
- And for reducing bacteria even in the deepest periodontal pockets and reducing periodontal pockets depth
- Extra-fine grain: 25 µm

Polishing powder 100g
9003032

HS Acclean Lemon Air Polishing Powder



- Sodium hydrogen carbonate prophy powder in an ultra-fresh and zingy lemon flavour.
- Made from a well-balanced selection of particles and rounded crystals of various sizes and shapes.
- Quickly remove plaque, stains, biofilm and colouration without causing pain to the patient or damaging the enamel.
- Neutralises any acids, and reduces bacteria and possible gum inflammation.

300g
9002099

Interdental Cleaning



HS Acclean Flosser

- Interdental flosser with pick helps to effectively remove plaque.
- Floss slides easily between teeth without breaking or shredding.

Dispenser Box
48 Packs of 3

Mint 1126517	Unscented 1126518
-----------------	----------------------



HS Acclean Kid's Flosser

- Colourfully assorted, sea creature character design flossers help to promote healthy oral hygiene for children.
- Ergonomically designed for easy handling and reach for back teeth.
- Floss slides comfortably between teeth to help remove plaque and food debris.

Dispenser Box
48 packs of 3: total 144 flossers
1126519



HS Acclean Patient Floss

- Acclean floss is a nylon, shred-resistant floss for easy, comfortable removal of plaque and debris between teeth.
- Acclean Floss Plus is coated in PTFE which makes it easier to slide between the teeth.

Floss 12yd x 72	Floss 200yd
Waxed 1557084	Waxed 1554658
Waxed Mint 1550833	Waxed Mint 1554783
Unwaxed 1559860	Unwaxed 1552880



HS Interdental Brush

- Interdental brush for on-the-go removal of plaque between teeth.
- Each brush is individually wrapped.

Tapered Ultrafine Pink, 2,8mm 1026324	Cylindrical Ultrafine Green, 2,7mm 1026325
Tapered Medium Purple, 3,2mm 1028793	Cylindrical Medium Blue, 3,6mm 1026327

HS Acclean Patient Tape

- Shred-resistant tape.

Tape 100yd
Waxed
1555728




Toothbrushes

HS Acclean Denture Brush

- Bristles on 2 sides, 1 side tapered and the other serrated flat trim.
- Assorted handle colours.



Denture Brush x 12
1026790



Acclean Gentle Massage Toothbrushes

- Soft rubber stimulators align bristles to gently massage and stimulate gums.
- Diamond-polished, rounded, multileveled bristle ends remove plaque, while rubber bristles massage along gum line.
- Assorted colors.

Toothbrushes x 72/Box
9004067



HS Disposable Toothbrushes with Paste

- Disposable toothbrushes preloaded with paste.
- Simply wet brush to activate.
- Individually wrapped.
- Assorted colours.
- Fluoride-free.

Pack x 100
9880707

Whitening Treatment & Accessories



HS Natural Elegance Bleach

- Carbamide peroxide-based tooth whitening gel.
- Non-drip consistency for easy application.
- Formulated to help reduce sensitivity.
- 10% and 16% formulations.
- 16% formula suitable for both day and night use.
- Mild mint flavour.
- Initial treatment to be completed by dentist or under direct supervision to ensure patient safety.

Home Kit 10%
Syringe 1.2ml x 8 plus Tips x 8
9008873

Home Kit 16%
Syringe 1.2ml x 8 plus Tips x 8
9008874




HS Whitening Tray Box

- Boxes for the storage and transport of bleaching trays, bite trays and similar splints.
- An assortment of 20 boxes in 5 different colours (blue, white, green, bright green and pink - 4 of each).
- Dimensions: 7.3 x 7.8 x 3cm.

Assorted x 20
9009611

Hand Instruments




Diagnostic Instruments	62
Periodontal Instruments	63
Restorative Instruments	65
Surgical Instruments	67
Orthodontic Instruments	68
Miscellaneous	68
Instrument Trays & Accessories	70
Pouches & Cassettes	70

Diagnostic Instruments

Henry Schein Air Collection

- A range of single and double ended instruments for operatory and periodontal procedures
- Unique and exclusive handle design
- Light as air with an average weight of 12 grams
- Stainless steel grade 304 handles & 420 tips
- Suitable for both washer disinfectors and autoclaves

Mirrors



HS Premium Mouth Mirrors

- Autoclavable stainless steel mirrors.


Front Surface	Magnifying	Plane
No 3 x 12 9000747	No 4 x 12 9000798	No 2 x 12 9880892
No 4 x 12 9000748	No 5 x 12 9000799	No 3 x 12 9000753
No 5 x 12 9000749		No 4 x 12 9000754
		No 5 x 12 9000755

Mirror Handles



HS Air Collection - Mirror Handle

Round
9008012



HS Resin Mirror Handles

- Ergonomic design.

Black 9001899	Blue 9001438	Green 9001905
------------------	-----------------	------------------



HS Mirror Handles


- Stainless steel mirror handles.

Octagonal Serrated S/S
9009470

Round Knurled S/S
9009471

Probes/Explorers

HS Air Collection - Probe Single Ended



Diagnostic Instruments

Probes/Explorers

HS Air Collection - Perio Probes

Double Ended



23/WHO
9008009



23/PR11
9008010



23/PR12
9008011

Single Ended



CPITN
9008017



PCP11
9008018



PCP12
9008019



PCP14
9008020



PCP15
9008021



Williams Round
9008022

Periodontal Instruments

Gracey Curettes

HS MaxiGrip Ergo Gracey



Gracey DE 00
9007827



Gracey DE 1/2
9007823



Gracey DE 7/8
9007824



Gracey DE 11/12
9007825



Gracey DE 13/14
9007826



Gracey DE 13/14ER
9007829



Columbia 4R/4L
9007832



Columbia DE 13/14
9007831

HS Air Collection - Gracey Curettes



1/2
9008047



3/4
9008050



5/6
9008051



7/8
9008052



9/10
9008108



11/12
9008048



13/14
9008049



15/16
9008109



17/18
9008110

Periodontal Instruments

Implant Scalers

HS Titanium Implant Scalers

- Blue-tipped titanium scalers.
- Designed for debriding implants and abutments.



Scalor Barnhart 5/6
9007561



Scalor 204S
9007562



Scalor H6/H7
9007563



Scalor 4R/4L
9007564



Scalor Langer 1/2
9007565

Scalers

HS Air Collection - Scalers

Single Ended



U15
9008059

Double Ended



204S
9008053



204SD
9008054



H6/7
9008055



S5/33
9008056



S6/S7
9008057



U15/30
9008060



Micro 301-302SI
9008100



Micro 301-302XSI
9008101

HS MaxiGrip Ergo Scalers



Barnhardt DE 1/2
9007928



Barnhardt DE 5/6
9007830



Younger Good DE 7/8
9007833



Scalor DE H6/7
9007834



Scalor DE H5/33
9007835



Scalor DE 204S
9007836



Scalor DE 204SD
9007837

Restorative Instruments

Amalgam Pluggers

HS Air Collection - Amalgam Pluggers



Small Smooth
9008064




Medium Smooth
9008063




Large Smooth
9008062

Carvers

HS Air Collection - Carvers



Hollenbach 1
9008076



Hollenbach 2
9008077



Hollenbach 3
9008107

Burnishers

HS Air Collection - Ball Burnishers



1L
9008068



X-Large 2.7/3.2mm
9008069



Large 2.0/2.5mm
9008070



Medium 1.2/1.6mm
9008071



Small 0.8/1mm
9008072

Anodised Composite Instruments

HS Air Collection - Composite Instruments - Anodised



F4
9008065



F5
9008066



F6
9008067

Restorative Instruments

Excavators

HS Air Collection - Excavators



123/124
9008037



125/126
9008038



127/128
9008039



129/130
9008040



131/132
9008041



141/142
9008042



153/154
9008043



Large
9008044



Medium
9008045



Small
9008046

Cement Spatulas

HS Air Collection - Cement Spatulas



Small
9008031



Medium
9008030



Large
9008029

Plastic Filling Instruments

HS Air Collection - Plastic Filling Instruments



1
9008094



4
9008085




5
9008086



6
9008087



18
9008095



49
9008096

Wax Instruments

HS Air Collection - LeCron Carver



LeCron
9008025

HS Wax Carver



LeCron Curved 17cm
9002006

HS Iwanson Measuring Device



Wax Measuring Gauge
9002266

Metal Measuring Gauge
9002265

Elevators

HS Periosteal Elevator - Williger

- For use during tooth extraction.



Williger 16cm
9001745

Extraction Forceps

HS Extracting Forceps

- High-quality stainless steel extraction forceps with superior grip.



Extracting Forceps Kit

Upper EF1, EF2, EF7, EF17, EF18, EF51A, EF67A
and Lower EF13, EF22, EF33A, EF79

9001053

HS Extracting Forceps

- High-quality stainless steel extraction forceps with superior grip.

Upper



EF1 Upper Incisal/Cuspid
9001045



EF2 Upper Incisal
9001046



EF7 Upper Bicuspid
9001047



EF17 Upper Molar Right
9001042



EF18 Upper Molar Left
9001043



EF51A Upper Root
9001044



EF67A Upper 3rd Molar
9001052

Lower



EF13 Lower Bicuspid/Cuspid
9001050



EF22 Lower Molar Left/Right
9001048



EF33A Lower Incisal Root
9001049



EF79 Lower 3rd Molar Left
9001051

Surgical Instruments

Forceps

HS Tissue Forceps

- For tissue handling during oral surgery.



Toothed 13cm
9002009



Toothed 18cm
9002010

Retractors

HS Tissue Retractor

Langenbeck



10 x 30mm
9002007



11 x 50mm
9002008

Orthodontic Instruments

Pliers

HS Mathieu Tie Pliers

Mathieu Tie Pliers 14cm
9002079



Miscellaneous

Articulating Forceps

HS Air Collection - Articulating Forceps

- Extremely lightweight, stainless steel handled instruments.

Miller
9008015



Miscellaneous

Calcium Hydroxide Placements

HS Air Collection - Calcium Hydroxide Placement Instrument



CaOH₂
9008023

Cord Packing Instruments

HS Air Collection - Cord Packing Instruments



Circle
9008032



Circle Serrated
9008033

Needle Holders

HS Needle Holder



Castroviejo 14cm
9008637



Castroviejo 17cm
9008638

Scissors

HS Scissors



Serrated Straight 11.5cm
9001997



Serrated Curved 11.5cm
9001998

Tweezers & Pliers

HS Air Collection - Tweezers

- Extremely lightweight, stainless steel handled instruments.



317 Serrated
9008013



Self-Locking-Serrated
9008016


HS Pliers

- Made from stainless steel to provide toughness and durability.
- Designed for various applications in the dental laboratory.



Serrated 13.5cm
9002060

Instrument Trays & Accessories



HS Aluminium Instrument Trays

- Solid or perforated aluminium trays.
- Not suitable for washer disinfectors.
- Size: 18 x 28cm.

Solid Trays

Silver
1005433
Blue
1008936
Green
9001274
Yellow
9001273

Perforated Trays

Silver
9001276
Blue
9001272
Green
9001270
Red
9001271

Tray Lid - Solid

Silver
1025760

Tray Lid - Perforated

Silver
9001433

Tray Insert

15 Instruments
9001434



HS Stainless Steel Instrument Trays

- Stainless steel instrument trays in standard 18cm x 28cm and small size 14cm x 18cm.
- Suitable for automated washer disinfectors.
- Instrument frames for 4, 8 or 16 instruments.

Trays 18 x 28cm

Base
9001431

Lid
9001429

Trays 18 x 28cm

Base
9001430

Lid
9001428

18 x 24cm

16 Instruments
9001482

Pouches & Cassettes



HS Instrument Cassettes

- Stainless steel locking instrument cassette tray with blue silicone holding rails.
- Fully vented design for optimal sterilisation.
- Upper and lower rails are interchangeable.
- Suitable for washer disinfectors and autoclaves.

5-Instrument Cassette

182.5 x 72.5 x 32mm
9008543

10-Instrument Cassette

182.5 x 145 x 32mm
9008544

16-Instrument Cassette

182.5 x 279 x 32mm
9008545

Rotary Instruments



Diamond Burs	72
Steel Burs	87
Tungsten Carbide Burs	88
Tungsten Carbide - Surgical Burs	93
Tungsten Carbide - Sterile Oral Surgery Burs	93
Bur Brushes	93

Diamond Burs

Maxima Burs


HS Maxima® Diamond Burs

Round



Round x 5

● FG 801-009 8782723	● FG 801-014 8782729	● FG 801-018 8783554	● FG 801-023 8782740	● FG 801-029 8783558
● FG 801-010 8782724	● FG 801-014 8782727	● FG 801-018 8782733	● FG 801-023 8782737	● FG 801-029 8783559
● FG 801-012 8782725	● FG 801-014 8782728	● FG 801-018 8782734	● FG 801-023 8782739	● FG 801-029 8782743
● FG 801-012 8782726	● FG 801-014 8783552	● FG 801-018 8783555	● FG 801-023 8782738	● FG 801-029 8782742
	● FG 801-016 8782732	● FG 801-021 8783556	● FG 801-023 8782741	
	● FG 801-016 8782730	● FG 801-021 8782735	● FG 801-025 8783757	
	● FG 801-016 8782731	● FG 801-021 8782736	● FG 801-025 8783758	
	● FG 801-016 8783553	● FG 801-021 8783557	● FG 801-025 9008867	



All Maxima® Diamonds
are Colour Coded for
easy identification!

Grit Guide

●

 Super Coarse

●

 Coarse

●

 Medium

●

 Fine

●

 Super Fine

○

 Ultra Fine

Round with Collar



Round with Collar x 5

● FG 802-010 8782747	● FG 802-014 8782751	● FG 802-018 8782753
● FG 802-012 8782749	● FG 802-016 8782752	

Inverted Cone



Inverted Cone x 5

● FG 805-009 8782755	● FG 805-012 8782758
● FG 805-010 8782756	● FG 805-016 8782761
● FG 805-012 8782757	

Inverted Cone with Collar



Inverted Cone with Collar x 5

● FG 806-010 8782766	● FG 807-012 8782771
● FG 806-012 8783759	● FG 807-016 8783605
● FG 806-012 8782767	● FG 807-016 8782772
● FG 806-014 8783654	● FG 807-018 8782773
● FG 806-014 8782768	
● FG 806-016 8782769	

Diamond Burs

Maxima Burs

Round End Taper

Round End Taper x 5



● FG 849-009 8783695	● FG 849-025 8783795
● FG 849-009 8783696	● FG 849-025 8783796
● FG 849-010 8783698	● FG 849-025 8783798
● FG 849-010 8783699	
● FG 849-010 8783700	

Round End Taper

Round End Taper x 5



● FG 855-012 8783634	● FG 855-014 8783637	● FG 855-025 8783704
● FG 855-012 8783635	● FG 855-014 8783638	● FG 855-025 8783705
● FG 855-012 8782882	● FG 855-014 8782883	● FG 855-025 8783706
● FG 855-012 8783636	● FG 855-014 8783639	

Round End Taper

Round End Taper x 5



● FG 850-010 8783764	● FG 850-016 8782874
● FG 850-010 8783765	● FG 850-016 8782872
● FG 850-010 8783766	● FG 850-016 8782873
● FG 850-012 8783626	● FG 850-016 8783628
● FG 850-012 8782868	● FG 850-018 8782876
● FG 850-012 87832867	● FG 850-018 8783629
● FG 850-014 8783627	● FG 850-018 8782875
● FG 850-014 8782869	● FG 850-018 8783630
● FG 850-014 8782871	
● FG 850-014 8782870	

Round End Taper

Round End Taper x 5




● FG 856-010 8783707	● FG 856-016 8782894	● FG 856-021 8783708
● FG 856-012 8783564	● FG 856-016 8782891	● FG 856-021 8783709
● FG 856-012 8782884	● FG 856-016 8782893	● FG 856-021 8783710
● FG 856-012 8782886	● FG 856-016 8782892	● FG 856-025 8783711
● FG 856-012 8782885	● FG 856-018 8782898	● FG 856-025 8783712
● FG 856-014 8782890	● FG 856-018 8782895	● FG 856-025 8783713
● FG 856-014 8782887	● FG 856-018 8782897	
● FG 856-014 8782889	● FG 856-018 8782896	
● FG 856-014 8782888		

Round End Taper -Long

Round End Taper x 5



● FG 856L-014 8783714	● FG 856L-016 8783718	● FG 856L-018 8783722	● FG 856L-020 8783767
● FG 856L-014 8783715	● FG 856L-016 8783719	● FG 856L-018 8783723	● FG 856L-020 8783768
● FG 856L-014 8783716	● FG 856L-016 8783720	● FG 856L-018 8783724	● FG 856L-020 8783769
● FG 856L-014 8783717	● FG 856L-016 8783721	● FG 856L-018 8783725	



All Maxima® Diamonds
are Colour Coded for
easy identification!

Grit Guide	
●	Super Coarse
●	Coarse
●	Medium
●	Fine
●	Super Fine
○	Ultra Fine

Diamond Burs

Maxima Burs

Flat End Taper

Flat End Taper x 5



● FG 845-014
8782836

● FG 845-014
8782837

● FG 845-016
8782838

Flat End Taper

Flat End Taper x 5



● FG 846-012
8782840

● FG 846-012
8782841

● FG 846-012
8783461

● FG 846-014
9008868

● FG 846-016
8782842

● FG 846-016
8782843

Flat End Taper

Flat End Taper x 5



● FG 847-010
8783686

● FG 847-010
8783687

● FG 847-012
8783794

● FG 847-012
8782844

● FG 847-012
8782846

● FG 847-012
8782845

● FG 847-014
8782848

● FG 847-014
8783562

● FG 847-014
8782847

● FG 847-014
8783563

● FG 847-016
8782852

● FG 847-016
8782849

● FG 847-016
8782851

● FG 847-016
8782850

● FG 847-018
8783688

● FG 847-018
8783689

● FG 847-018
8783690

Flat End Taper

Flat End Taper x 5



● FG 848-010
8783691

● FG 848-010
8783692

● FG 848-012
8783693

● FG 848-012
8783694

● FG 848-014
8783623

● FG 848-014
8782856

● FG 848-016
8782860

● FG 848-016
8782857

● FG 848-016
8782859

● FG 848-016
8782858

● FG 848-018
8783624

● FG 848-018
8782861

● FG 848-018
8782862

● FG 848-018
8783625

● FG 848-023
8782863

● FG 848-023
9008869

Modified Taper

Modified Taper x 5



● FG 847KR-016
8783622

● FG 847KR-016
8782853

● FG 847KR-016
8782855

● FG 847KR-016
8782854


Modified Taper

Modified Taper x 5



● FG 846KR-016
8783683

● FG 846KR-016
8783684



All Maxima® Diamonds
are Colour Coded for
easy identification!

Grit Guide

- Super Coarse
- Coarse
- Medium
- Fine
- Super Fine
- Ultra Fine

Diamond Burs

Maxima Burs

Round End Cylinder

Round End Cylinder x 5



● FG 838-009
8782827

Round End Cylinder

Round End Cylinder x 5



● FG 880-012
8783514

Round End Cylinder

Round End Cylinder x 5



● FG 881-012
8783516

● FG 881-012
8783517

● FG 881-014
8783518

Flat End Cylinder

Flat End Cylinder x 5



● FG 835-008
8782797

● FG 835-009
8782799

● FG 835-010
8783611

● FG 835-010
8782800

● FG 835-010
8782801

● FG 835-010
8783612

● FG 835-012
8783613

● FG 835-012
8782802

● FG 835-012
8782803

● FG 835-012
8783614

● FG 835-014
8782804

● FG 835-014
8782806

● FG 835-014
8782805

Flat End Cylinder

Flat End Cylinder x 5



● FG 836-010
8783668

● FG 836-010
8783669

● FG 836-012
8783615

● FG 836-012
8782809


● FG 836-012
8782810

● FG 836-012
8783616

● FG 836-014
8782813

● FG 836-014
8782811

● FG 836-014
8782812



All Maxima[®] Diamonds
are Colour Coded for
easy identification!

Grit Guide

- Super Coarse
- Coarse
- Medium
- Fine
- Super Fine
- Ultra Fine

Flat End Cylinder

Flat End Cylinder x 5



● FG 837-010
8783674

● FG 837-010
8783675

● FG 837-012
8783619

● FG 837-012
8782816

● FG 837-012
8782818

● FG 837-012
8782817

● FG 837-014
8782821

● FG 837-014
8782819

● FG 837-014
8782820

● FG 837-016
8783763

● FG 837-016
8782822

● FG 837-016
8782823

Diamond Burs

Maxima Burs

Modified Cylinder



Modified Cylinder x 5

● FG 835KR-010 8783760	● FG 835KR-012 8783661
● FG 835KR-010 8782808	● FG 835KR-012 8783662

Modified Cylinder



Modified Cylinder x 5

● FG 836KR-010 8783761	● FG 836KR-012 8783617	● FG 836KR-014 8783671
● FG 836KR-010 8783762	● FG 836KR-012 8782815	● FG 836KR-014 8783672

Modified Cylinder



Modified Cylinder x 5

● FG 837KR-012 8783620	● FG 837KR-014 8783789	● FG 837KR-016 8783791
● FG 837KR-012 8782824	● FG 837KR-014 8783790	● FG 837KR-016 8783792
● FG 837KR-012 8782825		

Modified Cylinder



Modified Cylinder x 5

● FG KS1 8783748
● FG KS2 8783749
● FG KS3 8783750

Bevelled Cylinder



Bevelled Cylinder x 5

● FG 877-009 8782960	● FG 877-012 8782962
● FG 877-010 8782961	● FG 877-012 8782963
● FG 877-010 8783781	● FG 877-012 8783782

Bevelled Cylinder



Bevelled Cylinder x 5

● FG 878-012 8783581	● FG 878-014 8782979	● FG 878-016 8783783
● FG 878-012 8782974	● FG 878-014 8783582	● FG 878-016 8782980
● FG 878-012 8782976	● FG 878-014 8782978	● FG 878-016 8783784
● FG 878-012 8782975	● FG 878-014 8782977	

Bevelled Cylinder



Bevelled Cylinder x 5

● FG 879-012 8783787	● FG 879-014 8783586	● FG 879-016 8783788
● FG 879-012 8782996	● FG 879-014 8782997	● FG 879-016 8782999
	● FG 879-014 8782998	● FG 879-016 9008872
	● FG 879-014 8783587	

Bevelled Cylinder



Bevelled Cylinder x 5

● FG 884-012 8783520
● FG 884-012 8783521
● FG 884-012 8783648

Diamond Burs

Maxima Burs

Bevelled Cylinder



Bevelled Cylinder x 5

● FG 885-010 8783802	● FG 885-012 8783804	● FG 885-014 8783820
	● FG 885-012 8783522	● FG 885-014 8783821
	● FG 885-012 8783523	● FG 885-014 8783822
	● FG 885-012 8783649	● FG 885-014 8783823

Bevelled Cylinder



Bevelled Cylinder x 5

● FG 886-010 8788840	● FG 886-014 8783524
● FG 886-012 8783742	● FG 886-014 8783525
● FG 886-012 8783743	
● FG 886-012 8783744	

Pointed Football



Pointed Football x 5


● FG 368-016 8783590	● FG 368-018 8783751	● FG 368-023 8782708
● FG 368-016 8782700	● FG 368-018 8783752	● FG 368-023 8782705
● FG 368-016 8782702	● FG 368-018 8783753	● FG 368-023 8782707
● FG 368-016 8782701	● FG 368-018 8783754	● FG 368-023 8782706
● FG 368-016 8782703	● FG 368-018 8783755	● FG 368-023 8782709
○ FG 368-016 8782704	○ FG 368-018 8783756	○ FG 368-023 8783591

Football



Football x 5

● FG 369-025 8783825
● FG 369-025 8782711



All Maxima® Diamonds
are Colour Coded for
easy identification!

Grit Guide

- Super Coarse
- Coarse
- Medium
- Fine
- Super Fine
- Ultra Fine

Football



Football x 5

● FG 379-012 8782712	● FG 379-016 8783592	● FG 379-018 8783596	● FG 379-023 8783826
● FG 379-014 8782714	● FG 379-016 87823593	● FG 379-018 8783597	● FG 379-023 8782718
● FG 379-014 8782713	● FG 379-016 8782716	● FG 379-018 8782717	● FG 379-023 8782720
	● FG 379-016 8782715	● FG 379-018 8783598	● FG 379-023 8782719
	● FG 379-016 8783594	● FG 379-018 8783599	● FG 379-023 8782721
	○ FG 379-016 8783595	○ FG 379-018 8783601	○ FG 379-023 8783602

Diamond Burs

Maxima Burs

Flame



Flame x 5

FG 860-010 8782921	FG 860-012 8782923	FG 860-014 8783575
FG 860-010 8782922	FG 860-012 8782925	FG 860-014 8782926
	FG 860-012 8782924	FG 860-014 8782927
	FG 860-012 8783574	FG 860-014 8783576
		FG 860-014 8783577

Flame



Flame x 5

FG 861-012 8783734	FG 861-014 8783774
FG 861-012 8783735	FG 861-014 8783775
FG 861-012 8783736	

Flame



Flame x 5

FG 862-009 8783777	FG 862-010 8783578	FG 862-012 8782933	FG 862-014 8782938	FG 862-016 8782941
FG 862-009 8783778	FG 862-010 8782929	FG 862-012 8782930	FG 862-014 8782935	FG 862-016 8782939
FG 862-009 8783779	FG 862-010 8782928	FG 862-012 8782932	FG 862-014 8782937	FG 862-016 8782940
	FG 862-010 8783579	FG 862-012 8782931	FG 862-014 8782936	FG 862-016 8783642
		FG 862-012 8782934	FG 862-014 8783641	FG 862-016 8782942
		FG 862-012 8783640		

Flame




Flame x 5

FG 863-010 8783737	FG 863-012 8783643	FG 863-014 8783780	FG 863-016 8782951
FG 863-010 8783738	FG 863-012 8782944	FG 863-014 8783739	FG 863-016 8782948
	FG 863-012 8782946	FG 863-014 8783740	FG 863-016 8782950
	FG 863-012 8782945	FG 863-014 8783741	FG 863-016 8782949
	FG 863-012 8782947		FG 863-016 8782952

Grit Guide

- Super Coarse
- Coarse
- Medium
- Fine
- Super Fine
- Ultra Fine



All Maxima® Diamonds are Colour Coded for easy identification!

Diamond Burs

Maxima Burs

Flame

Flame x 5

- FG 888-012
8783650
- FG 888-012
8783627
- FG 888-012
8783651
- FG 888-012
8783652


Flame

Flame x 5

- FG 889-009
8783528
- FG 889-009
8783530
- FG 889-009
8783529
- FG 889-009
8783653
- FG 889-010
8783531
- FG 889-010
8783533
- FG 889-010
8783532
- FG 889-010
8783534

Grit Guide

- Super Coarse
- Coarse
- Medium
- Fine
- Super Fine
- Ultra Fine



All Maxima® Diamonds
are Colour Coded for
easy identification!

Curettage

Curettage x 5

- FG 876K-012
8782958

Curettage

Curettage x 5

- FG 877K-012
8782964
- FG 877K-012
8782965
- FG 877K-014
8783580
- FG 877K-014
8782966
- FG 877K-014
8782968
- FG 877K-014
8782967
- FG 877K-016
8782972
- FG 877K-016
8782969
- FG 877K-016
8782971

Curettage

Curettage x 5

- FG 878K-012
8783583
- FG 878K-012
8782982
- FG 878K-012
8782983
- FG 878K-012
8783584
- FG 878K-014
8783585
- FG 878K-014
8782984
- FG 878K-014
8782986
- FG 878K-016
8782990
- FG 878K-016
8782987
- FG 878K-016
8782989
- FG 878K-016
8782988
- FG 878K-018
8782993
- FG 878K-018
8782991
- FG 878K-018
8782992
- FG 878K-021
8783785
- FG 878K-021
8782994
- FG 878K-021
8783786

Curettage

Curettage x 5

- FG 879K-014
8783647
- FG 879K-014
8783504
- FG 879K-014
8783505
- FG 879K-016
8783509
- FG 879K-016
8783506
- FG 879K-016
8783508
- FG 879K-016
8783507
- FG 879K-018
8783513
- FG 879K-018
8783510
- FG 879K-018
8783512
- FG 879K-018
8783511

Barrel

Barrel x 5

- FG 811-033
8783606
- FG 811-033
8783607
- FG 811-033
8783540
- FG 811-037
8783608
- FG 811-037
8783609
- FG 811-037
8783541

Diamond Burs

Maxima Burs

Needle



Needle x 5

FG 858-010 8783726	FG 858-012 8783565	FG 858-014 8782906
FG 858-010 8783727	FG 858-012 8782901	FG 858-014 8782903
FG 858-010 8783728	FG 858-012 8782900	FG 858-014 8782905
	FG 858-012 8782902	FG 858-014 8782904
		FG 858-014 8782907
		FG 858-014 8783566



All Maxima® Diamonds are Colour Coded for easy identification!

Grit Guide

- Super Coarse
- Coarse
- Medium
- Fine
- Super Fine
- Ultra Fine

Needle



Needle x 5

FG 859-010 8783770	FG 859-012 8783729	FG 859-014 8783771	FG 859-016 8783571	FG 859-018 8782917
FG 859-010 8782911	FG 859-012 8783730	FG 859-014 8782914	FG 859-016 8783772	FG 859-018 8782919
FG 859-010 8782910	FG 859-012 8783731	FG 859-014 8782913	FG 859-016 8783572	FG 859-018 8782918
FG 859-010 8782912	FG 859-012 8783732	FG 859-014 8782915	FG 859-016 8782916	FG 859-018 8782920
	FG 859-012 8783733	FG 859-014 8783569	FG 859-016 8783573	

Donut



Donut x 5

FG 909-037 8783745	FG 909-040 8783588
FG 909-037 8783746	FG 909-040 8783589
FG 909-037 8783747	FG 909-040 8783549
	FG 909-040 8783824

Xmas Tree



Xmas Tree x 5

FG 852-010 8783701	FG 852-012 8783631
FG 852-010 8783702	FG 852-012 8782879
FG 852-010 8783703	FG 852-012 8782878
	FG 852-012 8783632
	FG 852-012 8783633

Diamond Burs

Maxima Burs

Pear



Pear x 5

FG 830-009 8782782	FG 830-010 8782783	FG 830-012 8782785	FG 830-014 8782787	FG 830-021 8804396
FG 830-010 8782784	FG 830-012 8782786	FG 830-014 8782788		

Long Pear



Long Pear x 5

FG 830L-012 8782791	FG 830L-014 8783560	FG 830L-016 8782796
FG 830L-012 8782792	FG 830L-014 8782793	FG 830L-025 9008870
FG 830L-014 8782795		
FG 830L-014 8782794		

Depth Marker



Depth Marker

FG 828-026 8783658
FG 828-030 8783659
FG 834-021 8783546

Amalgam Buster



Amalgam Buster x 5

FG 813-012 8783655	FG 813-014 8783610	FG 813-016 8782776
FG 813-012 8783656	FG 813-014 8782775	

Tissue Protective End



Tissue Protective End x 5

FG 840-012 8783678
FG 840-014 8783680
FG 840-016 8783681

Interproximal



Interproximal x 5

FG 392-016 8782722
FG 392-016 8783603

Palatinal



Palatinal x 5

FG 899-031 8783548

End Cutting



End Cutting x 5

FG 839-012 8782832
FG 839-014 8783677
FG 839-016 8783793

Diamond Burs

Maxima Miniature Head Burs

HS Maxima® Miniature Burs



All Maxima® Diamonds are Colour Coded for easy identification!

Grit Guide

- Super Coarse
- Coarse
- Medium
- Fine
- Super Fine
- Ultra Fine

Flat End Taper



Flat End Taper x 5

● FG 846-012
9000919

Flat End Cylinder



Flat End Cylinder x 5

● FG 805-008
9000917

Football



Football x 5

● FG 379-009
9000915

X-Mas Tree



X-Mas Tree x 5

● FG 852-007
8783865

Inverted Cone



Inverted Cone x 5

● FG 805-010
9000916

Round End Cylinder



Round End Cylinder x 5

● FG 838-007
9000918

Turbo Diamond Burs

HS Maxima® Turbo Diamond Burs

Flat End Taper



Flat End Taper x 5

● FG 847-016
8783850

Flat End Cylinder



Flat End Cylinder x 5

● FG 837-014
8783845

Round End Taper



Round End Taper x 5

● FG 850-014
8783842

● FG 850-016
8783843

● FG 856-012
8783846

Bevelled Cylinder



Bevelled Cylinder x 5

● FG 878-012
8783847

● FG 879-012
8783852

● FG 879-012
8783851

● FG 879-014
8783853

Curettage



Curettage x 5

● FG 878K-014
8783848

Diamond Burs

Standard Diamond Burs

HS Diamond Burs

Round



Round

- | | | |
|-----------------------------|-----------------------------|-------------------------|
| ● 520 FG 801-010
9006640 | ● 708 FG 801-014
9006646 | ● FG 801-018
9006707 |
| ● 521 FG 801-012
9006667 | ● 522 FG 801-016
9006689 | ● RA 001-047
9009512 |
| ● 708 FG 801-014
9006642 | ● FG 801-016
9006715 | |

Round with Collar



Round with Collar

- | |
|-----------------------------|
| ● 525 FG 802-010
9006663 |
|-----------------------------|

Inverted Cone



Inverted Cone

- | |
|-----------------------------|
| ● 531 FG 805-010
9006664 |
| ● 531 FG 805-014
9006668 |

Inverted Cone with Collar



Inverted Cone

- | |
|-----------------------------|
| ● FG 807-014
9006662 |
| ● 535 FG 807-016
9006706 |

Round End Taper



Round End Taper

- | |
|-----------------------------|
| ● 501 FG 855-014
9006665 |
|-----------------------------|

Round End Taper



Round End Taper

- | |
|-----------------------------|
| ● 502 FG 850-014
9006666 |
|-----------------------------|

Flat End Taper



Flat End Taper

- | |
|-----------------------------|
| ● 550 FG 845-009
9006674 |
| ● 551 FG 845-012
9006675 |
| ● 552 FG 845-016
9006711 |

Flat End Taper



Flat End Taper

- | |
|--------------------------|
| ● FG 846W-016
9006686 |
|--------------------------|

Flat End Taper



Flat End Taper

- | |
|-----------------------------|
| ● 555 FG 847-012
9006676 |
|-----------------------------|

Flat End Taper



Flat End Taper

- | |
|-----------------------------|
| ● 556 FG 848-016
9006691 |
| ● 556 FG 848-016
9006716 |
| ● FG 848-018
9006682 |

Flat End Taper - Long



Flat End Taper - Long

- | |
|------------------------------|
| ● 557 FG 848L-012
9006677 |
| ● FG 848L-014
9006679 |

Diamond Burs

Standard Diamond Burs

Flat End Cylinder



Flat End Cylinder

● 546 FG 836-012
9006660

● 546 FG 836-014
9006650

Flat End Cylinder



Flat End Cylinder

● FG 835-009
9006656

● 541 FG 835-010
9006659

● FG 835-012
9006648

● 542 FG 835-014
9006671

● 543 FG 835-018
9006710

Round End Cylinder



Round End Cylinder

● 544 FG 838-009
9006672

● 545 FG 838-014
9006673

Flame



Flame

● FG 863-012
9006641

Flame



Flame

● 562 FG 862-012
9006658

● 562 FG 862-014
9006678

● 562 FG 862-016
9006695

● 562 RA 249-016
9009513

Flame



Flame

● FG 860-016
9006717

Needle



Needle

● FG 859-014
9006649

● 558 FG 859-018
9006712

Donut



Donut

● FG 909-037
9006719

Xmas Tree



Xmas Tree

● FG 852-012
9006680

Pear



Pear

● 526 FG 830-012
9006651

● 527 FG 830-014
9006645

Wheel



Wheel

● 512 FG 815-018
9006708

Egg



Egg

● FG 379-018
9006703

Diabolo



Diabolo

● 536 FG 813-012
9006669

Diabolo - Long



Diabolo - Long

● 538 FG 813L-016
9006709

Diamond Burs

Disposable Diamond Burs - Sterile

HS Sterile Single Use Diamond Burs

Round



- Round x 25
- 521 FG 801-010
9006728
 - 522 FG 801-016
9006729
 - 801 FG 801-023
9006730

Round End Taper



- Round End Taper x 25
- 502 FG 850-014
9006733

Round End Taper - Long



- Round End Taper - Long x 25
- FG 856L- 016
9006754
 - FG 856L-018
9006737

Round End Taper



- Round End Taper x 25
- | | |
|-----------------------------|-------------------------|
| ● 504 FG 856-012
9006786 | ● FG 856-016
9006735 |
| ● FG 856-014
9006734 | ● FG 856-018
9006736 |
| ● FG 856-014
9006786 | |

Flat End Taper



- Flat End Taper x 25
- FG 847-014
9006791
 - FG 847-016
9006747
 - FG 847-018
9006752

Flat End Taper



- Flat End Taper x 25
- FG 848-010
9006771
 - FG 848-016
9006788

Flat End Taper - Modified



- Flat End Taper - Modified x 25
- FG 847KR-016
9006756

Flat End Cylinder



- Flat End Cylinder x 25
- 541 FG 835-010
9006732
 - 541 FG 835-010
9006751
 - 542 FG 835-014
9006749

Bevelled Cylinder



- Bevelled Cylinder x 25
- FG 885-012
9006753
 - FG 885-016
9006750

Football



- Football x 25
- FG 257-023
9882123

Football



- Football x 25
- FG 368-018
9006743
 - FG 368-023
9006726
 - 651 FG 368-023
9006770
 - FG 368-023
9006775

Football



- Football x 25
- FG 379-023
9006727
 - FG 379-023
9006792

Diamond Burs

Disposable Diamond Burs - Sterile

Flame



- Flame x 25
- FG 862-012
9006740
 - 562 FG 862-014
9006741
 - FG 862-016
9006742

Needle



- Needle x 25
- FG 859-014
9006739
 - FG 859-014
9006781

Pear



- Pear x 25
- 526 FG 830-012
9006748

Pointed Cone



- Pointed Cone x 25
- FG 858-014
9006738

Inverted Cone



- Inverted Cone x 25
- FG 805-012
9006731
 - FG 805-016
9006724

Donut



- Donut x 25
- FG 909-040
9006745

Barrel



- Barrel x 25
- FG 811-033
9006725

Torpedo



- Torpedo x 25
- FG 878-012
9006757

Interproximal



- Interproximal x 25
- FG 392-016
9006785

Disposable Diamond Burs - Non-Sterile

HS Single Use Diamond Burs

- Non-sterile single-use diamond bur.

Football



- Football x 5
- FG 257 257-018
9008309

HS Steel Burs - Bulk Pack

Round



Round x 25

RA Size 1 008
9008367

RA Size 3 012
9008369

RA Size 5 016
9008371

RA Size 7 021
9008373

RA Size 9 025
9008377

RA Size 2 010
9008368

RA Size 4 014
9008370

RA Size 6 018
9008372

RA Size 8 023
9008374

RA Size 10 027
9008378

HS Steel Burs

Round



Round x 6

RA Size ½ 006
9008777

RA Size 5 016
9008773

RA Size 10 027
9008780

HP Size 11 031
9008798

RA Size 1 008
9008778

RA Size 6 018
9008774

RA Size 11 031
9008804

HP Size 12 033
9008799

RA Size 2 010
9008770

RA Size 7 021
9008775

RA Size 12 033
9008805

RA Size 3 012
9008771

RA Size 8 023
9008776

RA Size 4 014
9008772

RA Size 9 025
9008779

Inverted Cone



Inverted Cone x 6

RA Size ½ 006
9008781

RA Size 3 012
9008784

HP Size 2 010
9008786

RA Size 1 008
9008782

RA Size 4 014
9008785

RA Size 2 010
9008783

Flat Fissure Cross Cut



Flat Fissure Cross Cut x 6

RA Size 1 008
9008800

RA Size 3 012
9008802

RA Size 4 014
9008803

Cylinder



Cylinder x 6

RA Size 2 010
9008792

HP Size 1 008
9794084

HP Size 4 014
9794088

HP Size 2 010
9794085

HP Size 8 023
9794087

HP Size 3 012
9794086

Taper



Tapered x 6

RA Size 2 010
9008787

HP Size 2 010
9008788

HP Size 3 012
9008789

HP Size 5 016
9008790

Tungsten Carbide Burs

HS Tungsten Carbide Burs

Round



Round x 10	
FG Size ½ 005	FG Size 4 014
1007205	1004535
FG Size 1 008	FG Size 5 016
1004907	9888002
FG Size 2 010	FG Size 6 018
1000288	1003220
FG Size 3 012	FG Size 8 023
9888001	1006131

Round



Round x 10	
RA Size ½ 006	RA Size 4 014
1006373	1001847
RA Size 1 008	RA Size 5 016
1006319	1000664
RA Size 2 010	RA Size 6 018
1006411	1008765
RA Size 3 012	RA Size 8 023
9888003	1001401

Round



Round x 5	
HP Size 5 016	
5700261	
HP Size 6 018	
5700262	
HP Size 8 023	
5700263	

Pear



Pear x 10	
FG 330 Size 1 008	
1008426	
FG 331 Size 2 010	
1001120	
FG 332 Size 3 012	
1005616	

Flat Fissure Plain Cut



Flat Fissure Plain Cut x 10	
FG 57 Sire 2 021	
1002579	
FG 58 Size 3 012	
1000281	

Flat Fissure Cross Cut



Flat Fissure Cross Cut x 10	
FG 556 Size 1 008	FG 558 Size 3 012
1000898	1003104
FG 557 Size 2 010	FG 558L Size 3 012
1003307	1008157
FG 557L Size 2 010	
1009947	

Flat Fissure Cross Cut



Flat Fissure Cross Cut x 10	
RA 557 Size 2 010	
1003795	

Domed Fissure Plain Cut



Domed Fissure Plain Cut x 10	
FG 1157 Size 2 010	
1003098	
FG 1158 Size 3 012	
1005104	

Domed Fissure Cross Cut



Domed Fissure Cross Cut x 10	
FG 1557 Size 2 010	
1004372	
FG 1558 Size 3 012	
1009613	

Tapered Fissure Cross Cut

Tapered Fissure Cross Cut x 10



FG 699 Size 1 008
1009613

FG 700 Size 2 010
1007926

FG 700L Size 2 010
1003046

FG 701 Size 3 012
10045463

FG 701L Size 3 012
1007228

FG 702 Size 5 016
1009464

Tapered Fissure Cross Cut

Tapered Fissure Cross Cut x 5



HP 701 Size 3 012
5700266

Tapered Fissure Plain Cut

Tapered Fissure Plain Cut x 10



FG 169 Size 1 008
1001864

FG 171L Size 3 012
1002598

FG 170L Size 8 023
1003463

Amalgam Preparation

Amalgam Preparation x 10



FG 245 Size 8 023
1008848

Inverted Cone

Inverted Cone x 10



FG 34 Size 1 008
1008454

FG 37 Size 4 014
1009299

End Cutting

End Cutting x 10



FG 957 Size 2 010
1008216

FG 958 Size 3 012
1001158

Short Shank



Inverted Cone x 10

FG 35 Size 2 010
1002805

Tapered Fissure Plain Cut x 10

FG 170 Size 2 010
1009331

Tapered Fissure Cross Cut x 10

FG 700 Size 2 010
1002377

Flat Fissure Cross Cut x 10

FG 57 Size 2 010
1002975

FG 557 Size 2 010
1000273

FG 558 Size 3 012
1006142

Surgical



Surgical Round x 5

FG Size 2 010
1001344

FG Size 4 014
1009992

FG Size 6 018
1002734

FG Size 8 023
1006142

Surgical Fissure Plain Cut x 5

FG Size 3 012 058
1008404

Surgical Fissure Cross Cut x 5

FG Size 5 010 557
1008502

FG Size 3 012 558
1000282

Surgical Tapered Fissure Cross Cut x 5

FG Size 5 016 702
1000757

Crown Cutters

- Ideal for rapid reduction of all dental materials such as amalgam, porcelain, metals, and tooth structure.
- Sharp blade design and multiple crosscuts provide the perfect combination for crown removal and endodontic access.
- Single-piece construction eliminates bounce, chatter and vibration.



FG 331
9004176



FG 557
9004179



FG 245
9004178



FG 330
9004177



FG 1557
9004175



FG 1558
9004174

- One-piece burs made from 100% fine grade tungsten carbide, offering strength over many competitive products that only contain carbide heads.
- Cuts more effectively with greater blade contact on the tooth surface and reduced risk of chipping or breakage.



FG ½ 006
9004508



FG 33½ 006
9004511



FG 245 009
9004515



FG 330 008
9004516



FG 56 008
9004513



FG 556 008
9004518



FG 700 010
9004521



FG 1157 010
9004522



FG 1557 010
9004524

Tungsten Carbide Burs

Orthodontic Debonding Burs

HS Orthodontic Debonding Burs

- For debonding and adhesive removal.
- Controlled debonding without scratching or abraiding the enamel.

	8 Blade Burs x 5			
	FG 1172 016 1125568	FG 7406 018 1125565	RA 7006 018 1125573	FG 7642 010 1125574
	12 Blade Burs x 5			
	FG 7803 012 1125572	FG 7901 009 1125564	FG 7902 010 1125570	FG OF1 023 1125566
	18 Blade Burs x 5			
	FG 218 016 1125569	FG 118 S 016 1125561	FG 118 L 016 1125562	

Carbide Trimming & Finishing Burs

HS Gold Aesthetic Finishing Burs Safe-End EF Series

- 10-blade gold coloured bur.
- For optimum shaping and finishing of all restorations, particularly composite veneers.

	Safe End x 5
	FG EF3 008 9004302
	FG EF4 010 9004303
	FG EF6 014 9004301
	FG EF9 014 9004300

HS Gold Aesthetic Finishing Burs OF Series

- 12-blade gold coloured bur.
- For optimum shaping and finishing of all restorations, particularly composite veneers.

	OF Series x 5
	FG OF1 379-023 9004387
	FG OF2 247-009 9004388

Tungsten Carbide Burs

Carbide Trimming & Finishing Burs

HS Gold Trimming & Finishing Burs - 12 Bladed

- 12-bladed trimming and finishing burs with tungsten carbide head on a stainless steel shank.
- Easy and effective finishing of composite, amalgam and enamel.
- Various shapes for occlusal, lingual, labial and interproximal finishing.



Round x 5
FG 7006 018
9004355

FG 7008 023
9004356



Flame x 5
FG 7104 014
9004357

FG 7106 018
9004358

FG 7108 023
9004359



Flame Short x 5
FG 274 016
9004386



Cone x 5
FG 7606 018
9004366



Taper Long x 5
FG 7204 014
9004360

FG 7205 016
9004361

FG 7206 018
9004362



Egg x 5
FG 7404 014
9004363

FG 7406 018
9004364

FG 7408 023
9004365



Round End Taper x 5
FG 7642 010
9004384

FG 7664 014
9004385



Bullet x 5
FG 7801 009
9004371

FG 7802 010
9004372

FG 7803 012
9004373



Needle x 5
FG 7901 009
9004374

FG 7902 010
9004375

FG 7903 012
9004376



Taper Point x 5
FG 7610 008
9004367

FG 7611 010
9004368

FG 7612 014
9004369

FG 7613 014
9004370

HS Gold Fine Finishing Burs Safe-End EF Series

- 30-blade gold coloured bur.
- For final finishing of all dental materials.
- Ideal for use on ceramics.



FG 9006 x 5
9004310



FG 9714 x 5
9004308



FG 9214 x 5
9004311



FG 9903 x 5
9004309



FG 9406 x 5
9004306



FG 9904 x 5
9004307

Endo Access Burs

HS Endo Safe-End Access Burs

- Long, safe-ended access bur for opening the pulp chamber.
- Round, non-cutting tip prevents perforating both the pulp chamber floor and root canal walls.
- Removes debris cleanly and without vibration.



Endo Access x 5
Length 21mm
9007038

Tungsten Carbide - Surgical Burs

HS Tungsten Carbide Surgical Burs

- Two-piece bur with 100% carbide head on a stainless steel shank.



Surgical Round x 5
FG Size 2 010
1001344
FG Size 4 014
1009992
FG Size 6 018
1002734
FG Size 8 023
1006142



Surgical Fissure Plain Cut x 5
058 FG Size 3 012
1008404



Surgical Fissure Cross Cut x 5
557 FG Size 5 010
1008502
558 FG Size 3 012
1000282



Surgical Tapered Fissure Cross Cut x 5
702 FG Size 5 016
1000757

HS Tungsten Carbide Surgical Burs

- Surgical burs with 6 helical cutting edges.
- For sectioning impacted wisdom teeth, separating roots, and cutting horizontally at root level.



FG 151 23mm
9006299



FG 151L 28mm
9006300

Tungsten Carbide - Sterile Oral Surgery Burs

HS Oral Surgery Burs - Sterile

- Sterile single-use oral surgery burs.
- Individually wrapped.
- Sterilisation by gamma radiation.



Tapered Fissure
HP 702
9008379
RA 702
9008380




Round
HP 8
9008375
RA 8
9008376

Bur Brushes

Bur Brush Flat

- Nylon bur cleaning brush.
- HTM 01-05 compliant.

20 x 70mm
9793155



Hand Instrument Air Collection

From A Name You Can Trust.

Our Stainless Steel Air collection of single and double ended instruments are light as air with an average weight of 12 grams.

You can trust that Henry Schein Brand Hand Instruments provide unsurpassed quality—backed by the Henry Schein Seal of Excellence, your 100% guarantee of satisfaction.

Our Instrument Sets are especially designed for you:

**Air Diagnostic Instrument Set #9
900-8401**

- 1 Explorer #9
- 1 Mouth mirror handle with mirror
- 1 Colleague Tweezer

**Air Diagnostic Instrument Set #23
900-8184**

- 1 Probe SE #23
- 1 Mouth mirror handle with mirror
- 1 Colleague plier 317, serrated

See pages 61-70 for our selection of Instrument Products.



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!

Endodontics



Drills & Reamers	96
Barbed Broaches	96
K-Files	96
Hand Reamers	97
Hedstroem Files	97
NiTi Files	98
Finger Pluggers	98
Paste Fillers	98
GP Points	99
Paper Points	99
Greater Taper GP Points	100
Root Canal Sealers	100
Medicaments	100
Rubber Dam & Clamps	101
Endo Accessories	102

Drills & Reamers

HS Gates Glidden Drills

- Stainless steel blunt-ended tip to prevent perforation and aid better tactile sensitivity.
- For use with slow handpiece to widen the canal orifice.
- 32mm length.

Pack of 6

	32mm
Size 1	9880836
Size 2	9880837
Size 3	9880838
Size 4	9880839
Size 5	9880840
Size 6	9880841
Sizes 1-6	9880842

HS Peeso Reamers RA

- For enlargement of cervical part of root canal, preparation of root canal entrance, and pin-hole preparations.
- Stainless steel.
- 19mm working length.

Pack of 6

	32mm
Size 1	9880843
Size 2	9880844
Size 3	9880845
Size 4	9880846
Size 5	9880847
Size 6	9880848
Sizes 1-6	9880849

Barbed Broaches

HS Barbed Broaches - Sterile

- Sterile Barbed Broaches.

Pack of 10

	21mm
0 - XXXXFine	9008246
1 - XXXFine	9008247
2 - XXFine	9008248
3 - XFine	9008249
4 - Fine	9008250
5 - Medium	9008251
6 - Coarse	9008252

K-Files

Endoflex K-Files

- ISO standard stainless steel files with colour coded handle.

Pack of 6



21mm	25mm	31mm
6 1007112	6 1000512	15-40 1028480
45 1009397	45 1002032	45-80 1021620
8 1008298	8 1008509	
50 1009595	50 1001180	
10 1007076	10 1000227	
55 1006118	55 1008595	
15 1008332	15 1001109	
60 1007702	60 1008636	
20 1008457	20 1000386	
70 1007906	70 1008859	
25 1009024	25 1000571	
80 1007979	80 1008945	
30 1009072	30 1000616	
15-40 1009709	15-40 1001215	
35 1009254	35 1000689	
45-80 1004863	45-80 1005087	
40 1009336	40 1002878	

Hand Reamers

HS Endoflex Reamers

- ISO standard stainless steel files with colour coded handle.

Pack of 6



21mm	25mm
 10 1001317	 10 1002614
 8 1002426	 15 1002034
 15 1005631	 20 1002895
 20 1001862	 25 1003731
 25 1004759	 30 1003128
 30 1007587	 35 1001974
 35 1003811	 8 1002671
 40 1008507	 40 1003236
 45 1007245	 45 1001092
 50 1002481	 50 1003459
 55 1008353	 55 1009935
 60 1002202	 60 1009592
 70 1007242	 70 1002996
 80 1006682	 80 1001097
15-40 1002541	15-40 1001212
45-80 1008887	45-80 1003321

Hedstroem Files

HS Endoflex Hedstroem Files

- ISO standard stainless steel files with colour coded handle.

Pack of 6



21mm	25mm
 10 1004237	 10 1004240
 15 1003883	 15 1005139
 20 1004074	 20 1005318
 25 1000765	 25 1005362
 30 1004196	 30 1005429
 35 1004368	 35 1004085
 40 1004055	 40 1005678
 45 1000591	 45 1005716
 50 1004640	 50 1005767
 55 1003570	 55 1004210
 60 1005683	 60 1004290
 70 1009388	 70 1004453
 80 1003923	 80 1003784
15-40 1004722	15-40 1000069
45-80 1004137	45-80 1004596

NiTi Files

HS NT Files

- Made from nickel titanium alloy.
- Extremely flexible to follow the canal outline.
- Excellent memory.

Pack of 6



21mm	25mm
○ 15 9001188	○ 15 9001202
● 20 9001189	● 20 9001203
● 25 9001190	● 25 9001204
● 30 9001191	● 30 9001205
● 35 9001192	● 35 9001206
● 40 9001193	● 40 9001207
○ 45 9001194	○ 45 9001208
● 50 9001195	● 50 9001209
● 55 9001196	● 55 9001210
● 60 9001197	● 60 9001211
15-40 9001198	15-40 9001212

Finger Pluggers

HS Finger Pluggers

- Stainless steel colour coded plastic handles.

Pack of 4



21mm	25mm
● A 1008766	● A 1000995
● B 1009435	● B 1000379
● C 1003320	● C 1002022
● D 1009690	● D 1002187
Assorted 1003299	Assorted 1002619

Paste Fillers

Paste Filler with Spring

- Closed spring on shank absorbs pulse movement and gives breakage resistance.

Pack of 4



21mm	25mm
● 25 1 9882883	● 25 1 9880690
● 30 2 9882884	● 30 2 9880691
● 35 3 9882885	● 35 3 9880692
● 40 4 9882886	● 40 4 9880693
1-4 9882887	1-4 9880694

GP Points

HS Maxima® Gutta Percha Points

- Machine rolled gutta percha points.
- Flexible and pliable, yet firm enough to be inserted into the root canal.
- Packed in plastic drawer box.
- ISO sizes and colour coded tips.



Pack of 120

○ 15	● 35	15-40
9001230	9001234	9001238
● 20	● 40	45-80
9001231	9001235	9001239
● 25	○ 45	
9001232	9001236	
● 30	● 50	
9001233	9001237	

HS Gutta Percha Points

- GP points pack in vials.
- Colour coded tips for easy identification of ISO sizes.



Pack of 6 Vials

○ 15	● 35	45-80
9008974	9008978	9008993
● 20	● 40	90-140
9008975	9008979	9008994
● 25	○ 45	
9008976	9008980	
● 30	15-40	
9008977	9008992	

HS GP Auxiliary Points

- Non-ISO sized, machine-rolled auxiliary points.
- Flexible and pliable.
- Packed in a plastic drawer box.

Pack of 100

XXFine
9001241
XFine
9001242
Fine
9001243

Medium
9001244

Large
9001244



Paper Points

HS Maxima® Paper Points

- Colour coded tips packed in a plastic drawer.
- Sterile.




Pack of 200

○ 15	● 35	15-40
9001213	9001217	9001221
● 20	● 40	45-80
9001214	9001218	9001222
● 25	○ 45	90-140
9001215	9001219	9001223
● 30	● 50	
9001216	9001220	

HS Paper Points Vials

- Paper points in vials.
- Colour coded.
- Sterile.



Pack of 6

○ 15	● 35	● 55	● 100
9009013	9009017	9009021	9009026
● 20	● 40	● 60	15-40
9009014	9009018	9009022	9009031
● 25	○ 45	● 70	45-80
9009015	9009019	9009023	9009032
● 30	● 50	● 80	
9009016	9009020	9009024	

Greater Taper GP Points

HS Greater Taper Gutta Percha Points

- Tapered gutta percha points for use with matching 0.04 and 0.06 taper files, in conjunction with single-cone obturation techniques.
- Colour coded.
- Radiopaque.

Pack of 100

0.04

15

9008995

20

9008996

25

9008997

30

9008998

35

9008999

40

9009000

45

9009001

50

9009002

15-40

9009003

0.06

15

9009004

20

9009005

25

9009006

30

9009007

35

9009008

40

9009009

45

9009010

50

9009011

15-40

9009012



Medicaments



HS EDTA

- EDTA at 19% in a syringe format.
- For removal of inorganic substances by chelation.

EDTA Gel Syringe 7g x 2 plus 10 tips
9003700

HS Sodium Hypochlorite

- Aqueous sodium hypochlorite solution not exceeding 3.0%.
- Used for the irrigation, debridement and disinfection of root canals during and after instrumentation.


Solution 500ml
9008305



HS Utility Syringe

- Non-sterile, multipurpose curved tip syringe for endo irrigation or intraoral impression material delivery.

Curved Tip Syringe x 50
12cc
9006551



Root Canal Sealers



HS MTA Bioceramic Root Canal Sealer

- Bioceramic root canal sealer on an MTA base (mineral trioxide aggregate).
- Biocompatible for tissue recovery without reactions.
- High radiopacity.
- Aids in the quick recovery of bone and cementum formation.
- Working time 23 minutes, setting time 130 minutes.

Cartridge 9g plus Tips x 3
9880705

Base 18g and Catalyst 12g
9880706

HS Sterile Luer Lock Syringe

- 3-part disposable syringes.
- Smooth flow, tight, synthetic rubber piston.
- Latex-free, pyrogen-free and PVC-free.
- Single use, individually wrapped.
- Sterile.



Sterile Syringes x 100

3ml
9003759

5ml
9003760

10ml
9003761

HS Endodontic Blunt Needle

- Rounded, closed-end needle tips with side-port delivery.
- Tips eliminate hydraulic pressure and ensure a safer delivery of solutions when irrigating root canals and periodontal pockets.

Needles x 100

Blue 23G	Yellow 27G	Violet 30G
9006549	9006550	9006548



HS Endodontic Irrigating Needle

- Bendable side-vent irrigating tips designed for irrigating and cleansing post holes, canals, pockets, and fistulas.
- 100% latex free. Made from polypropylene and high grade stainless steel

Needles x 100

Blue 23G	Yellow 27G
9006546	9006545

Rubber Dam & Clamps

HS Rubber Dam Clamps

- Rubber dam clamps made from high-quality, stainless steel.
- Strong retention and excellent patient comfort.
- Matt finish to enable easier visibility.
- The assorted set comes complete with an autoclavable stand.

Assortment Kit

Autoclavable Stand with 11 Clamps, one each of 2, 2A, 7, 8A, 9, 12A, 13A, 14, 14A, W3, W8A
9004560



Size	Size 13A	
Wings	Clamps Winged	
Tooth Type	Lower Molar Serrated	
	9004554	
Size	Size 14	
Wings	Clamps Winged	
Tooth Type	Partially Erupted Molars	
	9004551	
Size	Size 12A	
Wings	Clamps Winged	
Tooth Type	Lower Molar Serrated	
	9004553	
Size	Size 2A	
Wings	Clamps Winged	
Tooth Type	Large Premolars	
	9004546	
Size	Size 14A	
Wings	Clamps Winged	
Tooth Type	Premolars	
	9004552	
Size	Size W8A	
Wings	Clamps Without Wings	
Tooth Type	Molars - Root	
	9004550	

Size	Size W3	
Wings	Clamps Without Wings	
Tooth Type	Lower Molars	
	9004549	
Size	Size 9	
Wings	Clamps Winged	
Tooth Type	Anterior	
	9004547	
Size	Size 2	
Wings	Clamps Winged	
Tooth Type	Premolars	
	9004545	
Size	Size 8A	
Wings	Clamps Winged	
Tooth Type	Molars - Root	
	9004555	
Size	Size 7	
Wings	Clamps Winged	
Tooth Type	Lower Molars	
	9004548	

Rubber Dam & Clamps

Dental Dam

HS Latex Rubber Dam

- Rubber dam made from natural latex.
- Powder-free.
- Mint flavoured.
- Green in colour.
- Size: 15cm x 15cm.



Rubber Dam x 36

Thin
9002650

Medium
9002651

Heavy
9002652

HS Non Latex Dam

- Latex-free dam.
- Powder-free.
- Blue in colour.
- Benefits from a stretchy elastic feel making it easier to use.
- Size: 15cm x 15cm.



Non-Latex Dam Blue
9003055

Dam Frames & Instruments



HS Rubber Dam Frame

- Stainless steel rubber dam frame.
- U shaped.

Frame
9006436

HS Rubber Dam Punch Ainsworth

- Stainless steel rubber dam punch for all types of rubber dam.

Ainsworth Type
1043375



Endo Accessories

HS Endo Cold Spray

- Chilling spray (medical device) suitable for icing pellets of cotton wool - for testing vitality of teeth.
- Also used for quick cooling of impression materials.
- Replacement for cryogenic cold spray.



Mint 200ml
9001383

Orange 200ml
9001384




HS Endo Sponges

- Foam sponges for round endo organisers.
- 5cm diameter.
- 0.5cm and 1cm thickness.
- The foam pads are autoclavable up to 135°C.

Sponges 0.5 x 5cm x 50
9002696

HS Endo Organiser Ring

- Endo organiser ring with foam insert.
- Fits neatly on finger.
- Autoclavable at 121° only.
- Integral ruler.



Endo Ring plus 50 Foam Refills

Ring - Light Green
9008384

Accessories

Foam Refills
Lavender x 25 and Grey x 25
9008385

HS Rubber Dam Clamp Forceps

- Stainless steel rubber dam clamp forceps in Brewer or lightweight pattern

Lightweight
1043536

Brewer Pattern
1043374



Small Equipment



Alginate Mixing Machines	104
Micromotors	104
Curing Lights	105
Handpiece Maintenance	106
X-Ray Viewers	106

Alginate Mixing Machines

HS Alginate Mixer

- Turns alginate into a bubble-free mix within seconds, ensuring a high-quality impression.
- Digital control panel timer can be set between 1 and 16 seconds.
- The mixer operates at a maximum of 2,900rpm.
- Includes multiplug, 4 mixing cups, 2 spatulas in different sizes, and a water dosing bottle.
- Dimensions: (W)21 x (H)29 x (D)25cm.
- Weight: 19kg.

Alginate Mixer
9003892



Micromotors

HS Micromotor L50

- Brushless micromotor with autocruise and forward/reverse button.
- Adjustable to 50,000rpm.
- Operated using foot or knee control.
- Weight: 3kg.
- Max torque: 7.8N/cm.

Micromotor L50
9002693



HS LED 1200 Cordless Curing Light

- Cordless curing light with an output of over 1200mW/cm².
- Digital control panel with audible signal for the last 3 seconds.
- High performance diode mechanism can polymerise an average 3mm thickness of material within 10 seconds on almost every resin-based composite in the market.
- Includes 1 x fibre-optic guide (8mm), 2 small orange tips and orange shield.
- Operates at a wavelength between 450 and 470nm.
- Long life 5W single blue LED bulb.
- Dimensions: 25cm x 3.5cm x 2.5cm
Charger unit: 11.5cm x 8cm x 7.5cm
Weight: 154g.

LED 1200
9006369



HS LED 1500 Cordless Curing Light

- Lightweight cordless curing light with output of over 1500mW/cm².
- Simple to operate control panel with a digital display and acoustic signal for the last 3 seconds.
- 3 different time adjustable modes: full mode, ramp mode, pulse mode.
- Including a fibre-optic 8mm guide, 2 small orange tips and orange shield.
- Operates at a wavelength between 450 and 470nm.
- Dimensions: 21 cm x 4.3cm x 4.3cm.
- Charger dimension: 19.2cm x 9.7cm x 6.5cm.
- Weight: 175g.

LED 1500
9006370



Handpiece Maintenance



HS Maxima® Oil Spray

- Universal spray for handpieces.
- Degreases, cleans and lubricates.

500ml
9001570



HS S-Spray

- Universal spray for the cleaning and care of Sirona turbines and handpieces.

250ml
9003292

X-Ray Viewers


HS XRV 30

- XRV 30 is an X-ray viewer with a filtered plexiglass screen.
- Scratch-resistant, coated stainless steel frame.
- Wall mountable for tabletop use.
- Screen dimensions - 300 x 165mm for viewing films up to 300 x 150mm.
- Uses 2 x 8W bulbs.

X-Ray Viewer
9003256



X-Ray & Imaging



Intraoral Films	108
Extraoral Films	108
X-Ray Protection	108
X-Ray Chemicals	109
X-Ray Accessories	110
Barrier Envelopes	110

Intraoral Films

HS DX-54 Child Periapical Film

- Child size periapical film for examination of the whole tooth and surrounding structures.
- D Speed.

Flexisoft x 100

22 x 35mm
9008387

HS DX-58 Adult Periapical Film

- Flexisoft D Speed Film for excellent image quality.

Flexisoft x 150

31 x 41mm
9008389

HS FX-58 Adult Periapical Film

- Flexisoft periapical X-ray film with high contrast.
- Easy-to-use packets.
- F Speed.

Flexisoft x 150

31 x 41mm
9008392



Extraoral Films

HS Pan-G Extraoral Film

- Panoramic films for green sensitive intensifying screens.
- Speed class depends on the intensifying screen used.
- As a guideline:-
Pan G + Fine screens = 100,
Pan G + Medium screens = 200,
Pan G + Regular = 400,
Pan G + Fast = 600.




Film x 100

12 x 30cm
9008173

15 x 30cm
9008174

18 x 24cm
9008175

X-Ray Protection



HS Thyroid Collar

- Thyroid collar to protect the thyroid glands.
- Shape provides a high level of comfort and protection.
- One size: adjusted by velcro band.

Blue
9881840



HS Concentrated Developer/Fixer

- Concentrated solution for manual X-ray processing.
- Dilute 500ml fixer or developer with 3 x 500ml water.
- Makes 2 litres of in use solution.

Developer 500ml
9002683

Fixer 500ml
9002684

HS Radiomat XR Developer and Fixer

- Radiomat XR concentrated developer and fixer.
- Suitable for automatic processors, for example Dürer Automat or Velopex.
- Each bottle contains 1.5 litres and makes 5 litres of solution.
- Use in appropriate tanks in processor.
- DO NOT MIX DEVELOPER AND FIXER. Follow the machine manufacturer's instructions for processing.
- Renew solutions in machine completely every three months for optimum function when using automatic regeneration and more frequently for low usage machines.
- For manual regeneration - top up twice weekly with 0.5 litres of the appropriate solution. This is a general rule for processing 7 films per day.

RadioMAT-XR 4 Pack

Developer 1.5L x 2 and Fixer 1.5L x 2
9002685

X-Ray Accessories

HS X-Ray Cleaning Films

- Pre-exposed coated cleaning films for automatic processors.
- Run the cleansing film through the system to clean rollers.

18 x 24cm x 100
9008390



Barrier Envelopes

HS Barrier Envelopes

- Phosphor plate envelopes to protect plates from exposure to light.
- Works with all brands of phosphor plate readers.
- No scissors required to open up.
- Easy tear-away seam.
- Easy-peel glue strip cover.
- Latex-free.

Black/Clear

Size 0 x 100
9009626

Size 1 x 100
9009627

Size 2 x 300
9009628



Laboratory Products



CAD/CAM Materials	112
Casting	112
Crown & Bridge	114
Duplication Material	114
Lab Putty	114
Model Plaster & Stone	115
Prosthetics	115
Finishing and Polishing	116
General Lab Items	117
Instruments	118
Rotary Instruments	121
Lab Small Equipment	121

CAD/CAM Materials

CAD/CAM Accessories

HS CAD/CAM Lab Finishing and Polishing

- CAD/CAM finishing and polishing kit HP designed for lab use.
- Diamond polishers and trimmers interspersed with diamond grit for adapting and correcting CAD/CAM restorations.
- Trimmers HP: speed 5000-10000rpm.
- Polishers HP: speed 7000-12000rpm.



Lab Finishing and Polishing Kit HP

Trimming Wheel and Cylinder in Coarse and Medium; and Polishing Wheel, Point and Unmounted Eclipse in Medium and Fine
9009510

Finishing and Trimming Refill HP

Coarse	Medium
Wheel 2661 9003351	Wheel 2660 9003353
Cylinder 2663 9003352	Cylinder 2662 9003354
Finishing and Polishing Refill HP	
Medium Unmounted	Fine Unmounted
Eclipse 2665 9003355	Eclipse 2669 9003358
Finishing and Polishing Refill HP	
Medium	Medium
Point 2668 9003356	Wheel 2666 9003357
Fine	Fine
Point 2664 9003359	Wheel 2667 9003360

HS CAD/CAM Chairside Finishing and Polishing

- CAD/CAM finishing and polishing RA instruments in fine, medium and coarse grits for quick and easy adjustment of CAD/CAM restorations, inlays, onlays, veneers and crowns.
- Diamond impregnated grinders and polishers that can be used in a wet or dry application to give a high lustre finish.
- Chairside use.
- Coarse RA: speed 5000-10000rpm.
- Medium, Fine RA: speed 7000-10000rpm.



Chairside Finishing and Polishing Kit RA

Trimming Wheel and Cylinder in Coarse and Medium; and Polishing Cup, Point and Unmounted Eclipse in Medium and Fine
9009511

Finishing and Trimming Refill RA

Coarse	Medium
Wheel 2661 9003500	Wheel 2660 9003502
Cylinder 2663 9003501	Cylinder 2662 9003503
Finishing and Polishing Refill RA	
Medium Unmounted	Fine Unmounted
Eclipse 2665 9003504	Eclipse 2669 9003507
Finishing and Polishing Refill RA	
Medium	Medium
Point 2668 9003505	Cup 2666 9003506
Fine	Fine
Point 2664 9003508	Cup 2667 9003509

Casting Accessories

HS Casting Rings

- Made of stainless steel.
- Height 55mm.

No 1 Ring 30mm
9002200

No 3 Ring 48mm
9002201

No 6 Ring 65mm
9002202



Casting Accessories



HS Profile Wax Round

- Tension-free wax wire.
- Malleable.
- Prevents accidental deformation prior to casting.
- High stability allows it to be waxed onto any object.
- Blue.

0.8mm 25g
9002238

1.0mm 25g
9002239



HS Waxwire

- Malleable, tension-free wax wire.
- Prevents accidental deformation prior to casting.
- High stability allows it to be waxed onto any object.

Hard Blue

3mm 250g
9002223
3.5mm 250g
9002224
4mm 250g
9002225
5mm 250g
9002226

Standard Green

3mm 250g
9002227
3.5mm 250g
9002228
4mm 250g
9002229
5mm 250g
9002230



HS Sprue Formers

- A range of formers for use with casting rings.

No 1 30mm
9002204

No 6 65mm
9002206

No 3 48mm
9002205

HS Electrolyt Solution

- Ready-to-use solution.
- Special additives for quick and effective electro-polishing.
- For use with copper, stainless steel or titanium electrodes.



1L
9009666

HS CAD/CAM Scan Spray

- CAD/CAM Scan Spray is designed to be used with CAD/CAM machines to make it possible to scan preparations.
- For intraoral use and for use on models.

Spray 50ml
9002655



Investment Material


HS Vest

- Investment material for crown and bridgework for all alloys.
- Suitable for precious, non-precious and pressable casting techniques.

Liquid 1L
9002957



Crown & Bridge



HS Dowel Pins

- Range of dowel pins and pins with sleeves in various sizes

Standard Size 1 x 100
9008139

Standard Size 3 x 100
9008141

Extended Tip Size 2 x 100
9008138

HS Die Spacer

- Paint-on acrylic resin.
- Fast drying.
- Different colours provide a different micron thickness.



Bottle 15ml

Transparent 3 Microns
9002408

Red 8 Microns
9002411

Blue 10 Microns
9002409

Silver 13 Microns
9002404

Thinner
9002403

Duplication Material



HS Z Dupe Silicone Flask

- Duplication flask for Z Dupe Silicone.
- Blue cuvette made of Makrolon plastic for silicone dublier technology.
- Adjustable.
- Anchoring by the conical cuff part.

Accessories

Duplication Flask
9008472

Lab Putty



HS Eurosil Max 2 Lab Putty

- Eurosil Max 2 Lab putty is a C-silicone for duplication of denture moulds, impression for copy dentures, insulation of teeth/gypsum in flask.
- Excellent reproduction with good final hardness.
- Superior breaking strength, no deformation in flask.
- Use with Eurosil Max 2 Catalyst.
- Hardness - 85 Shore A.
- Mixing time with gel: 30 sec.
- Total working time: 2:00 min.
- Setting time: 6:00 min.
- Mixing time with liquid: 30 sec.
- Total working time: 1:30 min.
- Setting time: 5:30 min.

Lab Putty

Blue 900ml
9008539

Blue 10kg
9008540

Catalyst Paste

Pink 60ml
9008538

Model Plaster & Stone

Separating Liquids

HS Plaster Separation Liquid

- Separating agent for all plasters.
- Ideal for sectioned models and the split cast technique.

1L
9002435
250ml with Spray Dispenser
9002436



Prosthetics

Abrasives

HS Pumice Mouse

- Pumice stone for cleaning dirty particles on instruments.
- Removes stains, discoloured areas and oily substances from skin.



Pumice Mouse
9002127

HS Pumice Powder

- Quality-grade pumice for pre-polishing of prostheses.

Medium
4.5kg
9003735
15kg
9003736



HS Resin Separation Liquid

- Separating liquid for use with resins/acrylics.

80ml
9002445
1L
9002473



Tray Material

HS Impression Tray Resin LC

- Light-cure tray material for the production of custom impression trays.
- Exhibits excellent characteristics such as handling, curing and quality consistency.
- Easily cured.
- Radiolucent (X-ray invisible).



Upper x 50
Green
9001489

Blue
9001490

Pink
9001491

Prosthetics

Bite Wax

HS Bitewax Bars

- Straight.
- Easily handled when warm.
- Maintains contour when cooled.



Lemon Soft 1850g
9001502

Lemon Soft 490g
9001503

Strawberry Hard 1850g
9001504

Strawberry Hard 490g
9001505

HS Bitewax Rims

- Curved.
- Easily handled when warm.
- Maintains contour when cooled.




Lemon Soft x 100
9001506

Strawberry Hard x 100
9001507

Casting Wax

HS Wax Clasps

- Made of high-quality casting wax.
- High flexibility and good adhesion.
- Easily removed from model.
- Store in a cool place.
- Green.



Molar x 200
9002233

Premolar x 200
9002234

Modelling Wax

HS Modelling Wax

- Standard modelling wax is medium hardness.
- Special modelling wax is medium hardness (slightly softer than standard).
- Very malleable.
- High beeswax content.
- Solidification point: Summer (56°C), Medium (55°C) and Winter (54°C).



Standard

Summer

1.25mm 500g
9001511

1.25mm 2500g
9001510

1.5mm 500g
9001517

1.5mm 2500g
9001516

Medium

1.25mm 500g
9001509

1.25mm 2500g
9001508

1.5mm 500g
9001515

1.5mm 2500g
9001514

Winter

1.25mm 500g
9001513

1.25mm 2500g
9001512

1.5mm 500g
9001519

1.5mm 2500g
9001518

Universal

1.25mm 500g
9002516

1.25mm 2500g
9002515

Special

Summer

1.25mm 500g
9001500

1.25mm 2500g
9001520

1.5mm 500g
9001521

1.5mm 2500g
9001499

Miscellaneous

HS Denture Finish Spearmint

- Reduces residual monomer deposits on the surface of dentures and imparts a pleasantly fresh taste.
- Dentures feel noticeably better in the mouth and are left looking sparkling clean.
- Can also be used on removable orthodontic appliances.



75ml
9002401

Finishing and Polishing

HS Polishing Brushes

Metal/Metal Centre
Goathair Soft ø 49mm
9001789

Plastic Centre 40mm
4-Row Straight ø 80mm
9001708



HS Polishing Pencil

- Black.

Pencil HP Hard x 12
9001701



HS Fabric Polishing Discs

- Polishing discs made from muslin, felt or chamois.



Chamois x 12

Unmounted 22mm
9001824

Felt x 12

HP 22mm
9001809

Muslin x 12

HP 22mm
9001811

Unmounted 22mm
9001825

General Lab Items

Glue

HS Instant Fix

- Cyanoacrylate adhesive.
- Rapid bonding or repair.

2g
9001841

10g
9001840

10g with Brush
9008886



Model Formers

HS Model Formers - Standard

- Made from high-grade rubber elastic material.




Standard

Size 1 Small 75 x 57 x 12mm
9001674

Size 1 Large 80 x 63 x 12mm
9001675

HS Model Formers - Split Cast

- Made from high-grade rubber elastic material.
- Shaped on the bottom to allow the possibility of split cast techniques.



Split Cast Technique

Size 2 Lower Arch 85 x 63 x 16mm
9001678

Size 2 Upper Arch 85 x 63 x 16mm
9001677

Mouth Guard & Splint Material

HS Thermo Guard

- Suitable for model duplication, mouthguards (sports) and nightguards.
- Soft, pliable rubber material.
- Does not bond to acrylic.
- Clear.



ø 120mm x 10

1.0mm 9001716	2.0mm 9001718	4.0mm 9001720
1.5mm 9001717	3.0mm 9001719	

HS Thermo Splint

- Suitable for periodontal splints, bite guards, retainer splints, invisible retainers, provisional crowns and bridges, temporary and fluoride splints.
- Hard elastic material.
- Highly transparent blanks.
- Will bond to acrylic.
- No predrying necessary.
- Clear.



ø 120mm x 10

0.5mm 9001749	1.0mm 9001751	2.0mm 9001753
0.75mm 9001750	1.5mm 9001752	

HS Thermo Bleach Round

- Tray material for model duplications, bracket transfer trays, bleaching trays and moulds for temporary crowns and bridges.
- Tough elastic material.
- Does not bond to acrylic.
- Clear.



ø 120mm x 10

0.5mm 9003738	1.0mm 9003740
------------------	------------------



General Lab Items



HS Coping & Spacer Foils

- Spacer foils for use when fabricating copings.
- Remove before embedding.

Spacer Foil Clear 0.10mm

42mm x 200
9001698

Coping Foil Opaque 0.60mm

42mm x 100
9001699



HS Rebasing Jig 3S

- Precise device with quick action collet nuts.
- Stainless steel metal plates.
- Brass columns and locking nuts.
- Collet nuts always ensure closure is in the correct position.
- Bite rising reduced to a minimum.
- Diameter: 115mm.
- Height: 100mm.

Jig
9001679

HS Silicone Spray

- Silicone spray, all-purpose lubricant and release agent.
- Prevents the bonding of 2 surfaces.

400ml
9002402



HS Bunsen Burner Gas

- Gas refill for various types of laboratory burners.

400ml
9002440



Instruments

Crown & Bridge



HS Specialty Brush For Wax

- For cleaning of wax patterns.

Brush
9002442



HS Brushes for Separation Liquid

- Brushes for the application of separation liquid.

Size 10
9002437

Size 12
9002438

Size 14
9002439

Instruments

Crown & Bridge



HS Ceramist Synthetic Brushes

- Ceramist brushes made from the finest gold-coloured synthetic hair.
- Seamless nickel ferrules.
- Black handle.

Size 00 9001553	Size 3 9001559	Size 7 9001563	Size 14 9001557
Size 0 9001552	Size 4 9001560	Size 8 9001564	
Size 1 9001554	Size 5 9001561	Size 10 9001555	
Size 2 9001558	Size 6 9001562	Size 12 9001556	



HS Ceramic Brushes

- A range of brushes for various laboratory applications.

Glazing Brush 9002406	Big Brush 9002410
Stain Brush 9002416	Ergonomic Handle No 4 9002417



HS Ceramic Instrument

- Fox mixing and condensation.

Ceramic Condenser
9002467




HS Red Sable Brush

- Brushes made from pure Kolinsky red sable hair.
- Seamless nickel ferrules.
- Black handle.

Size 00 9001754	Size 3 9001758	Size 7 9001762
Size 0 9001755	Size 4 9001759	Size 8 9001763
Size 1 9001756	Size 5 9001760	Size 9 9001764
Size 2 9001757	Size 6 9001761	Size 10 9001765

HS Diamond Tipped Forceps

- Stainless steel forceps with diamond tips.
- For the handling of crown and bridge work.



Diamond-Tipped Forceps
9002461



HS Diamond Tipped Jacket Grip Clamp

Curved 9002468	Straight 9002469
--------------------------	----------------------------



HS Ceramics Knife

- Ceramic knife with very thin blade for interdental cutting.

Ceramics Knife
9002463

Instruments

Crown & Bridge



HS Foil Scissors

- Curved

9002471
- Straight

9002472



HS Lab Tweezers

- Curved

9002465
- Straight

9002466

Lab Instruments



HS Crown and Bridge Veneer Instrument

- Colour coded instruments with easy grip handles.
- For use in veneering technique.

- Red

9002249
- Blue

9002250



HS Lab Instrument Set

- 10-piece instrument set in a leather case.

Leather Case
Containing: Iwansson Calipers, Mosquito Forceps Curved, Mosquito Forceps Straight, Foil Scissors Straight, Soldering Tweezers, College Tweezers, Wax Knife Small, Wax Knife Large, LeCron Carver and Zahle Carver

9002462

Wax Instruments

HS Colour Coded Thomas Instruments



Complete Set (Colour Coded)
Nos. 1, 2, 3, 4 and 5 in a Leather Case
9001741

- No 1 Gold

9001736
- No 3 Blue

9001738
- No 4 Red

9001739



HS PK Thomas Carver

- Set**
Carvers x 5 (No 1, No 2, No 3, No 4, No 5)
9002264
- No 1

9002257
- No 2

9002260
- No 3

9002261
- No 4

9002262
- No 5

9002263

Instruments


Wax Instruments



HS Wax Knife

Fahnenstock 13cm 9001650	Lessmann 13cm 9002270	Ermert 9002269
Fahnenstock 17cm 9001651	Lessmann 17cm 9002271	

HS Iwanson Measuring Device




Wax Measuring Gauge 9002266	Metal Measuring Gauge 9002265
--------------------------------	----------------------------------

Lab Small Equipment

HS Tray Cure

- Used for the polymerisation of light-cured tray materials.
- Can cure 2 trays at the same time.
- Digital timer can be set from 6 sec to 99 minutes.
- Automatically shuts down when chamber door is opened.
- 27W/460lm lamp has a lifespan of 5000 hours.

Tray Cure Unit
9002691



Rotary Instruments

Cutting & Separating Discs


HS Separating Discs

- Separating discs for use on precious, non-ferrous and chrome cobalt metals.



22 x 0.6mm x 100
9002885

Lab Equipment Accessories



HS Furnace Tong/Tweezers

Tong 43cm 9002247	Tweezer 30cm 9002459
-------------------------	----------------------------

Rotary Instruments

TC Cutters

HS Tungsten Carbide Cutters

Cross Cut 223

- Cross Cut Extra Coarse Black for plasters and acrylics.



Fig 79

● HP 194 070
9001937

Cross Cut 225

- Cross Cut Extra Coarse Black for plasters and acrylics.



Fig 251

● HP 274 060
9001936

Cross Cut 220

- Cross Cut Coarse Green for plasters and acrylics.



Fig 79

● HP 194 060
9001992



Fig 251

● HP 274 060
9001981

Cross Cut 190

- Cross Cut Standard Blue for plasters, acrylics, metal alloys and precious metals.



Fig 72

● HP 144 060
9001856



Fig 73

● HP 277 014
9001969

● HP 277 023
9001966



Fig 76

● HP 254 060
9001975



Fig 77

● HP 239 029
9001970

● HP 239 060
9001976



Fig 78

● HP 257 040
9001974

● HP 257 060
9001973



Fig 79

● HP 194 040
9001959

● HP 194 045
9001971

● HP 194 050
9001857

● HP 194 060
9001972



Fig 129

● HP 141 023
9001961



Fig 136

● HP 197 014
9001968



Fig 138

● HP 198 016
9001967



Fig 139

● HP 289 023
9001965



Fig 251

● HP 274 060
9001977



Fig 257

● HP 187 023
9001963



Fig 261

● HP 194 023
9001960



Fig 295

● HP 292 023
9001962



Fig 351

● HP 263 060
9001980

Rotary Instruments

TC Cutters

HS Tungsten Carbide Cutters

Cross Cut 140

- Cross Cut Fine Red for acrylics, metal alloys and precious metals.



Fig 73

● HP 277 014
9001938

● HP 277 023
9001939



Fig 77

● HP 237 023
9001940



Fig 77

● HP 239 060
9001955



Fig 78

● HP 257 012
9001941

● HP 257 060
9001957



Fig 79

● HP 194 031
9001953

● HP 194 040
9001944

● HP 194 045
9001954



Fig 129

● HP 141 023
9001951



Fig 137

● HP 225 023
9001942



Fig 138

● HP 198 016
9001952

● HP 198 023
9001949

● HP 198 040
9001860



Fig 250

● HP 277 040
9001861



Fig 251

● HP 274 060
9001958



Fig 257

● HP 187 023
9001943



Fig 257R

● HP 200 023
9001859



Fig 261

● HP 194 023
9001945



Fig 283

● HP 289 010
9001858



Fig 295

● HP 292 023
9001947



Fig 364R

● HP 144 023
9001946



Fig 351

● HP 263 060
9001956

Rotary Instruments

TC Cutters

HS Tungsten Carbide Cutters

Cross Cut 110

- Cross Cut Extra Fine Yellow for metal alloys, precious metals and ceramics.



Fig 77
● HP 237 023
9001987



Fig 79
● HP 194 040
9001983



Fig 139
● HP 289 023
9001986



Fig 364R
● HP 144 023
9001985

Cross Cut Smooth Cut 138

- Cross Cut Fine Smooth Cut Red for metal alloys.



Fig 79
● HP 194 040
9001862



Fig 79
● HP 194 070
9001934



Fig 251
● HP 274 060
9001935



Fig 261
● HP 194 023
9001933

Helical Cross Cut Fine 134

- Helical Cross Cut Fine Violet for acrylics, metal alloys, precious metals and titanium.



Fig 73
● HP 277 014
9001863



Fig 79
● HP 194 040
9001990



Fig 251
● HP 274 060
9001991



Fig 257R
● HP 200 023
9001867



Fig 351
● HP 198 040
9001868



Fig 79
● HP 194 045
9001930
● HP 194 060
9001932



Fig 364R
● HP 137 010
9001993
● HP 137 023
9001995
● HP 137 015
9001994

A

Abrasives 115

Addition Cured Silicone 4

Alginate Adhesive 8

Alginate – Krom..... 4

Alginate Mixer 104

Alginate Plus 4

Alloys & Accessories 40

Air/Water Syringes 10

Air Polishing Powder 58

Air Polishing Powder – Lemon..... 58

Amalgam Non Gamma 2 Caps 40

Amalgam Pluggers 65

Amalgam Polishing Kit..... 53

Anesthetics & Needles..... 1

Articulating Forceps 68

Articulating Paper 39

Aspirator Tips 10

Aspirator Tubes..... 10

Autoclavable Aspirator Tips..... 10

B

Ball Burnishers 65

Barbed Broaches..... 96

Barrier Envelopes..... 110

Barrier Products 18

Barrier Sleeves..... 18

Bibs & Aprons 13

Bibs & Aprons Accessories 13

Bite Registration 6

Bite Wax 116

Bitewax Bars..... 116

Bitewax Rims..... 116

Blood Pressure Monitors..... 32

Bond..... 39

Bouff ant Cap..... 32

Bowie Dick Test..... 27

Brushes for Separation Liquid..... 118

Brush for Wax 118

Bulk Fill Composites..... 40

Bunsen Burner Gas..... 118

Bur Brushes..... 93

Bur Disinfection 21

Burnishers..... 65

C

C-Fold Towels 12

CAD/CAM Accessories 112

CAD/CAM Chairside Finishing and Polishing..... 112

CAD/CAM Lab Finishing and Polishing 112

CAD/CAM Materials..... 112

CAD/CAM Scan Spray 113

Calcium Hydroxide Baseliners..... 34

Calcium Hydroxide Placements..... 69

Carbide Trimming & Finishing Burs..... 91

Caries Indicators 56

Carvers..... 65

Casting..... 112

Casting Rings..... 112

Casting Wax 116

Cement Mixing Tips - Brown 35

Cement Spatula 66

Cements & Liners..... 33

Ceramic Brushes..... 119

Ceramic Instrument 119

Ceramics Knife..... 119

Ceramist Synthetic Brushes 119

Clamps 101

Class 6 Indicator Tests Strips..... 27

Colour Coded Thomas Instruments 120

Compomers 41

Composite Instruments..... 65

Compule Dispenser 44

Concentrated Developer/Fixer..... 109

Condensation Silicone..... 5

Consumables..... 9

Contour Bibs 13

Coping & Spacer Foils..... 118

Cord Packing Instruments..... 69

Coreless Centrefeed Towel Rolls 12

Coreless Centrefeed Towel Rolls - Mini 12

Cotton Pellets..... 11

Cotton Products..... 11

Cotton Roll Dispenser 14

Cotton Rolls..... 11

Cotton Wool Balls..... 11

Cotton/Gauze Balls..... 11

Criterion Surgeons' Latex Powder-Free Gloves - Sterile..... 20

Criterion Vinyl Powder-Free Gloves 19

Crown & Bridge 47

Crown and Bridge Veneer Instrument..... 120

Curing Lights 105

Cutting & Separating Discs 121

D

Dam Frames & Instruments.....	102
Dental Dam.....	102
Dental Needles.....	2
Dentine Pins.....	46
Dentipins.....	46
Denture & Retainer Boxes.....	15
Denture Adjustment.....	52
Denture Box.....	15
Denture Brush.....	60
Denture Finish Spearmint.....	116
Desensitisers.....	56
Desensitising Mousse.....	56
Diagnostic Instruments.....	62
Diagnostic Instrument Sets.....	94
Diamond Burs.....	72
Diamond Polishing Paste.....	53
Diamond Tipped Forceps.....	119
Diamond Tipped Jacket Grip Clamp.....	119
Die Spacer.....	114
Digital Blood Pressure Unit.....	32
Digital Ear Thermometers.....	32
Dispensers.....	14
Dispensing Gun.....	6
Disposable Aspirator Tubes.....	10
Disposable Diamond Burs - Non-Sterile.....	86
Disposable Diamond Burs - Sterile.....	85
Disposable Impression Trays.....	7
Disposable Plastic Impression Trays.....	7
Disposable Prophy Angles – Oscillating.....	57
Disposable Pulp Kidney Dish.....	15
Disposable Saliva Ejectors.....	10
Disposable Scalpels - Sterile.....	30
Disposable Suction Filter.....	15
Disposable Toothbrushes with Paste.....	60
Disposable Tray Liners.....	12
Double Tray Paper.....	12
Dowel Pins.....	114
Drapes, Gowns & Accessories.....	31
Dri-Gards.....	13
Drills & Reamers.....	96
Drinking Cups.....	13
Dry Cleaning Cloth.....	12
Duplication Material.....	114
Dynamic Mixing Tips.....	6
DX-54 Child Periapical Film.....	108
DX-58 Adult Periapical Film.....	108

E

Earloop Face Masks.....	18
Economy Non-Sterile Latex Gloves Powder-Free.....	20
EDTA.....	100
Electrolyt Solution.....	113
Elevators.....	67
Endo Access Burs.....	92
Endo Accessories.....	102
Endo Cold Spray.....	102
Endo Organiser.....	102
Endo Organiser Ring.....	102
Endo Safe-End Access Burs.....	92
Endo Sponges.....	102
Endodontic Blunt Needle.....	101
Endodontic Irrigating Needle.....	101
Endodontics.....	95
Endoflex Hedstroem Files.....	97
Endoflex K-Files.....	96
Endoflex Reamers.....	97
Etch Gel.....	38
Etch Gel for Ceramics.....	38
EuroGauze Balls.....	11
EuroGauze Swabs - Non-Sterile.....	11
EuroSept Max Instrument Lemon.....	21
EuroSept Plus Alcohol-Free XL Wipes.....	24
EuroSept Plus Bur Disinfection.....	21
EuroSept Plus Hand Cream.....	27
EuroSept Plus Hand Cream Lotion.....	27
EuroSept Plus Hand Disinfection Liquid.....	26
EuroSept Plus Hand Disinfection Gel.....	26
EuroSept Plus Impression Disinfection.....	22
EuroSept Plus Instrument Disinfection.....	22
EuroSept Plus Instrument Lemon.....	21
EuroSept Plus Instrument Cleaner Enzymatic.....	22
EuroSept Plus Instrument Disinfection FV.....	22
EuroSept Plus Spittoon Disinfection.....	21
EuroSept Plus Maxima Spittoon Cleaner.....	21
EuroSept Plus Orange Oil Cleaner/Solvent.....	8, 25
EuroSept Plus Surface.....	23
EuroSept Plus Surface Alcohol-Free.....	23
EuroSept Plus Surface Celaner Limescale.....	21
EuroSept Plus Surface Disinfection FV.....	23
EuroSept Plus Surface Disinfection Alcohol Free Concentrate.....	23
EuroSept Plus Surface Lemon.....	24
EuroSept Plus Surface Non-Alcohol Concentrate.....	23
EuroSept Plus System Alcohol-Free Wet Wipes.....	25
EuroSept Plus Thermotabs.....	25
EuroSept Plus Tray Cleaner Concentrate.....	8, 25
EuroSept Plus Wash Lotion.....	26
EuroSept Plus Wash Lotion Sensitive.....	26
EuroSept Plus Wipes/XL Wipes.....	24

EuroSept Plus Wipes Alcohol-Free	24
EuroSept Plus Wash Lotion	26
EuroSept Plus Wash Lotion Sensitive	26
EuroSept Plus XL	24
Eurosil Max 2 Catalyst	5
Eurosil Max 2 Correction	5
Eurosil Max 2 Lab Putty	114
Eurosil Max 2 Putty	5
Evacuation	10
Excavators	66
Explorers/Probes	62
Extra Protection Long Nitrile Gloves	20
Extraction Forceps	67
Extraoral Films	108
Eye Protection	18

F

Fabric Polishing Discs	116
Face Masks - Earloop	18
Face Masks - Tie-On	19
Faceshield Visor Frames	18
Faceshield Visor Refills	18
Facial Tissues	12
Fibre Posts	46
Finger Pluggers	98
Finger Pulse Oximeters	32
Finishing & Polishing	51
Finishing & Polishing - Laboratory	116
Finishing and Polishing Discs	52
Finishing and Polishing Discs Extra Thin	52
Finishing and Polishing Strips	53
Finishing Burs	91
Flosser	59
Flosser - Kids	59
Floss Patient	59
Flowable Composite Tips	44
Fluoride Treatment	56
Fluoride Varnish	56
Foil Scissors	120
Forceps	68
Furnace Tong	121
FX-58 Adult Periapical Film	108

G

Gates Glidden Drills	96
Gauze	11
General Disposables	15

General Instruments	118
General Lab Items	117
Glass Ionomer Baseliners	34
Glass Ionomer Filling Materials	41
Glass Ionomer Luting Cements	34
Glove Dispenser	14
Gloves - Latex Powder-Free	20
Gloves - Latex Surgical	20
Gloves - Nitrile	20
Gloves - Nitrile Extra Protection Long	20
Gloves - Vinyl Powder-Free	19
Glue	117
Gold Aesthetic Finishing Burs OF Series	91
Gold Aesthetic Finishing Burs Safe-End EF Series	91
Gold Fine Finishing Burs Safe-End EF Series	92
Gold Trimming & Finishing Burs - 12 Bladed	92
GP Auxiliary Points	99
GP Points	99
Gracey Curettes	63
Greater Taper GP Points	100
Gutta Percha Points	99

H

Haemostatic Agents	48
Haemostatic Gel	48
Hand Cleaning & Disinfection	26
Hand Instruments	61
Hand Reamers	97
Hand Towel Dispensers	14
Hand Washes & Lotions	27
Handles	30
Handpiece Maintenance	106
Headrest Covers	13
Hedstroem Files	97
Helix Test	27
High Volume Evacuators	10
Hollenbach Carver	65
Hybrid Composites	40
Hypodermic Needles	2

I

Implant Cements	36
Implant Scalers	64
Impression Cleaning & Instrument Disinfection	22
Impression Materials	3
Impression Tray Resin LC	115
Impression Trays	7

Impression Trays with Integral Handle	7
Infection Control.....	17
Instant Fix.....	117
Instrument & Bur Cleaning & Disinfection.....	21
Instrument Cassettes.....	28, 70
Instrument Cleaner.....	21
Instrument Trays & Accessories.....	70
Instruments	61
Instruments - Lab 119	
Interdental Brush 59	
Interdental Wedges	44
Interfold Towels.....	12
Intraoral Films	108
Investment Material.....	113
Irrigation Syringes & Needles.....	100, 101
Iwanson Measuring Device	66, 121

J

Jumbo Etch.....	38
-----------------	----

K

K-Files.....	96
Krom-Alginate	4
Kid Flosser	59
Kidney Dish.....	15

L

Lab Equipment Accessories.....	121
Lab Instrument Set.....	120
Lab Instruments	118
Lab Instruments - Knives	119, 121
Lab Putty	114
Lab Small Equipment.....	121
Lab Tweezers.....	120
Laboratory Products.....	111
Latex Rubber Dam	102
LeCron Carver	66
LED 1200 Cordless Curing Light.....	105
LED 1500 Cordless Curing Light.....	105
LGI Plus Glass Ionomer Luting Cement.....	34
Light-Cured Composites.....	40
Limescale.....	21
Luting Cements	34

M

Masks	18
Mathieu Tie Pliers.....	68
Matrix Retainers & Bands	43
Matrix Strips.....	43
Max-Filltemp.....	36
MaxiGrip Ergo Gracey.....	63
MaxiGrip Ergo Scalars.....	64
Maxima Blood Pressure Monitor	32
Maxima Burs.....	72
Maxima Cotton Roll Dispenser.....	14
Maxima Cotton Rolls.....	11
Maxima Diamond Burs	72
Maxima Diamond Burs Amalgam Buster.....	81
Maxima Diamond Burs Barrel	79
Maxima Diamond Burs Bevelled Cylinder	76
Maxima Diamond Burs Curettage.....	79
Maxima Diamond Burs Depth Marker.....	81
Maxima Diamond Burs End Cutting	81
Maxima Diamond Burs Flame	78
Maxima Diamond Burs Flat End Cylinder	75
Maxima Diamond Burs Flat End Taper.....	74
Maxima Diamond Burs Inverted Cone.....	72
Maxima Diamond Burs Inverted Cone with Collar	72
Maxima Diamond Burs Long Pear.....	81
Maxima Diamond Burs Modified Cylinder.....	76
Maxima Diamond Burs Modified Taper	74
Maxima Diamond Burs Needle.....	80
Maxima Diamond Burs Pear	81
Maxima Diamond Burs Round	72
Maxima Diamond Burs Round End Cylinder	75
Maxima Diamond Burs Round End Taper.....	73
Maxima Diamond Burs Tissue Protective End	81
Maxima Diamond Burs Xmas Tree	80
Maxima Earloop Masks	18
Maxima Earloop Masks - Extra Protection.....	19
Maxima Earloop Masks - Sensitive.....	18
Maxima Earloop Masks with Eyeshield	19
Maxima Etch Gel.....	38
Maxima Fibre Posts	46
Maxima Gutta Percha Points.....	99
Maxima H-Files.....	16
Maxima K-Files	16
Maxima Miniature Burs.....	82
Maxima Natural Elegance 7th Generation Single Bond	39
Maxima Natural Elegance NH Composite.....	40
Maxima Natural Elegance NH Flow Composite.....	41
Maxima Needle Holders.....	69
Maxima Oil Spray.....	106
Maxima One Piece Carbide Burs.....	90
Maxima Paper Points	99

Maxima Tie-On Surgical Masks	19
Maxima Turbo Diamond Burs.....	82
MaxiTemp 1:1.....	48
MaxiTemp.....	48
Medicaments.....	100
Metal Trays.....	7
Microhybrid Composites.....	41
Microhybrid Flowable Composites.....	41
Micromotor L50	104
Micromotors	104
Mirror Handles.....	62
Mirrors	62
Miscellaneous.....	68, 116
Mixing Accessories	44
Mixing Bowl.....	44
Mixing Guns & Tips	6
Mixing Pads	36
Model Formers.....	117
Model Plaster & Stone.....	115
Modelling Wax	116
Mouth Guard & Splint Material.....	117
Mouth Mirrors.....	62
Mouthwash Alcohol-Free.....	56
Mouthwash & Tablets	56
Mouthwash Tablets	56
MTA Bioceramic Root Canal Sealer.....	100
Multi-Trays.....	7

N

Nanohybrid Composites	40
Nanohybrid Flowable Composites	41
Natural Elegance Bleach	60
Natural Elegance Bulk Fill Composite.....	40
Natural Elegance Flowable Composite.....	41
Natural Elegance Flowable Composite Tips	44
Natural Elegance Microhybrid Composite.....	41
Natural Elegance Sealant Kit	56
Natural Elegance Self-Etch Bond 6th Generation.....	39
Natural Elegance Single Coat Bond 5th Generation.....	39
Needles.....	2
Needle Holders	69
NiTi Files	98
Nitrile Gloves Powder-Free	20
Non Latex Dam.....	102
NT Files	98

O

Occlusion Spray	43
Oil Cleaner/Solvent.....	8, 25
Oral Hygiene & Whitening	55
Oral Surgery Burs - Sterile.....	93
Orange Solvents/Cleaners	8, 25
Orthodontic Debonding Burs	91
Orthodontic Retainer Boxes	15

P

Pan-G Extraoral Film.....	108
Paper Points.....	99
Paper Products	12
Part Vented Evacuator	10
Paste Fillers	98
Patient Apron.....	13
Patient Floss.....	59
Patient Monitoring	32
Patient Tape.....	59
Peeso Reamers RA	96
Perio Probes.....	63
Periodontal Instruments.....	63
Periosteal Elevators	67
Permanent Cements.....	35
PGA Rapid	30
PGA Violet.....	30
Pins & Posts.....	45
Pit & Fissure Sealants	56
PK Thomas Carver.....	120
Plaster Separation Liquid.....	115
Plastic Filling Instruments.....	66
Plastic Hub Needles	2
Plastic Tray Inserts	15
Plasticware, Dispensers & Containers	15
Pliers	69
Polishing Brushes	116
Polishing Discs.....	52
Polishing Pastes	53
Polishing Pencil.....	116
Polishing Points	53
Polishing Strips	53
Polycarbonate Crowns.....	49
Pouches & Cassettes	28
Preline Plus.....	34
Premium Mouth Mirrors.....	62
Premium Non-Sterile Latex Gloves Powder-Free	20
Pressure Indicators	43
Probes/Explorers	62
Profile Wax Round.....	113

Prophy Angles	57
Prophy Brushes	57
Prophy Brushes - Latex-Free	57
Prophy Cups	58
Prophy Cups - Latex-Free	58
Prophy Cups RA	58
Prophy Cups Screw Type	58
Prophy Paste	58
Prophy Powder	58
Prophylaxis	57
Prosthetics	115
Protective Sleeves	18
Pulse Oximeters	32
Pumice Mouse	115
Pumice Powder	115

R

Radiomat XR Developer and Fixer	109
Rebasing Jig 3S	118
Red Sable Brush	119
Reammers	46
Resin Cements	35
Resin Mirror Handles	62
Resin Separation Liquid	115
Restorative Instruments	65
Restorative Mixing Guns & Tips	44
Restoratives	37
Retainer Boxes	15
Retraction Cord	48
Retractors	68
RGI Glass Ionomer Baseline	34
RGI Glass Ionomer Fast 25 Capsules	42
RGI Glass Ionomer Filling Material	42
RGI Glass Ionomer XFast 20 Capsules	42
RGI Glass Ionomer Regular 3 5 Capsules	43
RGI Resin Reinforced Glass Ionomer	42
RGI Silver Reinforced Glass Ionomer	43
Root Canal Sealers	100
Rotary Instruments	71
Rotary NiTi Files	98
Round Tips	44
Round Tips Dispenser	14
Rubber Dam	102
Rubber Dam Clamp Forceps	102
Rubber Dam Clamps	101
Rubber Dam Frame	102
Rubber Dam Punch Ainsworth	102

S

S-Spray	106
Saliva Ejectors	10
Scalers	64
Scalpels	30
Scalpels Blades	30
Scalpels Handles	30
Scissors	69
Scrub Brush	28
Self-Seal Sterilisation Pouches	28
Separating Discs	121
Separating Liquids	115
Sharps Container	2
Shoe Covers - Latex-Free	32
Silicone Spray	118
Silk Black Sutures	31
Single Use Diamond Burs	86
Small Equipment	103
Sodium Hypochlorite	100
Solvents & Cleaners	8
Specialty Brush For Wax	118
Spittoon Cleaners	21
Sprue Formers	113
Stainless Steel Impression Trays	7
Stainless Steel Instrument Trays	70
Standard Diamond Burs	83
Steel Burs	87
Sterile Blades	30
Sterile Drape Kit	31
Sterile Luer Lock Syringe	100
Sterile Luer Syringe - 2 Part	31
Sterile Scalpels	30
Sterile Single Use Diamond Burs	86
Sterile Surgical Aspirators	10
Sterile Surgical Gowns	31
Sterireels	28
Stratosphere 2	40
Suction Filter	15
Surface Barriers	18
Surface Disinfection	23
Surface Limescale	21
Surgical	29
Surgical Aspirator Tip	10
Surgical Drape	31

Surgical Drape with Adhesive Strip32

Surgical Instruments..... 67

Surgical Masks 19

Sutures.....30

Syringe Tips with Metal Insert - 3 in 110

Syringe Tips with Plastic Insert - 3 in 110

Syringe Tips without Insert - 3 in 110

Syringes..... 2

System Wipes25

T

T-Mixer 1:1 Tips..... 6

T-Mixer Dispensing Gun 6

TC Cutters.....122

Temporary Cements 35

Temporary Crown & Bridge Materials.....48

Temporary Crowns 49

Temporary Filling Materials 36

Temporary Root Canal Filling Material.....34

Testing.....27

Thermo Bleach Round..... 117

Thermo Guard117

Thermo Splint.....117

Thyroid Collar108

Tie-On Face Mask.....19

Tissue Forceps..... 68

Tissue Retractor68

Titanium Implant Scalars..... 64

Tofflemire Bands..... 43

Tofflemire Retainer Bands..... 43

Toilet Rolls12

Toothbrushes60

Toothbrushes Gentle Massage.....60

Towel Roll Dispensers..... 14

Transparent Mixing Pads..... 36

Tray Adhesives..... 8

Tray Cleaners.....8, 25

Tray Cure121

Tray Liners.....12

Tray Material.....115

Tray Paper12

Trays & Inserts15

Tungsten Carbide - Sterile Oral Surgery93

Tungsten Carbide - Surgical Burs93

Tungsten Carbide Burs 88

Tungsten Carbide Crown Cutter Burs90

Tungsten Carbide Cutters.....90

Turbo Diamond Burs 82

Tweezers & Pliers69

U

Universal Silicone Polishing Kit.....53

Utility Syringe 100

V

Vest..... 113

VPS Hydro 380..... 5

VPS Hydro 380 Accesories..... 6

VPS Premium Hydro Bite..... 6

VPS Hydro Putty..... 4

VPS Hydro Wash Light Body..... 4

VPS-CADCAM..... 6

W

Wash Lotions.....26

Washer Disinfection..... 25

Wax Bite Rims116

Wax Carver66

Wax Clasps116

Wax Instruments 66, 120

Wax Knife.....121

Waxwire113

Whitening Tray Box..... 60

Whitening Treatment & Accessories60

Wipes.....24

X

X-Ray & Imaging..... 107

X-Ray Accessories..... 110

X-Ray Chemicals.....109

X-Ray Cleaning Films 110

X-Ray Developer and Fixer.....109

X-Ray Protection.....108

X-Ray Viewers.....106

XRV 30.....106

Z

Z Dupe Silicone Flask 114

Notes

Notes

**BECAUSE SO MANY RELY ON YOU...
YOU CAN RELY ON US.**

Europe and South East Asia: **+49 1522 1748061 • david.burckhardt@henryschein.de**
Middle East: **+971 50 4813292 • ghassan.nasser@henryschein.com**