

# PRODUCT OVERVIEW

## EDGEGLIDEPATH™

4 files in a box



The EdgeGlidePath™ heat-treated FireWire NiTi file features the flexibility and strength that EdgeEndo is known for, replacing multiple instruments to create a glidepath faster with far less effort and with a 1mm maximum flute diameter.

|                                                                                   |  | Tip size | Length |         |         |         |
|-----------------------------------------------------------------------------------|--|----------|--------|---------|---------|---------|
|  |  | 19       | 19mm   | 21mm    | 25mm    | 31mm    |
|                                                                                   |  |          |        | 9887892 | 9887893 | 9887894 |

## EDGETAPER™

6 files in a pack



EdgeTaper™ files are non-heat-treated rotary files, available in 21mm, 25mm and 31mm lengths.

|                                                                                     |                                    |          | Length  |         |         |         |
|-------------------------------------------------------------------------------------|------------------------------------|----------|---------|---------|---------|---------|
|                                                                                     |                                    | Tip size | 19mm    | 21mm    | 25mm    | 31mm    |
|    | SX                                 | 19       | 9887938 |         |         |         |
|    | S1                                 | 17       |         | 9887932 | 9887933 | 9887934 |
|    | S2                                 | 20       |         | 9887935 | 9887936 | 9887937 |
|   | F1                                 | 20       |         | 9887917 | 9887918 | 9887919 |
|  | F2                                 | 25       |         | 9887920 | 9887921 | 9887922 |
|  | F3                                 | 30       |         | 9887923 | 9887924 | 9887925 |
|  | F4                                 | 40       |         | 9887926 | 9887927 | 9887928 |
|  | F5                                 | 50       |         | 9887929 | 9887930 | 9887931 |
|  | Assorted<br>SX, S1, S2, F1, F2, F3 |          |         | 9888005 | 9888006 | 9888035 |

## EDGETAPER PLATINUM™

6 files in a pack

The EdgeTaper Platinum™ files are heat-treated rotary files, available in 21mm, 25mm and 31mm lengths. EdgeTaper Platinum™ features our heat-treated Firewire™ NiTi which increases its resistance to cyclical fatigue.

|                                                                                     |                                    |          | Length  |         |         |         |
|-------------------------------------------------------------------------------------|------------------------------------|----------|---------|---------|---------|---------|
|                                                                                     |                                    | Tip size | 19mm    | 21mm    | 25mm    | 31mm    |
|  | SX                                 | 19       | 9887916 |         |         |         |
|  | S1                                 | 17       |         | 9887910 | 9887911 | 9887912 |
|  | S2                                 | 20       |         | 9887913 | 9887914 | 9887915 |
|  | F1                                 | 20       |         | 9887895 | 9887896 | 9887897 |
|  | F2                                 | 25       |         | 9887898 | 9887899 | 9887900 |
|  | F3                                 | 30       |         | 9887901 | 9887902 | 9887903 |
|  | F4                                 | 40       |         | 9887904 | 9887905 | 9887906 |
|  | F5                                 | 50       |         | 9887907 | 9887908 | 9887909 |
|  | Assorted<br>SX, S1, S2, F1, F2, F3 |          |         | 9888000 | 9888004 | 9888034 |





# PRODUCT OVERVIEW

| EDGEFILL™ FOR EDGETAPER™                                                          |    |     |          | 6 pcs in a box |
|-----------------------------------------------------------------------------------|----|-----|----------|----------------|
| Plastic carrier thermal obturator for use with EdgeTaper™ and EdgeTaper Platinum™ |    | ISO | Order No |                |
|  | F1 | 20  | 9887974  |                |
|  | F2 | 25  | 9887976  |                |
|  | F3 | 30  | 9887978  |                |
|  | F4 | 40  | 9887980  |                |
|  | F5 | 50  | 9887982  |                |





| EDGECORE™ FOR EDGETAPER™                                                                 |    |     |          | 5 pcs & 1 verifier in a box |
|------------------------------------------------------------------------------------------|----|-----|----------|-----------------------------|
| Pink gutta percha core thermal obturator for use with EdgeTaper™ and EdgeTaper Platinum™ |    | ISO | Order No |                             |
|         | F1 | 20  | 9887973  |                             |
|         | F2 | 25  | 9887975  |                             |
|       | F3 | 30  | 9887977  |                             |
|       | F4 | 40  | 9887979  |                             |
|       | F5 | 50  | 9887981  |                             |



| EDGETAPER™ GUTTA POINTS                                                             |    |     |          | 60 pcs in a box |
|-------------------------------------------------------------------------------------|----|-----|----------|-----------------|
| Gutta percha points compatible with EdgeTaper™ and EdgeTaper Platinum™              |    | ISO | Order No |                 |
|  | F1 | 20  | 9887986  |                 |
|  | F2 | 25  | 9887987  |                 |
|  | F3 | 30  | 9887988  |                 |
|  | F4 | 40  | 9887989  |                 |
|  | F5 | 50  | 9887990  |                 |



| EDGETAPER™ PAPER POINTS                                                             |    |     |          | 60 pcs in a box |
|-------------------------------------------------------------------------------------|----|-----|----------|-----------------|
| Paper points compatible with EdgeTaper™ and EdgeTaper Platinum™                     |    | ISO | Order No |                 |
|  | F1 | 20  | 9887994  |                 |
|  | F2 | 25  | 9887995  |                 |
|  | F3 | 30  | 9887996  |                 |
|  | F4 | 40  | 9887997  |                 |
|  | F5 | 50  | 9887998  |                 |

