



***Minimally Invasive Dentistry
with Proven, Patient-Friendly
Results***

waterlase*MDX™

Waterlase is the all-tissue laser trusted by the most dentists worldwide

Over 28 million patients all over the world have experienced minimally invasive Waterlase Dentistry™. The Waterlase MDX is an cost-effective way to add Waterlase Dentistry to your practice. Waterlase **YSGG** technology has been proven to provide the best clinical & patient reported outcome in perio procedures with minimum discomfort, swelling and bleeding, backed by a landmark, randomized, multi-center, blinded study. Perform **80+ FDA cleared indications** including fast and effective crown/veneer removal.

	Laser Medium
	SureFire™ Trunk fibers
	Maximum Power
	Maximum Energy
	Graphical User Interface
	Fractional Handpiece
	Colors

Waterlase iPlus	Waterlase MDX300
YSGG	YSGG
✓	✓
10 W	8 W
600 mJ	300 mJ
✓	✗
✓	✗
Platinum Racing Red	Metallic Only

INTUITIVE INTERFACE




Easy to use interface with presets make operation efficient and quick to learn!

MAXIMUM FLEXIBILITY



Extremely light and flexible fiber optic cable help alleviate user fatigue.

MAXIMUM VERSATILITY



The ergonomic, contra-angle Waterlase Gold Handpiece is designed for maximum versatility and access, compared to bulky articulating arms.

MDX SPECIFICATIONS

Wavelength: Er,Cr:YSGG, 2780 nm

Power: 0.1 to 8.0 W

Pulse Repetition Rates: 10-50 Hz Clinician Selectable

Pulse Energy: 0-300mJ

Laser Classification: 4

Operating Voltage: 100-230 VAC

Dimensions: 10.5 x 19 x 32 in (27 x 48 x 91 cm)

Weight: 75lbs (34 kg)

Visit biolase.com for more information.

BIOLASE
27042 Towne Centre Drive, Suite 270
Foothill Ranch, CA 92610 USA
+1 949.361.1200 + 833.BIOLASE
biolase.com

EU REPRESENTATIVE
MT Promedt Consulting GmbH
Ernst-Heckel-Straße 7
66386 St. Ingbert, Germany
+49 6894 581020 mt-procons.com

Rx Only **MD**

BIOLASE
LEADER IN DENTAL LASERS™

biolase.com