

Smart and economic choice!

MEG-ENGINE ECO



Touch once for automatic deep cleanse of contra angle!

Automatic cleansing function enables the user to maintain the handpiece performance for a longer period of time

Change procedure with one touch of the screen

Implant and surgical procedure parameters can be freely set with the touch of a button!



Built-in Fiber Optic LED!

Perform your operation with a built-in fiber optic LED Experience the difference with MEG-ENGINE ECO's two-level lighting system!

Ultra light handle!

Provides a lighter, comfortable handle grip than previous models—more durable, but 30% less handle weight due to use of new materials.



	Other brands	MEG-ENGINE ECO
Maximum decibel	56.81dB	44.48dB
Rise in temperature (After 5 minute use)	18°C	8°C


4th generation motor

Perform your operation with a quieter 4th generation motor that does not overheat!

MEG-ENGINE ECO components



Handpiece



Bracket



Micro Motor



Handpiece holder



Foot pedal



Motor



Irrigation Tube (5ea)
*disposable

Technical data

Coolant flow rate at 100%	110 ml/min (at least)
Voltage	250 V - F 3.15 AL
Maximum Torque	70 Ncm
Motor Speed range	300 - 40,000 rpm
Dimensions (mm)	230mm(W) x 262mm(L) x 125mm(H)
Weight (kg)	1.9
Operation control	Touch Screen
Optic system	O
PRE-SET Surgical Protocol	O
Data Storage Mode	X

Product Name	Ref.C
Set MEG-ENGINE ECO, 230V	COXO-ENG-ECO
C6-22 contra angle 20:1	COXO-CX235C6
M1310 motor with cable 1.8m	COXO-ENG-MOTOR
Multifunctional Foot Control	COXO-ENG-PEDAL
Irrigation tubing set 2.2m disp osable set(5p cs) :1회용	COXO-ENG-TUBE
MEG ENGINE ECO 20: 1 contra angle handpiece cartridge	COXO-ENG-HP-CART
MEG ENGINE ECO 20: 1 contra angle handpiece shaft	COXO-ENG-HP-SHAFT
MEG ENGINE ECO power cord	COXO-ENG-POWER
MEG ENGINE ECO fuse	COXO-ENG-FUSE
MEG ENGINE ECO motherboard circuit	COXO-ENG-CIRC

* Products with a REF number can be purchased separately; however, the bracket and handpiece holder cannot be purchased separately.