

**DEVELOPED AND
RECOMMENDED BY
MEDICAL SPECIALISTS**

Lumoral is a CE-marked medical device developed by a team of Finnish dentists, doctors and researchers. This innovative treatment has undergone extensive clinical research and testing, consistently delivering outstanding results.



Manufacturer: Koite Health Ltd, Karaportti 5, 02610, Espoo, Finland
+358 50 306 5452 info@lumoral.com

Lumoral is a CE marked medical device – CE 0598

Nikinmaa S, Moilanen N, Sorsa T, Rantala J, Alapulli H, Kotiranta A, Auvinen P, Kankuri E, Meurman JH, Pättilä T. Indocyanine Green-Assisted and LED-Light-Activated Antibacterial Photodynamic Therapy Reduces Dental Plaque. Dent J (Basel). 2021 May 3;9(5):52. doi: 10.3390/dj9050052. PMID: 34063662; PMCID: PMC8147628.

Nikinmaa S, Alapulli H, Auvinen P, Vaara M, Rantala J, Kankuri E, et al. (2020) Dual-light photodynamic therapy administered daily provides a sustained antibacterial effect on biofilm and prevents Streptococcus mutans adaptation. PLoS ONE 15(5): e0232775. <https://doi.org/10.1371/journal.pone.0232775>

Nikinmaa S.; Podonyi, A.; Raivio, P.; Meurman, J.; Sorsa, T.; Rantala, J.; Kankuri, E.; Tauriainen, T.; Pättilä, T. Daily Administered Dual-Light Photodynamic Therapy Provides a Sustained Antibacterial Effect on Staphylococcus aureus. Antibiotics 2021, 10, 1240. <https://doi.org/10.3390/antibiotics10101240>

Henttilä, J. ., Laakamaa, N., Sorsa, T., Meurman, J., Välimaa, H., Nikinmaa, S., Kankuri, E., Tauriainen, T., & Pättilä, T. (2021). Dual-Light Photodynamic Therapy Effectively Eliminates Streptococcus Oral Biofilms: Dual-light therapy, Photodynamic therapy, Streptococcus oralis, soft-tissue infections, alternative antimicrobial therapy. . Journal of Pharmacy & Pharmaceutical Sciences, 24, 484–487. <https://doi.org/10.18433/jpps32084>

Dörfer C, Benz C, Aida J, Campard G. The relationship of oral health with general health and NCDs: a brief review Int Dent J 2017 Sep 67 Suppl 2:14–18. PMID: 29023744; DOI: 10.1111/rdj.12360

Suunterveyden merkitys yleisterveydelle | Hammaslääkäriliitto (hammaslaakariliitto.fi)

Lumoral.

**Clean teeth,
healthy gums and
fresh breath**



lumoral.com

lumoral.com

lumoral.com

ELEVATE YOUR ORAL HEALTH TO NEW STANDARDS WITH LUMORAL

Lumoral is the first medical solution that selectively targets its antibacterial effect on dental plaque, protecting the healthy bacterial balance in the mouth. The treatment reduces inflammation in the gums and prevents decay and tartar build-up.

When used regularly as an addition to traditional oral hygiene, the antibacterial method has been shown to improve oral health and periodontal treatment outcomes.

ORAL HEALTH'S CRITICAL ROLE IN OVERALL HEALTH

Oral health has a huge impact on overall health. Untreated oral inflammation can, at its worst, put you at risk of cardiovascular disease, as bacteria in your mouth can easily spread to the rest of your body. You will benefit from enhanced oral health self-care, especially if you have pre-existing dental problems or a general medical condition (e.g. diabetes, rheumatic diseases, cardiovascular diseases).

LUMORAL'S TARGETED ACTION ELIMINATES PLAQUE BACTERIA BY OVER 99.99%



PREVENT ORAL DISEASES CAUSED BY BACTERIA

Lumoral treatment

- ✓ reduces tartar build-up
- ✓ prevents tooth decay
- ✓ prevents gingivitis
- ✓ treats irritated and inflamed gums
- ✓ freshens breath
- ✓ maintains a healthy bacterial balance in the mouth
- ✓ removes stains

RESULTS IN AS LITTLE AS 10 MINUTES TWICE A WEEK



1. SWISH

Swish your mouth with the minty Lumorinse® mouth rinse for 60 seconds. The mouth rinse's light-sensitive molecule adheres to the plaque surface.



2. ACTIVATE

Place the mouthpiece in your mouth and apply light with Lumoral for 10 minutes. The light activates the mouth rinse, which triggers the antibacterial effect of the treatment.



3. BRUSH

Brush your teeth and clean the interdental spaces carefully. The treatment inactivates plaque bacteria, leaving your teeth amazingly clean and smooth!

**You will notice a difference in the cleanliness
of your teeth after the first use.**